

ASSEMBLY INSTRUCTION

External Joystick / External Joystick Lite



Description

Article number

External Joystick - Mounting on armrest, opposite side as Omni display

☐

1016504

External Joystick Lite - Mounting on armrest opposite side as Omni display

☐

1016511



						
Included in kit	1015013	1015014	00355.0477	00355.0471	1016718	1016716
1016504	1		1	1	1	1
1016511		1	1	1	1	1



EN

1. Assemble the Omni screen bracket.
2. Assemble the swing-away controller bracket for the Omni screen.
3. Assemble swing-away bracket for the joystick on the opposite side.
4. Fasten the black and red buddy buttons to the armrest.
5. Set up the cabling.

 In order to set up the system, see R-net Omni technical manual.

NO



1. Monter Omni-skjermbraketten
2. Monter den utsvingbare styreboksbraketten til Omni skjerm.
3. Monter den utsvingbare braketten til joysticken på motsatt side.
4. Fest de svarte og røde Buddy-knappene til armlenet.
5. Sett opp kablene.


 Se den tekniske manualen for R-net Omni for oppsett av systemet.

SE

1. Montera Omni-skärmens fäste.
2. Montera undansvängbara controller fästet för Omni-skärm.
3. Montera undansvängbara fästet för manöverboxen på motsatt sida.
4. Fäst de svarta och röda buddy-kontakterna till armstödet.
5. Installera kablaget.

 För att installera systemet, se R-net Omni teknisk manual.





Sunrise Medical HCM B.V.
Vossenbeemd 104
5705 CL Helmond
The Netherlands

T: +31 (0)492 593 888
customerservice@sunrisemedical.nl
www.SunriseMedical.eu