




Youngster 3



Fit for Adventures



Youngster 3

The lightest folding wheelchair
with growth adjustability!

Zippie® Youngster 3™:

The new lightweight Youngster 3 offers integrated seat width and seat depth growth. The new frame design gives your wheelchair a stylish and unique look.

On top of that the Youngster 3 offers new designed, innovative and lightweight options using latest technology.





Patented 20 mm seat width growth via cross brace for growing user (without new parts).



New lightweight fork and frame for maximum adjustability. Improved look and better manoeuvrability.



40 mm seat depth growth via seat tubes for growing user (without new parts).

- Lightweight package:

Lightweight fork, lightweight seat upholstery, lightweight anti-tips, Safari package

Explore your world! Have fun, enjoy playing with your friends and be active

Options

The new Youngster 3™ offers you and your parents a wide range of options for more comfort, freedom, fun and safety.



New selection of back upholsteries with a wide range of different versions and an improved adjustability.



Lightweight seat upholstery for easy handling and transport.

22 frame colours!

Select a colour for each component. Cushions are always black




Modern sideguard design in combination with the new Safari wheel lock for more safety.

Spokeguards

There are 32 fantastic spokeguard designs you can choose from. Here you can see 4 designs! See our spoke guard catalogue for the whole range.



Technical specifications



Frame colour : 22 colours available (see order form)

Weight of heaviest part: 4.4 kg

Upholstery colour: Black / Vented (black & silver) / Alcantara (grey & brown)

Product weight: from 8.5 kg

Max user weight: 85 kg

Front wheel diameter: 4" / 5" / 6"

Growth: Seat width: + 20 mm without parts
Seat depth: + 40 mm without parts

Rear wheel diameter: 20" / 22" / 24"

Frame Colours



It is best to select colours using the metal colour plates at your dealers and to use the order form to finalise your desired colour scheme

IMPORTANT NOTE

Due to a recent update, the colours shown may not reflect current availability. For the most up-to-date selection, please refer to the latest order form or contact your local Sunrise Medical representative.

Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.de

For further information on the full specification, options, and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.



04.2025



www.SunriseMedical.eu