

# SUPERIOR CLINICAL SEATING



# SUPERIOR CLINICAL SEATING



Achieving the optimal seated position, providing postural stability and effective pressure redistribution while respecting the individual anatomic dimensions of each patient – that is what JAY® is about.

Compromises in clinical use are not acceptable. Effectiveness and easy of use are the top priorities. For this reason, Jay® develops products that address these challenges. Every single aspect is thought through, down to the last detail. A brand that is able to adjust to every body type.

The result is a technological masterpiece: “Jay®” combining stability, effective postural and pressure management, easy handling and comfort to the highest degree.





Cushion	JAY® J3	JAY® J2
<b>POSITIONING</b>	<b>The ultimate modular seating solution with Air or Fluid technology</b>	<b>Stability and positioning with skin integrity</b>
<b>PURPOSE</b>	The ultimate modular seating solution. Extreme skin protection and positioning.	Designed for client with moderate to aggressive postural needs and / or high risk of skin breakdown.
<b>SKIN INTEGRITY</b>	■ ■ ■ ■	■ ■ ■ ■
<b>SHEAR</b>	■ ■ ■ ■	■ ■ ■ ■
<b>STABILITY</b>	■ ■ ■ ■	■ ■ ■ ■
<b>POSITIONING</b>	■ ■ ■ ■	■ ■ ■ ■
<b>SEATING TOLERANCE</b>	■ ■ ■ ■	■ ■ ■ ■
<b>SIZES</b>		
<b>WIDTH (CM)</b>	30 - 60 cm	36 - 61 cm
<b>DEPTH (CM)</b>	30 - 60 cm	41 - 51 cm
<b>HEIGHT FRONT (CM)</b>	9.5 cm / 11 cm (Standard / Deep)	7.6 cm
<b>HEIGHT REAR (CM)</b>	9.5 cm / 11 cm (Standard / Deep)	7.6 cm
<b>USER WEIGHT (KG)</b>	150 kg / 227 kg	150 kg
<b>PRODUCT WEIGHT (MIN)</b>	Fluid: 1.7 kg / Air: 1.1 kg	2,0 kg
<b>PRODUCT WEIGHT (MAX)</b>	Fluid: 4.5 kg / Air: 2.8 kg	5.4 kg
<b>CUSHION TYPE</b>		
<b>OPEN CELL FOAM</b>	-	-
<b>OPEN CELL FOAM+ FLUID</b>	-	-
<b>CLOSED CELL FOAM+ RETICULATE FOAM</b>	-	-
<b>CLOSED CELL FOAM+ FLUID</b>	-	Yes
<b>CLOSED CELL FOAM+ FLUID/ AIR</b>	Yes	-
<b>COVER TYPE</b>		
<b>INCONTINENCE</b>	Yes	Yes
<b>AIR EXCHANGE / MICROCLIMATIC</b>	Yes	Yes
<b>GUARANTEE / WARRANTY</b>		
<b>CUSHION</b>	2 years	2 years

The score given is based on a subjective assessment of the cushions performance in relation to the rest of the JAY® cushion range





## JAY® J2 DC

### Extreme skin integrity with positioning

Designed for clients with extreme skin risk and may require postural stability and / or positioning.



36 - 61 cm  
36 - 51 cm  
10.2 cm  
10.2 cm

150 kg  
2.3 kg  
5.4 kg

-  
-  
-  
Yes  
-

Yes  
Yes

2 years

## JAY® GS

### Paediatric stability, positioning and growth

The ideal cushion for the pediatric clients with high positioning needs who is a moderate to high risk for skin breakdown.



26 / 31 / 36 / 41 / 46 cm  
33 / 38 / 44 / 49 / 51 cm  
7 cm  
7 cm

113 kg  
0.8 kg  
3 kg

-  
-  
-  
Yes  
-


Yes  
No

2 years

## JAY® J2 Plus

### Pressure and shear reduction for the larger user

Pressure and shear reduction for the larger user



51 - 67 cm  
46 - 56 cm  
10.2 cm  
10.2 cm

295 kg  
3.2 kg  
5.0 kg

-  
-  
-  
Yes  
-

Yes  
Yes

2 years



Cushion	JAY® Care	JAY® Active
---------	-----------	-------------

<b>POSITIONING</b>	<b>Elderly care and posture</b>	<b>Low profile pressure and shear reduction</b>
<b>PURPOSE</b>	The ideal cushion for elderly patients with posterior pelvic tilt with moderate skin integrity and postural support needs. Features a longer seat well.	Low profile cushion, designed for the active users who requires high level pressure redistribution and stability.
<b>SKIN INTEGRITY</b>	■ ■ ■	■ ■ ■ ■
<b>SHEAR</b>	■ ■ ■ ■ ■	■ ■ ■ ■ ■
<b>STABILITY</b>	■ ■ ■	■ ■ ■ ■
<b>POSITIONING</b>	■ ■ ■	■ ■ ■
<b>SEATING TOLERANCE</b>	■ ■ ■	■ ■ ■ ■
<b>SIZES</b>		
<b>WIDTH (CM)</b>	41 / 46 / 51 cm	36 / 41 / 46 / 51 cm
<b>DEPTH (CM)</b>	45 / 50 cm	41 / 46 / 51 cm
<b>HEIGHT FRONT (CM)</b>	9.6 cm	5 cm
<b>HEIGHT REAR (CM)</b>	9.6 cm	7 cm
<b>USER WEIGHT (KG)</b>	113 kg	113 kg
<b>PRODUCT WEIGHT (MIN)</b>	3,3 kg	2.2 kg
<b>PRODUCT WEIGHT (MAX)</b>	4.1 kg	3.8 kg
<b>CUSHION TYPE</b>		
<b>OPEN CELL FOAM</b>	-	-
<b>OPEN CELL FOAM+ FLUID</b>	Yes	Yes
<b>CLOSED CELL FOAM+ RETICULATE FOAM</b>	-	-
<b>CLOSED CELL FOAM+ FLUID</b>	-	-
<b>CLOSED CELL FOAM+ FLUID/ AIR</b>	-	-
<b>COVER TYPE</b>		
<b>INCONTINENCE</b>	Yes	Yes
<b>AIR EXCHANGE / MICROCLIMATIC</b>	No	Yes
<b>GUARANTEE / WARRANTY</b>		
<b>CUSHION</b>	2 years	2 years

The score given is based on a subjective assessment of the cushions performance in relation to the rest of the JAY® cushion range



## JAY® Xtreme

**Lightweight cushion providing moderate pressure redistribution and shear protection for the active user**

Lightweight, moderate skin protection for active users: Stability and skin protection.



## JAY® Easy Fluid

**Skin integrity with stability**

Skin integrity with stability, and positioning for the active user.



## JAY® Lite

**Extremely Lightweight, breathable and skin protection cushion for the active user**

Extremely Lightweight, breathable and skin protection cushion for the active user.



36 / 38 / 41 / 44 / 46 / 51 cm  
36 / 38 / 41 / 44 / 46 / 51 cm  
5.5 cm  
5.5 cm

25 - 60 cm  
25 - 60 cm  
6.3 cm  
6.3 cm

35 - 60 cm  
35 - 60 cm  
8 cm  
7.5 cm

150 kg  
2.2 kg  
3.8 kg

150 kg  
1.8 kg  
6 kg

125 kg  
0.55 kg  
1.62 kg

-  
Yes  
-  
-  
-

-  
Yes  
-  
-  
-

Yes  
No  
No

Yes  
Yes

Yes  
Yes

No  
Yes

2 years

2 years

2 years



Cushion	JAY® Zip Cushion	JAY® Easy Visco
<b>POSITIONING</b>	<b>Designed for pediatrics</b>	<b>High level comfort and stability</b>
<b>PURPOSE</b>	Designed to meet the pediatric client's unique needs.	Lightweight cushions with high level comfort and stability.
<b>SKIN INTEGRITY</b>	■ ■ ■	■ ■ ■
<b>SHEAR</b>	■ ■	■ ■
<b>STABILITY</b>	■ ■ ■ ■	■ ■ ■ ■
<b>POSITIONING</b>	■ ■ ■	■ ■ ■
<b>SEATING TOLERANCE</b>	■ ■ ■ ■	■ ■ ■ ■
<b>SIZES</b>		
<b>WIDTH (CM)</b>	20 - 40 cm	25 - 60 cm
<b>DEPTH (CM)</b>	20 - 46 cm	25 - 60 cm
<b>HEIGHT FRONT (CM)</b>	5 cm	6.3 cm
<b>HEIGHT REAR (CM)</b>	4,5 cm	6.3 cm
<b>USER WEIGHT (KG)</b>	75 kg	150 kg
<b>PRODUCT WEIGHT (MIN)</b>	0.3 kg	0.50 kg
<b>PRODUCT WEIGHT (MAX)</b>	1.25 kg	2.70 kg
<b>CUSHION TYPE</b>		
<b>OPEN CELL FOAM</b>	Yes	Yes
<b>OPEN CELL FOAM+ FLUID</b>	-	-
<b>CLOSED CELL FOAM+ RETICULATE FOAM</b>	-	-
<b>CLOSED CELL FOAM+ FLUID</b>	-	-
<b>CLOSED CELL FOAM+ FLUID/ AIR</b>	-	-
<b>COVER TYPE</b>		
<b>INCONTINENCE</b>	Yes (inner cover)	Yes
<b>AIR EXCHANGE / MICROCLIMATIC</b>	Yes (outer Cover)	Yes
<b>GUARANTEE / WARRANTY</b>		
<b>CUSHION</b>	2 years	2 years

The score given is based on a subjective assessment of the cushions performance in relation to the rest of the JAY® cushion range





## JAY® Soft Combi P

### Long life stability

Designed for the client who has symmetrical posture and low risk of skin breakdown needing low to moderate postural support.



25 - 60 cm  
25 - 60 cm  
5.9 cm  
4.4 cm

150 kg  
0.34 kg  
1.98 kg

Yes  
-  
-  
-  
-

Yes  
No

5 years

## JAY® Basic

### Long life comfort

Designed for the symmetrical client with low skin risk who requires comfort and mild postural support.



25 - 60 cm  
25 - 60 cm  
4.8 cm  
4.7 cm

113 kg  
0.3 kg  
1.73 kg

Yes  
-  
-  
-  
-

Yes  
No

2 years



**Backrest**

**JAY® J3 back**

**POSITIONING**

**Simply - Effective - Seating**

**PURPOSE**

Anthropometrically designed for support and positioning. Wide range of shapes and sizes to match user's clinical needs.

**SIZES**

WIDTH (CM) 31 / 36 / 41 / 46 / 51 cm  
 DEPTH (CM) SC 5 cm / MC 8 cm / DC & PDC 15 cm  
 HEIGHT (CM) LT (17 / 20 / 24 cm) / MT (30 / 34 / 38 cm)  
 UT (42 / 46 / 50 cm) / SH (53 / 57 / 61 cm)

**USER WEIGHT (KG)**

136 kg

**COVER TYPE**

**INCONTINENCE**

No

**AIR EXCHANGE / MICROCLIMATIC**

Yes

**GUARANTEE / WARRANTY**

2 years





### JAY® Easy back

### JAY® Zip back

#### Easy to prescribe, to order and to use

#### Designed for pediatrics

The Jay Easy backrest is designed to offer simple and comfortable support.

Lightweight backrest designed specifically for kids anatomical measurements.

36 / 41 / 46 / 51 cm  
5 cm & 15 cm  
38 / 46 / 53 cm

20 / 25 / 30 / 35 cm  
2,5 / 3,8 / 5 / 6,3 cm  
15 / 20 / 25 / 30 cm

136 kg

75 kg

No  
Yes

Yes (inner cover)  
Yes (outer cover)

2 years

2 years





**Sunrise Medical GmbH**

Kahlbachring 2-4  
69254 Malsch/HD  
Tel.: +49 (0) 72 53 / 9 80-0  
Fax: +49 (0) 72 53 / 9 80-111  
[www.SunriseMedical.com](http://www.SunriseMedical.com)