



# ATTITUDE JUNIOR

WHERE ADVENTURES BEGIN



# THE RIGHT ATTITUDE FOR DISCOVERERS

THE FIRST BIKE OF CHOICE TO OPEN  
UP A NEW WORLD OF EXPLORATION  
AND FREEDOM.

---

We want kids to love riding just  
as much as we do. That's  
why our Junior bikes ride  
better, last longer,  
and carry kids safely  
to more adventures.

// Empulse Attitude Junior  
with Zippie Youngster 3





**ATTITUDE**  
*Junior*

# CYCLING IS IDEAL FOR KIDS

- n It contributes to natural growth development
- n Raises self-confidence through participation and independence
- n Enhances coordination and body awareness
- n Strengthens the cardiovascular system






# READY TO GO IN SECONDS

The Attitude Junior's clever docking mechanism attaches to your wheelchair in a matter of seconds...and it's just as quick to remove. Ideal for explorers who don't have time to lose when they want to change between cycling or wheelchair use.

After undocking only two small pins (170 g) remain on the wheelchair.

The docking mechanism on the Attitude Junior bike is so simple, it just locks securely to the front of the wheelchair – without any need for parts or cluttered wiring around the knees. So there's nothing coming between them and the comfort that they have come to expect from their wheelchair.



- 1 Lock in the docking pins.
- 2 Swing the upper docking connector downwards.
- 3 Lock the connection by pulling the red lever.





# IDEAL FOR KIDS ON THE GROW!

Go together – grow together. Easy depth- and angle adjustment of the Attitude Junior bike ensures an ergonomic well-set crank position (depth up to 100 mm, angle range 28°). This is important for an energy efficient driving experience and a must have for growing kids.





# WHERE ADVENTURES BEGIN

---



The Attitude Junior is packed with features that will have your child riding with control, confidence and comfort. 16" Front wheels, the 8 speed hub-gearing and horizontal handles provide an exceptional driving performance which is perfect for all terrains. Plus, the powerful backpedal brake allows a smooth speed control. Light-weight and easily manoeuvrable, it is approved for fun anywhere that it goes.





# ATTITUDE

## Junior



### // LET IT ROLL

The 8 speed hub-gearing is perfect for all terrains. The powerful backpedal brake allows for smoother speed control.



### // HORIZONTAL HANDLES

The horizontal push handle position is ideal for children and supports the paediatric arm ergonomics.



### // ERGONOMIC HANDLE

The ergonomic handle helps to manoeuvre the bike easily with just one hand while undocking or storing it away.



### // PARKING BRAKE

Easy operating additional parking brake with lock-in function as standard on all Attitudes.



# TECHNICAL DATA

## ATTITUDE JUNIOR // Add-on wheelchair bike for kids

Max. user weight	100 kg
Front frame width (adj. + 2 cm)	260 - 460 mm
Weight of parts at chair	170 g
Total bike weight	14,5 kg
Wheel size	16"
Gearing shifters	Grip shift
Gearing speeds	8 speed hub gearing

V-crank width	240 mm
Crank length	110 mm / 130 mm / 150 mm
Crank depth adjustment	100 mm
Crank angle adjustment	28°
Brake	Backpedal brake
Parking brake	yes
Frame colour	22

### ⚠ IMPORTANT NOTE

Due to a recent update, the colours shown may not reflect current availability. For the most up-to-date selection, please refer to the latest order form or contact your local Sunrise Medical representative.



### Sunrise Medical GmbH

Kahlbachring 2-4  
69254 Malsch/HD  
Tel.: +49 (0) 72 53 / 9 80-0  
Fax: +49 (0) 72 53 / 9 80-222  
E-Mail: [international@sunrisemedical.de](mailto:international@sunrisemedical.de)  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.eu](http://www.SunriseMedical.eu).



04/2025

