

Combination Agreement

The below mentioned products are in principle appropriate to be combined with the below mentioned **Otto Bock** products provided that the combination is made expertly and according to the restrictions of both manufacturers. The restrictions can be asked at the manufacturers.

SUNRISE MEDICAL product	Otto Bock product
Quickie Wheeldrive	Motus CS Motus CV Avantgarde CV Avantgarde CS Avantgarde 4 Avantgarde VR Avantgarde T Start M2 Effect Start M5 Comfort Start4 M2 Ventus Motus

Please note that mounting and limitations has to be effected according to the Sunrise Medical and Otto Bock mounting instructions by authorized and skilled staff only.

The following limitations must be considered:

- The above listed combinations can only be used up to a speed of max 6kph
- The combination for Motus, Ventus and Avantgarde (including its variations) is released for castor sizes of 5,5", 6" and 7" only.
- The combination with the Start series is released for castor sizes 7" and 8" only

The maximum values and restrictions (e.g. user weight) mentioned in the user manual of both manufacturers must be observed.

Machining, bending, welding or bracing on any safety relevant components is not allowed.

The combination agreement has been made in two copies, one for each party. Authorized representatives of the parties have executed this agreement.

SUNRISE MEDICAL Europe
 Kahlbachring 2-4
 69254 Malsch
 Germany

Otto Bock Mobility Solutions GmbH
 Lindenstraße 13
 07426 Königsee-Rottenbach
 Germany

Malsch, 30th October 2017



Gerd Weis
 (Senior test engineer)

3. November 2017

