

## Combination Agreement

The below mentioned products are in principle appropriate to be combined with the below mentioned **Meyra** products provided that the combination is made expertly and according to the restrictions of both manufacturers. The restrictions can be asked at the manufacturers.

| SUNRISE MEDICAL product | Meyra product  |
|-------------------------|--|
| Quickie Wheeldrive      | Hurricane<br>ZX1, ZX3<br>X-Star<br>Avanti, Avanti Pro<br>Eurochair <sup>2</sup> 2.750; 2.850<br>Eurochair 1.750, 1.850<br>Service 3.600<br>Smart <sup>F</sup> , Smart <sup>S</sup><br>X1 |

Please note that mounting and limitations has to be effected according to the Sunrise Medical and MEYRA mounting instructions by authorized and skilled staff only.

The maximum values and restrictions (e.g. user weight or total weight) mentioned in the user manual of both manufacturers must be observed.

Machining, bending, welding or bracing on any safety relevant components is not allowed. The minimum steering wheel size is 142mm.


The damages wich stand in the direct connection with the wheeldrive are taken over from sunrise medical.

The combination agreement has been made in two copies, one for each party. Authorized representatives of the parties have executed this agreement.

**SUNRISE MEDICAL Europe**  
Kahlbachring 2-4  
69254 Malsch  
Germany

**MEYRA GmbH**  
Meyra Ring 2  
32689 Kalletal-Kalldorf  
Germany

Malsch, 8<sup>th</sup> May 2017



Gerd Weis  
(Senior test engineer)



Piotr Baczynski  
(General Manager)