





Nakkeputer




Beskrivelse


Artikkelnummer

A Nakkepute Ergo		Skifergrått stoff	☐	12775512
A Nakkepute Ergo		Sort hud	☐	12785512
B Nakkepute liten		Skifergrått stoff	☐	12778008
B Nakkepute liten		Sort hud	☐	12788008
C Nakkepute stor	Vatt	Skifergrått stoff	☐	12778006
C Nakkepute stor	Vatt	Sort hud	☐	12788006
D Nakkepute stor	Sacco kuler	Skifergrått stoff	☐	12778007
D Nakkepute stor	Sacco kuler	Sort hud	☐	12788007


NOVA




1x**A**




1x**B**




1x**C**



1x**D**




1



2

Tilpasning ③

Nakkepute Liten ③ kan snus opp ned ved å tre stoffet med lodd gjennom smøystolen på baksiden av puten, som vist i figur 3.





Sunrise Medical AS
Delitoppen 3
1540 Vestby
Norway
T: +47 66963800
E: post@sunrisemedical.no
www.SunriseMedical.no