

Q500® M/F/R/H

with SEDEO® PRO



Rev:04/2024
Valid:11/2024


ORDER FORM

Errors and omissions excepted.

| Billing Address Or Customer-ID | | | Shipping Address | Identical To Billing Address | |
|--|-------------|-------------|-------------------------------------|------------------------------|--|
| C.-ID: | | | C.-ID: | | |
| Name: | | | Name: | | |
| Street: | | | Street: | | |
| Town/PC: | | | Town/PC: | | |
| E-Mail: | | | E-Mail: | | |
| Phone: | | | Phone: | | |
| Order Information | | | User Weight | | |
| Order | Calculation | Built-4-Me® | User Weight Capacity: 160 kg | | |
| Purchase Order No.: | | | User Weight: kg | | |
| Art. No.: Q500 M: XXX, Q500 F: XXX, Q500 R: XXX, Q500 H: XXX | | | Reimbursement Code: xxx | | |
| Order Date: | Day | Month | Year | Base Price € XX.XXX,- | |

Client Measurements

| | |
|----------------------------|----|
| a) Top Of Head To Seat | mm |
| b) Top Of Shoulder To Seat | mm |
| c) Underarm To Seat | mm |
| d) Elbow To Seat | mm |
| e) Top Of Pelvis To Seat | mm |
| f) Shoulder Width | mm |
| g) Trunk Width | mm |
| i) Hip Width | mm |
| j) Actual Sitting Depth | mm |
| m) Foot Length | mm |
| n) Lower Leg Length | mm |
| o) Actual Seat Width | mm |



Seating Specifications

| | | | |
|---|--------------------------|-----------------------------------|--|
| All Bases Prepared For DAHL Docking Station | | Seat Type | SEDEO PRO |
| Typical WC Weight | Without Batteries 140 kg | Seat Size Range | Width 380 - 560 mm, Depth 400 - 560 mm |
| | With Batteries 165 kg | Seat To Floor Height Range | 415 mm - 490 mm |

MWD Base Specifications

| | | | |
|-----------------------------|--|------------------------|----------------------------|
| User Weight Capacity | 160 kg | Drive Wheels | 13" Drive Wheels (Black) |
| Base Size | Width Min. 620 mm | Tires / Castors | 7" Front / 7" Rear (Black) |
| Drive Control | VR2, Option R-Net | Charger | 8 A Dual Mode Charger |
| Driving Range | Up To 40 km | Battery Size | 50, 60 And 80 Ah Batteries |
| Drive Motors | 2-Pole 6 kph, 4-Pole 6/8/10 And 12.5 kph | Colour Accents | 7 Colour Accents |
| Ground Clearance | 70 mm | Frame Colour | Black |

FWD Base Specifications

| | | | |
|-----------------------------|--|------------------------|--------------------------------|
| User Weight Capacity | 160 kg | Drive Wheels | 13" / 14" Drive Wheels (Black) |
| Base Size | Width Min. 610 mm | Tires / Castors | 9" Rear (Black), Option 10" |
| Drive Control | VR2, Option R-net | Charger | 8 A Dual Mode Charger |
| Driving Range | Up To 40 km | Battery Size | 50, 60 And 80 Ah Batteries |
| Drive Motors | 2-Pole 6 kph, 4-Pole 6/8/10 And 12.5 kph | Colour Accents | 7 Colour Accents |
| Ground Clearance | 65 mm | Frame Colour | Black |

RWD Base Specifications

| | | | |
|-----------------------------|--|------------------------|--------------------------------|
| User Weight Capacity | 160 kg | Drive Wheels | 13" / 14" Drive Wheels (Black) |
| Base Size | Width Min. 610 mm | Tires / Castors | 9" Front (Black), Option 10" |
| Drive Control | VR2, Option R-net | Charger | 8 A Dual Mode Charger |
| Driving Range | Up To 40 km | Battery Size | 50, 60 And 80 Ah Batteries |
| Drive Motors | 2-Pole 6 kph, 4-Pole 6/8/10 And 12.5 kph | Colour Accents | 7 Colour Accents |
| Ground Clearance | 65 mm | Frame Colour | Black |

HYBRID Base Specifications

| | | | |
|-----------------------------|--|------------------------|--------------------------------|
| User Weight Capacity | 160 kg | Drive Wheels | 13" / 14" Drive Wheels (Black) |
| Base Size | Width Min. 620 mm | Tires / Castors | 9" Front / 4" Rear (Black) |
| Drive Control | VR2, Option R-net | Charger | 8 A Dual Mode Charger |
| Driving Range | Up To 40 km | Battery Size | 50, 60 And 80 Ah Batteries |
| Drive Motors | 4-Pole 6 kph, Options 8, 10 And 12.5 kph | Colour Accents | 7 Colour Accents |
| Ground Clearance | 65 mm | Frame Colour | Black |

Base

Q500 Drive Bases

| | | € |
|-----------|-----------------------------------|---|
| QOP000023 | Q500 M (MWD Base Chair) | – |
| QOP000024 | Q500 F (FWD Base Chair) | – |
| QOP000029 | Q500 R (RWD Base Chair) | – |
| QOP000030 | Q500 H (HYBRID Base Chair) | – |

Suspension Setting (MWD Only)

| | | € |
|-----------|--|-------|
| QOP150001 | Standard Suspension | – |
| QOP150002 | Firm Suspension Recommended with QOP130009 Seat elevator/tilt module, Recommended with Powered tilt (QOP130004/QOP130012/ QOP130007 for clients over 75 kg user weight | XXX,- |

Motor Packages

| | | € |
|--------------------------------|---|-------|
| MWD Base | | |
| QOP180011 | 2-Pole Motor, 6 kph / 4 mph, Standard | XXX,- |
| MWD/FWD/RWD/HYBRID Base | | |
| QOP180012 | 2-Pole Motor, 6 kph / 4 mph, Performance Only With Drive Wheel 14" | XXX,- |
| QOP180009 | 4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed) | XXX,- |
| QOP180010 | 4-Pole Motor, 8 kph / 5 mph | XXX,- |
| QOP180007 | 4-Pole Motor, 10 kph / 6.5 mph If FWD Selected: Highly Recommended With Traction Control System | XXX,- |
| QOP180008 | 4-Pole Motor, 12.5 kph / 8 mph If FWD Selected: Only In Combination With Traction Control System | XXX,- |

Tracking

| | | € |
|-----------|---|-------|
| QOP110009 | Traction Control System For R-Net (GYRO Technology) | XXX,- |
| QOP110014 | Omit Traction Control System For R-Net | XXX,- |

Drive Wheels 

| | | € |
|--------------------------------|--|-------|
| MWD/FWD/RWD/HYBRID Base | | |
| QOP070065 | Drive Wheel, 13", Pneumatic, Knobbly, Black | XXX,- |
| QOP070069 | Drive Wheel, 13", Solid, Knobbly, Black | XXX,- |
| QOP070034 | Drive wheel, 14", Pneumatic, Knobbly, Black | XXX,- |
| QOP070033 | Drive wheel, 14", Solid, Knobbly, Black | XXX,- |
| FWD/RWD Base | | |
| QOP070032 | Drive Wheel, 14", Pneumatic, V-Type, Black | XXX,- |
| QOP070031 | Drive Wheel, 14", Solid, V-Type, Black Not Available With HYBRID | XXX,- |
| QOP070091 | Drive Wheel, 14", Pneumatic, Trelleborg, Grey Only on FWD; Only With QOP080065 Castor Wheel, 9", Solid, Grey; Only with QOP180007 4-Pole Motor 10 kph; Only for Norway, Only For NL Supplied Products | XXX,- |

Castor Wheels 

| | | € |
|------------------------|---|-------|
| MWD Base | | |
| QOP080213 | Castor Wheel, 7", Solid, Black | XXX,- |
| FWD/HYBRID Base | | |
| QOP080210 | Castor Wheel, 9", Solid, Black | XXX,- |
| QOP080211 | Castor Wheel, 9", Pneumatic, Black | XXX,- |
| QOP080065 | Castor Wheel, 9", Solid, Grey Only on FWD; Only with QOP070091 Drive Wheel 14" Pneumatic Trelleborg Grey; Only for Norway, Only For NL Supplied Products | |
| RWD Base | | |
| QOP080210 | Castor Wheel, 9", Solid, Black | XXX,- |
| QOP080211 | Castor Wheel, 9", Pneumatic, Black | XXX,- |
| QOP080058 | Castor Wheel, 10", Solid, Black Only With Drive Wheel 14" | XXX,- |
| QOP080212 | Castor Wheel, 10", Pneumatic, Black Only With Drive Wheel 14" | XXX,- |

Mudguards (FWD/RWD/HYBRID Only)  €

| | | |
|-----------|---|-------|
| QOP010007 | Mudguard, Drive Wheel Standard Included on HYBRID Base | XXX,- |
| QOP010008 | Mudguard, Castor Wheel, 9" | XXX,- |

Kerb Climber (RWD Only)  €

| | | |
|-----------|--------------|-------|
| QOP090038 | Kerb Climber | XXX,- |
|-----------|--------------|-------|

Colour Accents €

| | | |
|-----------|-----------------------------|-------|
| QOP140005 | Apple Green (RAL 110 70 77) | XXX,- |
| QOP140004 | Black (RAL 9005) | XXX,- |
| QOP140003 | Blue (RAL 260 40 45) | XXX,- |
| QOP140006 | Orange (RAL 060 60 80) | XXX,- |
| QOP140002 | Red (RAL 040 40 67) | XXX,- |
| QOP140007 | Silver (RAL 9006) | XXX,- |
| QOP140001 | White (RAL 9010) | XXX,- |

LED Lights & Indicators  €

| | | |
|-----------|--|-------|
| QOP110201 | LED Light And Indicator QOP110030 Lighting Module VR2 Or QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | XXX,- |
| QOP110202 | LED Fender Light QOP110030 Lighting Module VR2 or QOP110032 Actuator Control & Lights Module required – Not Road Approved And Without Indicators; On Q500 M Only With QOP040123 Armrest, Reclining, Flip Back or QOP040070 Armrest, Reclining, Multi Adjustable, Flip Back or QOP040073 Armrest, Reclining, Multi Adjustable, Flip Back, With Lock Function; Not With QMP130012 Posterior Tilt, Powered, -5° - 25°; On Q500 F and R Only With QOP010007 Mudguard, Drive Wheel; Not With SH415 | XXX,- |
| QOP110209 | Omit Light And Indicator | XXX,- |

Seat Heights Measured From Front Of Seat Pan To Floor. Obtain Quote To Verify Seat To Floor Height €

| | | |
|-------|---|---|
| SH415 | 415 mm Only With 60 Ah Batteries, Not With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered | - |
| SH425 | 425 mm Only With 60 Ah Batteries, Only With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered | - |
| SH44 | 440 mm Not With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered | - |
| SH45 | 450 mm Only With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered | - |
| SH465 | 465 mm | - |
| SH49 | 490 mm | - |

Batteries & Chargers  €

| Batteries | | |
|------------------|--|-------|
| QOP160003 | Battery, 50 Ah, Gel, Maintenance-Free | XXX,- |
| QOP160001 | Battery, 60 Ah, Gel, Maintenance-Free | XXX,- |
| QOP160005 | Battery, 80 Ah, Gel, Maintenance-Free Only With Battery Charger 220 V / 10 A To Reduce Reloading Time | XXX,- |
| QOP160033 | Battery, 60 Ah (C20), Gel, Maintenance-Free, Set, Type Exide Only for NL Helmond supplied products | XXX,- |
| QOP160034 | Battery, 70 Ah (C20), Gel, Maintenance-Free, Set, Type Exide Only for NL Helmond supplied products | XXX,- |
| QOP160009 | Omit Batteries | XXX,- |
| Chargers | | |
| QOP160101 | Charger, Automatic, 220 V / 8 A | XXX,- |
| QOP160100 | Charger, Automatic, 220 V / 8 A, UK | XXX,- |
| QOP160102 | Charger, Classic, 220 V / 8 A, IP54 | XXX,- |
| QOP160103 | Charger, Automatic, 220 V / 10 A, UK | XXX,- |
| QOP160104 | Charger, Automatic, 220 V / 10 A | XXX,- |
| QOP160029 | Charger, Automatic, Vitron Impulse25 smart (IXP4) 220 V / 8 A Only for Norway | XXX,- |
| QOP160109 | Omit Charger | XXX,- |
| QOP160035 | Charger, Automatic, IVRA, 220V / 8A Only for NL Helmond supplied products | XXX,- |

Powered Seating

Power Seating Options

€

| | | |
|-----------|--|---|
| QOP130011 | Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°) | – |
| QOP130004 | Posterior Tilt, Powered, 0° - 30° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | – |
| QOP130012 | Posterior Tilt, Powered, -5° - 25° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required, Not With Seat Height 415 Or 440 mm | – |
| QOP130007 | Posterior Tilt, Powered, 0° - 50° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | – |
| QOP130005 | Posterior Tilt, Powered, -5° - 45° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required, Not With Seat Height 415 Or 440 mm | – |
| QOP130009 | Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | – |
| QOP130006 | Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered Not With QOP110036 Remote Joystick, VR2; In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | – |

C-Me Function C-Me Function Require QMP130009 / QMP130006. Not available with QMP030084 /QMP030085 Backrest, Manual Recline

€

| | | |
|-----------|---|-------|
| QOP110034 | C-Me Function, Enabled 4.5 kph Maximum Speed When Lifted (Tilt Angle <9°), Not Available With 50° Lift/Tilt In Combination With 90°-150° Power Recline | – |
| QOP110035 | C-Me Function, Disabled 1.5 kph Maximum Speed When Lifted | XXX,- |

Backrest Frame Type

€

| | | |
|-----------|---|-------|
| QOP030083 | Backrest, Tool Angle Adjustable, 90° - 120° | – |
| QOP030084 | Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right Not On Same Joystick Side | XXX,- |
| QOP030085 | Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left Not On Same Joystick Side | XXX,- |
| QOP030086 | Backrest, Powered Recline, 90° - 130° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | XXX,- |
| QOP030087 | Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required | XXX,- |
| QOP030229 | Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear with Relax function Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required; Only With QMP130006 Seat Elevator (300mm) Tilt Module (Tilt 0° - 50°), Powered | XXX,- |

Backrest Frame Option

€

| | | |
|-----------|------------------------------------|-------|
| QOP030098 | Backrest, Flip Down, Quick Release | XXX,- |
|-----------|------------------------------------|-------|

Back Cushion Type

€

SEDEO PRO (See Section A)

| | | |
|-----------|-------------------------|---|
| QOP030104 | SEDEO PRO, Back Cushion | – |
|-----------|-------------------------|---|

JAY J3 (See Section B)

| | | |
|-----------|----------------------|-------|
| QOP030054 | JAY J3, Back Cushion | XXX,- |
|-----------|----------------------|-------|

Interfaces

| | | |
|-----------|--|---|
| QOP030106 | SEDEO PRO, Interface Only | – |
| QOP030107 | SEDEO PRO, Interface Only For Back Cushion, Flexible | – |
| QOP030055 | JAY J3, Interface Only | – |
| QOP030209 | Omit Back Cushion | – |

A – SEDEO PRO Back Size

€

Back Cushion Width

| | | |
|------|---|---|
| BW38 | 380 mm, For Effective Back Width 380 mm | – |
| BW42 | 420 mm, For Effective Back Width 380 - 420 mm | – |
| BW46 | 460 mm, For Effective Back Width 420 - 460 mm | – |
| BW50 | 500 mm, For Effective Back Width 460 - 500 mm | – |
| BW54 | 540 mm, For Effective Back Width 500 - 540 mm | – |
| BW56 | 560 mm, For Effective Back Width 520 - 560 mm | – |

Back Cushion Height

| | | |
|------|--|---|
| BH44 | 440 mm, For Effective Back Height 440 - 500 mm Only With SEP030116 SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps & QOP030117 SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adj Straps | – |
| BH50 | 500 mm, For Effective Back Height 500 - 560 mm | – |
| BH56 | 560 mm, For Effective Back Height 560 - 620 mm | – |

A – SEDEO PRO Back Cushion

€

| | | |
|-----------|--|-------|
| QOP030109 | SEDEO PRO, Back Cushion, Profiled | XXX,- |
| QOP030113 | SEDEO PRO, Back Cushion, Profiled, Adjustable Straps | XXX,- |
| QOP030114 | SEDEO PRO, Back Cushion, Deep Contour, Passive | XXX,- |
| QOP030115 | SEDEO PRO, Back Cushion, Flexible, Deep Contour, Passive, Adjustable Straps Only Available For NL | XXX,- |
| QOP030116 | SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps | XXX,- |
| QOP030117 | SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps | XXX,- |
| | SEDEO PRO, Back Cushion, Adaptive (Choose B4Me Requirements Separately) | XXX,- |

A – SEDEO PRO Back Cushion Cover

€

| | | |
|-----------|--|-------|
| QOP030059 | Back Cushion Cover, Incontinence Protection | – |
| QOP030108 | Back Cushion Cover, Comfort | XXX,- |
| QOP030230 | Back Cushion Cover, Comfort, Double zipper only for QOP030116 Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps and QOP030117 Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps. Only For Norway, Only For NL Supplied Products | XXX,- |

A – SEDEO Lateral Thoracic Supports Only With QOP030109 SEDEO PRO, Back Cushion, Profiled & QOP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps

€

| | | |
|-----------|--|-------|
| QOP030607 | SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right | XXX,- |
| QOP030608 | SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left | XXX,- |

B – JAY J3 Back Size

€

| | | |
|-----------|--------------------------|-------|
| JAY030111 | JAY J3 Back Size, 410 mm | XXX,- |
| JAY030112 | JAY J3 Back Size, 460 mm | XXX,- |
| JAY030113 | JAY J3 Back Size, 510 mm | XXX,- |

B – JAY J3 Back, Mid Deep Contour

€

| | | |
|-----------|---|-------|
| JAY030042 | JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm | XXX,- |
|-----------|---|-------|

B – JAY J3 Back, Deep Contour

€

| | | |
|-----------|---|-------|
| JAY030021 | JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm | XXX,- |
|-----------|---|-------|

Headrest Pad

€

| | | |
|-----------|---|-------|
| QOP030430 | SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm Only Use With QOP030433 or QMP030434 SEDEO Headrest Mount | XXX,- |
| QOP030431 | SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm Only Use With QOP030433 or QMP030434 SEDEO Headrest Mount | XXX,- |
| QOP030432 | SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm Only Use With QOP030433 or QMP030434 SEDEO Headrest Mount | XXX,- |
| WHT030402 | Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm Only Use With WHT030408 / WHT030409 / WHT030411 Whitmyer Headrest Mount | XXX,- |
| WHT030404 | Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm Only Use With WHT030408 / WHT030409 / WHT030411 Whitmyer Headrest Mount | XXX,- |

Headrest Mount Choose A Whitmyer Headrest Mount And Headrest Receiver In Case You Would Like To Add A Head Control In The Special Control Section

| | | |
|-----------|--|-------|
| QOP030433 | SEDEO, Headrest Mount | XXX,- |
| QOP030434 | SEDEO, Headrest Mount, With Quick Release | XXX,- |
| WHT030408 | Whitmyer, Headrest Mount, LINX ² | XXX,- |
| WHT030409 | Whitmyer, Headrest Mount, AXYS | XXX,- |
| WHT030411 | Whitmyer, Headrest Mount, COBRA XTRA Flip Back | XXX,- |

Headrest Receiver

| | | |
|-----------|---|-------|
| QOP030119 | Headrest Receiver, Back Cushion Mount Only With QOP030109 SEDEO PRO, Back Cushion, Profiled & QOP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps | XXX,- |
| QOP030120 | Headrest Receiver, Backrest Frame Mount | XXX,- |

Push Handles

| | | |
|-----------|---------------------------------|-------|
| QOP030204 | Push Handles, Height-Adjustable | XXX,- |
|-----------|---------------------------------|-------|

Seat Rail

| | | |
|-----------|--|-------|
| QOP020085 | Seat Rail, Short, For Seat Depth 400 - 560 mm | - |
| QOP020069 | Seat Rail, Long, For Seat Depth 460 - 560 mm Not With Seat Width 380 - 420 mm | XXX,- |

Seat Cushion Type

SEDEO (See Section C)

| | | |
|-----------|---------------------|---|
| QOP020027 | SEDEO, Seat Cushion | - |
|-----------|---------------------|---|

JAY (See Section D)

| | | |
|-----------|-------------------|-------|
| QOP020515 | JAY, Seat Cushion | XXX,- |
|-----------|-------------------|-------|

Omit (Select Seat Width And Depth in Section C)

| | | |
|-----------|-------------------|-------|
| QOP020110 | Omit Seat Cushion | XXX,- |
|-----------|-------------------|-------|

C – SEDEO Seat Size

Seat Cushion Width

| | | |
|------|---|---|
| SW38 | 380 mm, For Effective Seat Width 380 mm | - |
| SW42 | 420 mm, For Effective Seat Width 380 - 420 mm | - |
| SW46 | 460 mm, For Effective Seat Width 420 - 460 mm | - |
| SW50 | 500 mm, For Effective Seat Width 460 - 500 mm | - |
| SW54 | 540 mm, For Effective Seat Width 500 - 540 mm | - |
| SW56 | 560 mm, For Effective Seat Width 520 - 560 mm | - |

Seat Cushion Depth

| | | |
|------|--|---|
| SD40 | 400 mm, Comes With Cushion Size 460 mm Not with QOP090128 or QOP090129 Positioning belt Short or Long | - |
| SD42 | 420 mm, Comes With Cushion Size 460 mm | - |
| SD44 | 440 mm, Comes With Cushion Size 460 mm | - |
| SD46 | 460 mm, Comes With Cushion Size 500 mm | - |
| SD48 | 480 mm, Comes With Cushion Size 500 mm | - |
| SD50 | 500 mm, Comes With Cushion Size 500 mm | - |
| SD52 | 520 mm, Comes With Cushion Size 560 mm | - |
| SD54 | 540 mm, Comes With Cushion Size 560 mm | - |
| SD56 | 560 mm, Comes With Cushion Size 560 mm | - |

C – SEDEO Seat Cushions  €

| | | |
|-----------|---|-----|
| QOP020073 | SEDEO, Seat Cushion, Profiled, Height 70 mm | XXX |
| QOP020074 | SEDEO, Seat Cushion, Memory, Height 70 mm Only With QOP020518 Seat Cushion Cover, Incontinence Protection, Only For NL | XXX |
| QOP020075 | SEDEO, Seat Cushion, Support, Height 90 mm | XXX |
| QOP020076 | SEDEO, Seat Cushion, Support Memory, Height 90 mm | XXX |
| QOP020077 | SEDEO, Seat Cushion, Form, Height 110 mm Only With QOP020518 Seat Cushion Cover, Incontinence Protection, Only For NL | XXX |
| QOP020078 | SEDEO, Seat Cushion, Stabilizer, Height 90 mm | XXX |
| | SEDEO, Seat Cushion, Adaptive (Choose B4Me Requirements Separately) | XXX |

C – SEDEO Seat Cushion Cover €

| | | |
|-----------|---|-----|
| QOP020518 | Seat Cushion Cover, Incontinence Protection | XXX |
| QOP020080 | Seat Cushion Cover, Comfort | XXX |

C – SEDEO Seat Boards €

| | | |
|-----------|---|-----|
| QOP020082 | SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear | XXX |
| QOP020083 | SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear | XXX |
| QOP020084 | SEDEO, Seat Pan, 75 mm, Depth-Adjustable, Cushions Guide At Rear Only For NL | XXX |

D – JAY Seat Size €**Seat Cushion Width**

| | | |
|------|--|---|
| SW40 | 400 mm | – |
| SW42 | 420 mm | – |
| SW44 | 440 mm | – |
| SW46 | 460 mm | – |
| SW48 | 480 mm | – |
| SW50 | 500 mm, Comes With Cushion Size 500 mm | – |
| SW52 | 520 mm, Comes With Cushion Size 500 mm | – |
| SW54 | 540 mm, Comes With Cushion Size 560 mm | – |
| SW56 | 560 mm, Comes With Cushion Size 560 mm | – |


Seat Cushion Depth

| | | |
|------|---------------------------------------|---|
| SD40 | 400 mm | – |
| SD42 | 420 mm | – |
| SD44 | 440 mm | – |
| SD46 | 460 mm | – |
| SD48 | 480 mm | – |
| SD50 | 500 mm | – |
| SD52 | 520 mm | – |
| SD54 | 540 mm Only With Seat Width 560 mm | – |
| SD56 | 560 mm Only With Seat Width 560 mm | – |

D – JAY Seat Cushions  €

| | | |
|-----------|--|-------|
| JAY020006 | JAY Easy Visco With Flat Base | XXX,- |
| JAY020011 | JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm | XXX,- |
| JAY020013 | JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm | XXX,- |
| JAY020015 | JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm | XXX,- |
| JAY020017 | JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm | XXX,- |

| Knee & Thigh Supports | | € |
|-----------------------|--|-------|
| QOP020342 | Knee Support, Lateral, 120 x 120 mm, Swing Away, Right Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger Receivers Only, QOP090122 Transfer Bar | XXX,- |
| QOP020341 | Knee Support, Lateral, 120 x 120 mm, Swing Away, Left Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger Receivers Only, QOP090122 Transfer Bar | XXX,- |
| QOP020113 | Thigh Support, Lateral, Small, 120 x 120 mm, Removable, Right Not With QOP040035 Side Guard, Wedge Shaped, Padded | XXX,- |
| QOP020114 | Thigh Support, Lateral, Small, 120 x 120 mm, Removable, Left Not With QOP040035 Side Guard, Wedge Shaped, Padded | XXX,- |
| QOP020121 | Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Right Not With QOP040035 Side Guard, Wedge Shaped, Padded | XXX,- |
| QOP020122 | Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Left Not With QOP040035 Side Guard, Wedge Shaped, Padded | XXX,- |

| Armrest Types  | | € |
|---|---|-------|
| QOP040019 | Armrest, Single Post, Flip Back In case of MVD base, not with Seat height 415/seat height 425 & Tilt, Powered, 0° - 30°/ Tilt, Powered, 0° - 50° | - |
| QOP040021 | Armrest, Single Post, Removable In case of MVD base, not with Seat height 415/seat height 425 & Tilt, Powered, 0° - 30°/ Tilt, Powered, 0° - 50° | XXX,- |
| QOP040022 | Armrest, Single Post, SEDEO LITE Arm Pad, Integral Foam, 380 mm Right/Left Included; Side Guard Included; Only With QOP110082 Joystick Mounting, Parallel, Swing Away; Not With QOP090040 Tray Table, Width-Adjustable, Swing Away & QOP090046 Tray Table, Half-Side, Swing Away; Not With Special Controls Mounted On Armrest; Not With QOP090120 Bag For Armrest | XXX,- |
| QOP040123 | Armrest, Reclining, Flip Back | XXX,- |
| QOP040070 | Armrest, Reclining, Multi Adjustable, Flip Back | XXX,- |
| QOP040073 | Armrest, Reclining, Multi Adjustable, Flip Back, With Lock Function Only For Norway, Only For NL Supplied Products | XXX,- |

| Arm Pads  | | € |
|--|--|-------|
| QOP040023 | Arm Pad, Integral Foam, 300 mm, Right | XXX,- |
| QOP040024 | Arm Pad, Integral Foam, 300 mm, Left | XXX,- |
| QOP040025 | Arm Pad, Integral Foam, 400 mm, Right | - |
| QOP040026 | Arm Pad, Integral Foam, 400 mm, Left | - |
| QOP040027 | Arm Pad, Integral Foam, 450 mm, Right With QOP040019 Flip Back Armrest, Seat Width +20 mm | XXX,- |
| QOP040028 | Arm Pad, Integral Foam, 450 mm, Left With QOP040019 Flip Back Armrest, Seat Width +20 mm | XXX,- |

| Comfort Arm Pads | | |
|------------------|--|-------|
| QOP040204 | Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right | XXX,- |
| QOP040203 | Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left | XXX,- |
| QOP040074 | Arm Pad, Comfort, Padded, Stretch Rubber, 300 mm, Right Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040075 | Arm Pad, Comfort, Padded, Stretch Rubber, 300 mm, Left Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040076 | Arm Pad, Comfort, Padded, Stretch Rubber, 400 mm, Right Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040077 | Arm Pad, Comfort, Padded, Stretch Rubber, 400 mm, Left Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040078 | Arm Pad, Comfort, Padded, Removable Cover, 300 mm, Right Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040079 | Arm Pad, Comfort, Padded, Removable Cover, 300 mm, Left Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040080 | Arm Pad, Comfort, Padded, Removable Cover, 400 mm, Right Only For Norway, Only For NL Supplied Products | XXX,- |
| QOP040081 | Arm Pad, Comfort, Padded, Removable Cover, 400 mm, Left Only For Norway, Only For NL Supplied Products | XXX,- |

| Hemiplegia Arm Pads | | |
|---------------------|---|-------|
| QOP040031 | Arm Pad, Hemiplegia, Integral Foam, Right Not With QOP110052 Joystick Mount, Right, QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away, QOP090110 Rear View Mirror, Right; not with QOP110362 Hand Control Mounting Location Right; not with QOP110119/QOP110123 OMNI2 Module Mount Right | XXX,- |
| QOP040032 | Arm Pad, Hemiplegia, Integral Foam, Left Not With QOP110053 Joystick Mount, Right, QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away, QOP090118 Rear View Mirror, Left; not with QOP110361 Hand Control Mounting Location Left; not with QOP110120/QOP110124 OMNI2 Module Mount Left | XXX,- |

Side Guard Pads For Armrest €

| Side Guard Pads For Armrest, Single Post, Flip Back / Removable | | |
|---|--|-------|
| QOP040033 | Side Guard Not With QOP040123 Armrest, Reclining, Flip Back | XXX,- |
| QOP040034 | Side Guard, Integral Foam Seat Width -20 mm, Not With QOP040123 Armrest, Reclining, Flip Back | XXX,- |
| QOP040035 | Side Guard, Wedge Shaped, Padded Seat Width -80 mm, Only With QOP040021 Armrest, Single Post, Removable | XXX,- |
| Side Guard Pads For Armrest, Reclining, Flip Back | | |
| QOP040036 | Side Guard, Mounted On Seat Rail Bracket Only With QOP040123 Armrest, Reclining, Flip Back | XXX,- |

Centre Mount Legrest €

| QOP050086 | Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard In Combination With Lift Or Tilt, Choose QOP050087 Centre Calf Support, Integral Foam Or QOP050092 Calf Pad, Integral Foam, Pair; Limited Knee Angle Possibilities With RWD/HWD | XXX,- |
|-------------|---|-------|
| Pad Options | | |
| QOP050087 | Centre Calf Support, Integral Foam | XXX,- |
| QOP050092 | Calf Pad, Integral Foam, Pair | XXX,- |

Fixed Centre Mount Legrest €

| QOP050140 | Fixed Centre Mount Legrest Not with RWD; In Combination With Lift Or Tilt Choose QPP050092 Calf Pad, Integral Foam, Pair; Only For UK Supplied Products | XXX,- |
|-------------|--|-------|
| Footplates | | |
| QOP050072 | Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm | XXX,- |
| Pad Options | | |
| QOP050092 | Calf Pad, Integral Foam, Pair | XXX,- |

Power Centre Mount Legrest €

| QOP050076 | Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension In Combination With Lift Or Tilt, Choose QOP050092 Calf Pad, Integral Foam, Pair; Not with RWD | XXX,- |
|-------------|--|-------|
| QOP050071 | Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor Not With VR2: QOP110036, QOP110037, QOP110040; In Case Of R-Net: QOP110032 required; In Combination With Lift Or Tilt, Choose QOP050092 Calf Pad, Integral Foam, Pair; Not with RWD | XXX,- |
| QMP050182 | Centre Mount Legrest, Powered, 90° - 180°, 135 mm Extension, Short, LLL 300mm In Combination With Lift Or Tilt, Choose QOP050092 Calf Pad, Integral Foam, Pair; Only For Norway, Only For NL Supplied Products; Not with RWD | |
| Footplates | | |
| QOP050072 | Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm | XXX,- |
| QOP050073 | Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm | XXX,- |
| Pad Options | | |
| QOP050092 | Calf Pad, Integral Foam, Pair | XXX,- |

SEDEO PRO Swing Away Legrest €

| QOP050015 | SEDEO PRO, Legrest, 70°, Swing Away | XXX,- |
|------------------|--|-------|
| QOP050016 | SEDEO PRO, Legrest, Manual Elevating, Swing Away | XXX,- |
| QOP050017 | SEDEO PRO, Legrest, Power Elevating, Swing Away | XXX,- |
| QOP050183 | SEDEO PRO, Legrest, Power Elevating, Swing Away, Short, LLL 300 mm | XXX,- |
| QOP050000 | Omit Legrest | XXX,- |
| Hanger Receivers | | |
| QOP050170 | SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thickness 70 mm And Higher | XXX,- |
| QOP050171 | SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thickness 50 mm And Lower | XXX,- |
| QOP050080 | SEDEO PRO, Hanger Receivers Only, High | XXX,- |
| QOP050172 | SEDEO PRO, Hanger Receivers Only, Low | XXX,- |

Footplates

| | | |
|-----------|--|-------|
| QOP050088 | Footplate, For Seat Width 385 - 400 mm | XXX,- |
| QOP050089 | Footplate, For Seat Width 420 - 460 mm | XXX,- |
| QOP050090 | Footplate, For Seat Width 480 - 520 mm | XXX,- |
| QOP050091 | Footplate, For Seat Width 540 - 560 mm | XXX,- |
| QOP050081 | Footplates, Coupled Choose Also A Footplate; Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away | XXX,- |

Support Options 

| | | |
|-----------|--|-------|
| QOP050082 | Calf Support, Integral Foam, Depth-Adjustable, Pair Not On RWD With QOP050015 SEDEO PRO, Legrest, 70°, Swing Away | XXX,- |
| QOP050083 | Calf Support, Padded, Depth-Adjustable, Pair Not On RWD With QOP050015 SEDEO PRO, Legrest, 70°, Swing Away | XXX,- |
| QOP050084 | Heel Loop | XXX,- |
| QOP050085 | Calf Strap Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away | - |
| QOP050200 | Knee Pad, Integral Foam | XXX,- |

Amputee Support

| | | |
|-----------|--|-------|
| QOP050142 | Amputee Support, Integral Foam, Right | XXX,- |
| QOP050143 | Amputee Support, Integral Foam, Left | XXX,- |
| QOP050013 | Amputee Support, Padded, Right Extra On A Set Of Swing Away Legrests, Attached Lose | XXX,- |
| QOP050014 | Amputee Support, Padded, Left Extra On A Set Of Swing Away Legrests, Attached Lose | XXX,- |

Positioning Supports & Accessories

€

| | | |
|-----------|---|-------|
| QOP090128 | Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm | - |
| QOP090129 | Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm | XXX,- |
| QOP090110 | Rear View Mirror, Right Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away | XXX,- |
| QOP090118 | Rear View Mirror, Left Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away | XXX,- |
| QOP090009 | Crutch Holder Not With QOP030086 Powered Recline, 90° - 130°, QOP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear | XXX,- |
| QOP090119 | Brackets For Bag On Backrest Frame Not With QOP030086 Powered Recline, 90° - 130°, QOP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear | XXX,- |
| QOP090120 | Bag For Armrest Not With QOP040123 Armrest, Reclining, Flip Back | XXX,- |
| QOP090122 | Transfer Bar Not With QOP020342/QMP020341 Knee Support, Lateral, Swing Away, Right/Left & QOP050015-17 SEDEO PRO, Legrest, Swing Away & QOP050080 SEDEO PRO, Hanger Receivers Only | XXX,- |
| QOP030605 | Abduction Pommel Not With (Power) Centre Mount Legrests | XXX,- |
| QOP090205 | Dahl Lockplate Mounted Only For Norway, Only For NL Supplied Products | |

Tray Table (Mounted Opposite Joystick)

| | | |
|-----------|--|-------|
| QOP090040 | Tray Table, Width-Adjustable, Swing Away | XXX,- |
| QOP090046 | Tray Table, Half-Side, Swing Away Not with QOP090126 raised Edges, Acrylic For Tray Table | XXX,- |

Tray Table Accessories

| | | |
|-----------|--------------------------------------|-------|
| QOP090126 | Raised Edges, Acrylic For Tray Table | XXX,- |
|-----------|--------------------------------------|-------|

Joystick / Electronics**Remote Joystick – VR2** 

€

| | | |
|-----------|--|-------|
| QOP110036 | Remote Joystick, VR2 | - |
| QOP110037 | Remote Joystick, VR2, 2 Actuators | XXX,- |
| QOP110040 | Remote Joystick, VR2, 2 Actuators And Lights | XXX,- |

Joystick Mounting  €

| | | |
|-----------|---|-------|
| QOP110050 | Joystick Mounting, Fixed | – |
| QOP110056 | Joystick Mounting, Parallel, Swing Away, Capsulated | XXX,- |
| QOP110082 | Joystick Mounting, Parallel, Swing Away Only With QOP040022 Armrest, Single Post, SEDEO LITE | XXX,- |

Joystick Mounting Location €

| | | |
|-----------|-----------------------------------|---|
| QOP110052 | Joystick Mounting Location, Right | – |
| QOP110053 | Joystick Mounting Location, Left | – |

Joystick Knobs (Optional)  €

| | | |
|-----------|--|-------|
| QOP090505 | Joystick Knob, Golf Ball, Ø 42 mm | XXX,- |
| QOP090506 | Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP090507 | Joystick Knob, Soft Ball, Ø 70 mm | XXX,- |
| QOP090508 | Joystick Knob, Ball, Ø 35 mm | XXX,- |
| QOP090045 | Joystick Knob, Allround, Ball, Ø 48 mm | XXX,- |
| QOP090509 | Joystick Knob, Chin Cup, Ø 37 mm | XXX,- |
| QOP090512 | Joystick Knob, Stick, 80 x 15 mm | XXX,- |
| QOP090511 | Joystick Knob, Mushroom, Ø 45 mm | XXX,- |

Electronic Programmers – VR2  €

| | | |
|-----------|--------------------------------------|-------|
| QOP090065 | Mobility Programmer With Dongle, VR2 | XXX,- |
|-----------|--------------------------------------|-------|

Power Module – VR2 €

| | | |
|-----------|--|-------|
| QOP110002 | Power Module, VR2, 90 A Power Legrests Not Available In Combination With C-Me | XXX,- |
|-----------|--|-------|

Lighting Module – VR2 €

| | | |
|-----------|---|-------|
| QOP110030 | Lighting Module, VR2 If LED Light And Indicator Is Chosen, Then This Is Required | XXX,- |
|-----------|---|-------|

Attendant Controls €

| | | |
|-----------|--|-------|
| QOP110021 | Attendant Control, Dual | XXX,- |
| QOP110020 | Joystick Mounted As Attendant Control Choose Also Remote Joystick | XXX,- |

Attendant Mounting Location

| | | |
|-----------|------------------------------------|---|
| QOP110025 | Attendant Mounting Location, Right | – |
| QOP110026 | Attendant Mounting Location, Left | – |

4-Actuator Control Box (Direct Actuator Control) €

| | | |
|-----------|--|-------|
| QOP110054 | 4-Actuator Control Box Attached To Armrest, Right Not Available With C-Me | XXX,- |
| QOP110055 | 4-Actuator Control Box Attached To Armrest, Left Not Available With C-Me | XXX,- |

Remote Joystick – R-Net  €

| | | |
|-----------|---|-------|
| QOP110015 | Remote Joystick, R-Net, LED, Actuators | – |
| QOP110016 | Remote Joystick, R-Net, LED, Actuators And Lights | – |
| QOP110017 | Remote Joystick, R-Net, Colour LCD | – |
| QOP110018 | Remote Joystick, R-Net, Colour LCD, Heavy Duty | XXX,- |
| QOP110013 | Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included | XXX,- |
| QOP110023 | Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included | XXX,- |
| QOP110019 | Omit Remote Joystick, R-Net | XXX,- |

Joystick Mounting  €

| | | |
|-----------|---|-------|
| QOP110050 | Joystick Mounting, Fixed | – |
| QOP110056 | Joystick Mounting, Parallel, Swing Away, Capsulated | XXX,- |
| QOP110082 | Joystick Mounting, Parallel, Swing Away Only With QOL040022 Armrest, Single Post, SEDEO LITE | XXX,- |

Joystick Mounting Location

€

| | | |
|-----------|-----------------------------------|---|
| QOP110052 | Joystick Mounting Location, Right | – |
| QOP110053 | Joystick Mounting Location, Left | – |

Display Language Selection

€

| | | | | | | | |
|-----------|---------|-----------|------------|-----------|---------|-----------|--------|
| QOP119000 | German | QOP119004 | Portuguese | QOP119009 | Swedish | QOP119013 | Slovak |
| QOP119001 | English | QOP119005 | Italian | QOP119010 | Greek | QOP119014 | Polish |
| QOP119002 | French | QOP119006 | Dutch | QOP119011 | Finnish | QOP119015 | Danish |
| QOP119003 | Spanish | QOP119008 | Norwegian | QOP119012 | Czech | | |

Joystick Knobs (Optional)

€

| | | |
|-----------|---------------------------------------|-------|
| QOP090505 | Joystick Knob, Golf Ball, Ø 42 mm | XXX,- |
| QOP090506 | Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP090507 | Joystick Knob, Soft Ball, Ø 70 mm | XXX,- |
| QOP090508 | Joystick Knob, Ball, Ø 35 mm | XXX,- |
| QOP090045 | Joystick Knob, Allround, Ball Ø 48 mm | XXX,- |
| QOP090509 | Joystick Knob, Chin Cup, Ø 37 mm | XXX,- |
| QOP090512 | Joystick Knob, Stick, 80 x 15 mm | XXX,- |
| QOP090511 | Joystick Knob, Mushroom, Ø 45 mm | XXX,- |

Electronic Programmers – R-Net

€

| | | |
|-----------|----------------------------------|-------|
| QOP090066 | PC Programmer With Dongle, R-Net | XXX,- |
|-----------|----------------------------------|-------|

Power Modules – R-Net

€

| | | |
|-----------|----------------------------|-------|
| QOP110003 | Power Module, R-Net, 90 A | XXX,- |
| QOP110004 | Power Module, R-Net, 120 A | XXX,- |

Seating / Lighting Module Control – R-Net

€

| | | |
|-----------|---|-------|
| QOP110032 | Actuator Control & Lights Module, R-Net, ISM-6L | XXX,- |
|-----------|---|-------|

Attendant Controls

€

| | | |
|-----------|--|-------|
| QOP110021 | Attendant Control, Dual | XXX,- |
| QOP110020 | Joystick Mounted As Attendant Control Choose Also Remote Joystick | XXX,- |

Attendant Mounting Location

| | | |
|-----------|------------------------------------|---|
| QOP110025 | Attendant Mounting Location, Right | – |
|-----------|------------------------------------|---|

Table Top Controls

€

| | | |
|-----------|---|-------|
| QOP110143 | Table Top Control, LCD, Centre Tray Mounted | XXX,- |
|-----------|---|-------|

Tray Table Mounting Location

| | | |
|-----------|-------------------------------------|---|
| QOP110080 | Tray Table Mounting Location, Left | – |
| QOP110081 | Tray Table Mounting Location, Right | – |

#1 Ctrl+5 Switch Box Options

€

| | | |
|-----------|---|-------|
| QOP110043 | #1 Ctrl+5 Switch Box With Free Assignable Functions If combined with attendant controll, Ctr/+5 will be mounted to the armrest | XXX,- |
|-----------|---|-------|

#1 Ctrl+5 Switch Box Mounting Location

| | | |
|-----------|--|---|
| QOP110180 | #1 Ctrl+5 Switch Box Mounting Location, Behind Joystick Not Possible With QOP110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away) | – |
| QOP110181 | #1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick | – |

#1 Ctrl+5 Switch Type

| | | |
|-----------|--|---|
| QOP110175 | #1 Ctrl+5 Switch Type With Push Buttons | – |
| QOP110176 | #1 Ctrl+5 Switch Type With Toggle Switches | – |

Direct Access Only For Norway, Only For NL Supplied Products

| | | |
|-----------|----------------------------------|-------|
| QOP110174 | Direct Access, Behind Joystick | XXX,- |
| QOP110188 | Direct Access, Opposite Joystick | XXX,- |

Special Controls**OMNI2 Input Control Module – R-Net** Choose OMNI2 Mount And Input Device in Sections Below And Choose Switches in Section E  €

| | | |
|-----------|---|-------|
| QOP110116 | OMNI2 Input Control Module, R-Net, Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker | XXX,- |
| QOP110117 | OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker | XXX,- |
| QOP110118 | OMNI2 Input Control Module, R-Net, Scanning Mode, Traction Control System (GYRO Technology), Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker | XXX,- |
| QOP110008 | Programming Service Not Discountable | XXX,- |

OMNI2 Module Mount If Used With Standard Joystick Selection, It Must Be Mounted Opposite Side €

| | | |
|-----------|--|-------|
| QOP110119 | OMNI2 Module, Fixed Mount, Right | – |
| QOP110120 | OMNI2 Module, Fixed Mount, Left | – |
| QOP110123 | OMNI2 Module, Swing Away Mount, Right Not With QOP090110 Rear View Mirror Right | XXX,- |
| QOP110124 | OMNI2 Module, Swing Away Mount, Left Not With QOP090118 Rear View Mirror Left | XXX,- |

Hand Control For Switch Selection See Section E €**Proportional Input Devices**

| | | |
|-----------|--|-------|
| QOP110301 | Hand Control, Switch-It MicroGuide | XXX,- |
| QOP110302 | Hand Control, Switch-It MicroPilot | XXX,- |
| QOP110304 | Hand Control, Switch-It VersaGuide | XXX,- |
| QOP110305 | Hand Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP110306 | Hand Control, Allround Joystick | XXX,- |
| QOP110307 | Hand Control, Allround Joystick Light Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm | XXX,- |
| QOP110308 | Hand Control, Allround Joystick Heavy Duty | XXX,- |
| QOP110309 | Hand Control, Allround Joystick, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module | XXX,- |
| QOP110310 | Hand Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module | XXX,- |
| QOP110311 | Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module | XXX,- |

Proportional Input Devices, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker

| | | |
|-----------|--|-------|
| QOP110044 | Hand Control, Switch-It MicroGuide | XXX,- |
| QOP110045 | Hand Control, Switch-It MicroPilot | XXX,- |
| QOP110047 | Hand Control, Switch-It VersaGuide | XXX,- |
| QOP110048 | Hand Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP110049 | Hand Control, Allround Joystick | XXX,- |
| QOP110059 | Hand Control, Allround Joystick Light Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm | XXX,- |
| QOP110060 | Hand Control, Allround Joystick Heavy Duty | XXX,- |

Hand Control Mounting Platforms €

| | | |
|-----------|--|-------|
| QOP110342 | Hand Control Mounting, Link-It To The Space Disk, Hand Tray Only For Micro Series Joysticks | XXX,- |
|-----------|--|-------|

Hand Control Mounting 

| | | |
|-----------|---|-------|
| QOP110350 | Hand Control Mounting, Link-It To The Armrest, Fixed | XXX,- |
| QOP110353 | Hand Control Mounting, Link-It To The Armrest, Swing Away | XXX,- |

Hand Control Mounting Location

| | | |
|-----------|---------------------------------------|-------|
| QOP110362 | Hand Control Mounting Location, Right | XXX,- |
| QOP110361 | Hand Control Mounting Location, Left | XXX,- |

Head Control For Switch Selection See Section E. Choose A Whitmyer Headrest Mount On Previous Pages 

€

Proportional Input Devices

| | | |
|--|---|-------|
| QOP110370 | Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads | XXX,- |
| QOP110371 | Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads | XXX,- |
| Proportional Input Devices, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker | | |
| QOP110061 | Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads | XXX,- |
| QOP110062 | Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads | XXX,- |

Chin Control 

€

Proportional Input Devices

| | | |
|-----------|---|-------|
| QOP110470 | Chin Control, Switch-It MicroGuide | XXX,- |
| QOP110471 | Chin Control, Switch-It MicroPilot | XXX,- |
| QOP110472 | Chin Control, Switch-It VersaGuide | XXX,- |
| QOP110473 | Chin Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP110474 | Chin Control, Allround Joystick | XXX,- |
| QOP110475 | Chin Control, Allround Joystick Light Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP110476 | Chin Control, Allround Joystick, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module | XXX,- |
| QOP110477 | Chin Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 Module Not Required, Obligated to Have Master R-Net Module; Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |

Proportional Input Devices, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker

| | | |
|-----------|--|-------|
| QOP110070 | Chin Control, Switch-It MicroGuide | XXX,- |
| QOP110071 | Chin Control, Switch-It MicroPilot | XXX,- |
| QOP110072 | Chin Control, Switch-It VersaGuide | XXX,- |
| QOP110073 | Chin Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |
| QOP110074 | Chin Control, Allround Joystick | XXX,- |
| QOP110075 | Chin Control, Allround Joystick Light Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide | XXX,- |

Chin Control Mounting 

| | | |
|-----------|--|-------|
| QOP110480 | Chin Control Mounting, Link-It To Chin Harness, Small If You Need Switches Mounted On The Harness, You Can Only Use Allround Joystick | XXX,- |
| QOP110483 | Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right | XXX,- |
| QOP110484 | Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left | XXX,- |
| QOP110485 | Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right | XXX,- |
| QOP110486 | Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left | XXX,- |
| QOP110487 | Chin Control Mounting, Link-It To Headrest, Swing Away, Left | XXX,- |
| QOP110488 | Chin Control Mounting, Link-It To Headrest, Swing Away, Right | XXX,- |

Link-It Extension Kit

| | | |
|-----------|--|-------|
| QOP110590 | Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips | XXX,- |
| QOP110591 | Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips | XXX,- |

Sip 'N Puff For Switch Selection See Section E

€

| | | |
|---|--|-------|
| QOP110510 | Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest Headrest To Be Selected | XXX,- |
| Sip 'N Puff, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker | | |
| QOP110076 | Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest Headrest To Be Selected | XXX,- |

E – Switches Select 1 Switch And Mounting Position Accordingly. For A Second Switch Please Use #2 Switch Selection Part.
Single Function For Special Input Devices 

€

#1 Switch Selection

| | | |
|-----------|-----------------------------------|-------|
| QOP110611 | #1 Single Switch, Ø 64 mm, Black | XXX,- |
| QOP110612 | #1 Single Switch, Ø 64 mm, Red | XXX,- |
| QOP110613 | #1 Single Switch, Ø 64 mm, Green | XXX,- |
| QOP110523 | #1 Egg Switch – Black | XXX,- |
| QOP110524 | #1 MicroLite Switch – Black | XXX,- |
| QOP110651 | #1 Single Switch, Ø 30 mm, Black | XXX,- |
| QOP110654 | #1 Single Switch, Ø 30 mm, Blue | XXX,- |
| QOP110653 | #1 Single Switch, Ø 30 mm, Green | XXX,- |
| QOP110652 | #1 Single Switch, Ø 30 mm, Red | XXX,- |
| QOP110655 | #1 Single Switch, Ø 30 mm, Yellow | XXX,- |
| QOP110534 | #1 Twister Switch – Black | XXX,- |
| QOP110535 | #1 Twister Switch – Green | XXX,- |
| QOP110536 | #1 Twister Switch – Red | XXX,- |

#1 Switch Mounting 

| | | |
|-----------|--|-------|
| QOP110560 | #1 Switch Mounting, Link-It To The Hand Control Not With HD Allround Joystick Or Space Disk Mount | XXX,- |
| QOP110566 | #1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) Armrest Rail Mount | XXX,- |
| QOP110561 | #1 Switch Mounting, Link-It To The Headrest | XXX,- |
| QOP110563 | #1 Switch Mounting, Link-It To The Seat Rail Operated By Knee, If Two Switches Selected, One Mounted On The Right Side, One On The Left | XXX,- |
| QOP110565 | #1 Switch Mounting, Link-It To Chin Boom No Left / Right Switch Selection Required | XXX,- |
| QOP110569 | #1 Omit Switch Mounting Switches Attached Lose | XXX,- |

#1 Switch Mounting Location

| | | |
|-----------|------------------------------------|-------|
| QOP110572 | #1 Switch Mounting Location, Right | XXX,- |
| QOP110571 | #1 Switch Mounting Location, Left | XXX,- |

#1 Satellite Switches Only For Allround Joystick, Maximum 2 Selectable

| | | |
|-----------|-----------------------------|-------|
| QOP110537 | #1 Satellite Switch – Black | XXX,- |
| QOP110538 | #1 Satellite Switch – Green | XXX,- |
| QOP110539 | #1 Satellite Switch – Red | XXX,- |

#2 Switch Selection

| | | |
|-----------|-----------------------------------|-------|
| QOP110631 | #2 Single Switch, Ø 64 mm, Black | XXX,- |
| QOP110632 | #2 Single Switch, Ø 64 mm, Red | XXX,- |
| QOP110633 | #2 Single Switch, Ø 64 mm, Green | XXX,- |
| QOP110088 | #2 Egg Switch – Black | XXX,- |
| QOP110089 | #2 MicroLite Switch – Black | XXX,- |
| QOP110671 | #2 Single Switch, Ø 30 mm, Black | XXX,- |
| QOP110674 | #2 Single Switch, Ø 30 mm, Blue | XXX,- |
| QOP110673 | #2 Single Switch, Ø 30 mm, Green | XXX,- |
| QOP110672 | #2 Single Switch, Ø 30 mm, Red | XXX,- |
| QOP110675 | #2 Single Switch, Ø 30 mm, Yellow | XXX,- |
| QOP110095 | #2 Twister Switch – Black | XXX,- |
| QOP110096 | #2 Twister Switch – Green | XXX,- |
| QOP110097 | #2 Twister Switch – Red | XXX,- |

#2 Switch Mounting 

| | | |
|-----------|--|-------|
| QOP110098 | #2 Switch Mounting, Link-It To The Hand Control Not With HD Allround Joystick Or Space Disk Mount | XXX,- |
| QOP110099 | #2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) Armrest Rail Mount | XXX,- |
| QOP110100 | #2 Switch Mounting, Link-It To The Headrest | XXX,- |
| QOP110104 | #2 Switch Mounting, Link-It To The Seat Rail Knee, If 2 Switches Selected Only Left And Right Mount | XXX,- |
| QOP110105 | #2 Switch Mounting, Link-It To Chin Boom No Left / Right Switch Selection Required | XXX,- |
| QOP110106 | #2 Omit Switch Mounting Switches Attached Lose | XXX,- |

#2 Switch Mounting Location

| | | |
|-----------|------------------------------------|-------|
| QOP110108 | #2 Switch Mounting Location, Right | XXX,- |
| QOP110107 | #2 Switch Mounting Location, Left | XXX,- |

#2 Satellite Switches Only For Allround Joystick, Maximum 2 Selectable

| | | |
|-----------|-----------------------------|-------|
| QOP110111 | #2 Satellite Switch – Black | XXX,- |
| QOP110112 | #2 Satellite Switch – Green | XXX,- |
| QOP110113 | #2 Satellite Switch – Red | XXX,- |

Multi-Function Switches Used In Combination With Assignable Button Programming

| | | |
|-----------|---|-------|
| QOP110550 | Mini 2 Button (Stereo Plug) Attached Lose, When Ordered. No Link-It Fixation Possible | XXX,- |
| QOP110551 | Mini 2 Button (Dual Mono Plug) Attached Lose, When Ordered. No Link-It Fixation Possible | XXX,- |
| QOP110552 | Mini 4 Button InLine (DB9 Connector) Attached Lose, When Ordered. No Link-It Fixation Possible | XXX,- |
| QOP110554 | 4-Way Toggle Switched Joystick | XXX,- |

Accessories

| | | |
|-----------|---|-------|
| QOP110580 | Mono Extension Cable (3.5 mm Jack Socket) (2 m) | XXX,- |
| QOP110581 | Stereo Extension Cable (3.5 mm Jack Socket) (2 m) | XXX,- |
| QOP110582 | Extension Cable DB9 (1 m) | XXX,- |
| QOP110583 | Y-Cable, 2 Female Mono To 1 Male Stereo | XXX,- |

Safety Devices

| | | |
|-----------|---|-------|
| QOP090147 | Remote Stop Switch (9 m Activation Range, Fit Via OBC Port) Not With VR2 | XXX,- |
|-----------|---|-------|

Drive Wheels & Castor Wheels & Mudguards & Kerb Climber

◀ BACK



// QOP070069



// QOP080213



// QOP010007



// QOP010008



// QOP090038

LED Lights & Chargers

◀ BACK



// QOP110209



// QOP110209



// QOP160101



// QOP160104

SEDEO PRO Back Cushions

◀ BACK



// QOP030109



// QOP030113



// QOP030114



// QOP030115



// QOP030116



// QOP030117

JAY J3 Back

◀ **BACK**



// JAY030012



// JAY030013



// JAY030045



// JAY030014

| Seat Width in mm | Shallow Contour | Mid Contour | Mid Deep Contour / Extra Deep | Deep Contour / Extra Deep |
|--------------------------------|-----------------|-------------|----------------------------------|------------------------------|
| | JAY030012 | JAY030013 | JAY030045/JAY030050 | JAY030014/JAY030051 |
| Maximum User Trunk Width in mm | | | | |
| 360 | 330 | 295 | 280 | 280 |
| 410 | 380 | 345 | 330 | 330 |
| 460 | 430 | 395 | 380 | 380 |
| 510 | 480 | 445 | 430 | 430 |

JAY J3 Back, Spine Align Fit Kit & Positioning Boxes

◀ **BACK**



// JAY090060-JAY090034



// JAY090059-JAY090031



// JAY090035



// JAY090036

JAY J3 Back, Chest Strap

◀ **BACK**



// JAY090040 - JAY090042

Headrest Pads & Headrest Mounts & Push Handles

← BACK




// QOP030430 & QOP030431



// WHT030401-WH030403



// WHT030404-WHT030406



// QOP030433



// WHT030407



// WHT030408



// WHT030409



// WHT030410



// WHT030411



// QOP030204

SEDEO Seat Cushions

← BACK



// QOP020073



// QOP020074



// QOP020075



// QOP020076



// QOP020077



// QOP020078

JAY Seat Cushions

◀ BACK



// JAY020005



// JAY020003



// JAY020006



// JAY020010



// JAY020011



// JAY020018



// JAY020019



// JAY020020

Armrest Types & Arm Pads

◀ BACK



// QOP040019



// QOP040021



// QOP040123




// QOP040027



// QOP040023



// QOP040204



// QOP040070

Fixed Centre Mount Legrest & Centre Mount Legrest

◀ BACK



// QOP050086



// QOP050140

SEDEO PRO Swing Away Legrest & Support Options & Power Centre Mount Legrest

◀ BACK



// QOP050016



// QOP050076

Remote Joystick– R-Net & Joystick Mounting

◀ BACK



// QOP110017



// QOP110013



// QOP110018



// QOP110056

Remote Joystick – VR2

◀ BACK



// QOP110036



// QOP110037

Joystick Knobs

◀ BACK



// QOP090505



// QOP090045



// QOP090507



// QOP090508



// QOP090509



// QOP090512



// QOP090511



// QOP090506

Electronic Programmers – R-Net

◀ BACK



// QOP090066

Power Modules – R-Net & Seating / Lighting Module Control – R-Net

◀ BACK



// QOP110004

Attendant Controls

◀ BACK



// QOP110021

#1 Ctrl+5- Switch Box Options & OMNI2 Input Control Module – R-Net

◀ BACK



// QOP110175



// QOP110176



// QOP090068



// QOP110116

Head Control

◀ BACK



// QOP110401 - QOP110406



// QOP110370 // QOP110371

Chin Control & Hand Control

◀ BACK



// QOP110470 // QOP110301



// QOP110471 // QOP110302



// QOP110472 // QOP110473 //
QOP110304 // QOP110305



// QOP110306 - QOP110309
// QOP110474 - QOP110477



// QOP110475 // QOP110477
// QOP110307 // QOP110310

Table Top Controls

◀ BACK



// QOP110143

Switches

◀ BACK



// QOP110651 - QOP110655
// QOP110671 - QOP110675



// QOP110523



// QOP110524



// QOP110611 - QOP110613
// QOP110631 - QOP110633



// QOP110537 - QOP110539



// QOP110534 - QOP110536



// QOP110551



// QOP110552

Switch Mounting

◀ BACK



// QOP110526 // QOP110566



// QOP110566 // QOP110526
// QOP110529



// QOP110560 // QOP110572
// QOP110013 // QOP110056
// QOP110529 // QOP110526



// QOP110563

Hand Control Mounting

◀ **BACK**



// QOP110353 // QOP110302
// QOP110362



// QOP110353 // QOP110302
// QOP110362 // QOP110342



// QOP110302



// QOP110353 // QOP110361
// QOP110524 // QOP110526
// QOP110560



// QOP110353 // QOP110306
// QOP110362 // QOP110538
// QOP110539

Chin Control Mounting

◀ **BACK**



// QOP110480 // QOP110470
// QOP110471



// QOP110483 // QOP110484



// QOP110483 // QOP110484



// QOP110485 // QOP110486



// QOP110485 // QOP110474



// QOP110487 // QOP110471
// QOP110526 // QOP110529
// QOP110561 // QOP110572
// QOP110571



// QOP110487 // QOP110474

Comments For Customisations


 built4me

Q500 M/F/R/H with SEDEO PRO

€

| | |
|-------------------------------|------------------|
| Base Price | 00,000.00 |
| Total Of Additional Equipment | |
| Packaging And Transport | |
| Exclusive VAT | |
| VAT | |
| Total Price | |

Date

Signature

Day Month Year

**Sunrise Medical**

Ul. Walczaka 114a,66 400,
 Gorsaw, Gorzów Wlkp. Poland
 Phone: + 48 95 737 25 14
 E-Mail: cs.pl@sunrise-medical.pl
 www.SunriseMedical.eu