



OCTANE

SUB4

TITANIUM
TECHNOLOGY:
A NEW LEVEL OF DURABILITY



*Excluding wheels and cushion

OCTANE **SUB4**

THE WORLDS FIRST MADE TO MEASURE
TITANIUM WHEELCHAIR GUARANTEED
UNDER 4KG*

The Octane Sub4 elevates the daily wheelchair to a whole new level of performance and engineering down to the finest detail.

We took a classic design and redefined it even further by combining Sub4 technology with a hand-crafted made to measure titanium grade 9 frame. Resulting, in an ultra-lightweight and rigid wheelchair that is perfectly designed to fit you and your lifestyle.

TITANIUM GRADE 9:

- STRENGTH TO WEIGHT RATIO
- SMOOTH RIDING
- NO CORROSION
- DURABILITY





Frame finish options: polished or brushed



Sub4 technology – guaranteed under 4 kg*

KEY FEATURES

- Made to measure
- Less than 4 kg guaranteed (excluding wheels and cushion)
- Rigid and lightweight
- Titanium grade 9 frame
- Low profile castor arm design
- Airtech Sub4 technology
- Low profile cross braces
- Optimised single tube design
- Max user weight 110 kg



Low profile castor arm design



Single post axle

OCTANE **SUB4**

TECHNICAL DATA

OCTANE SUB4 // Titanium Wheelchair

| | | | |
|--------------------------------|---------------------------|--|---|
| Max. user weight | 110 kg | Total length | 880 mm |
| | | Total height | 900 mm |
| Frame version | Fully welded open frame | Max. pivot width | 880 mm |
| Frame inset | Custom-made | Max. safe slope | 7° (for use of parking brakes) |
| Frame angle | Custom-made | Front wheel | 3" / 4" / 5" |
| Seat width | 250 - 440 mm | Rear wheel | 24" / 25" |
| Seat depth | 250 - 440 mm | Wheelchair weight | From 6.2 kg |
| Seat height front | 380 - 520 mm | Dismanteled weight | From 3.6 kg |
| Seat height rear | 360 - 500 mm | Mass of heaviest part | 2.1 kg |
| Backrest height | 100 - 400 mm | Removable parts | Rear wheels (for transport only) |
| Backrest angle | Custom-made | Colours | - Frame finish options: Brushed frame or polished frame - 5 accent colours (decal on frame), 6 backrest upholstery stitching colours, 9 Spinerger colour options |
| Centre of gravity (COG) | 25 - 125 mm | Intended use as seat in motor vehicle | Not released as seat in motor vehicles |
| Wheel camber | 1° / 2° / 3° | Transport without user | Fix the wheelchair safe in the vehicle (e.g. car or plane) |
| Total width | Max. 680 mm (SW + 240 mm) | Intended use and environment | For people with limited mobility (children & adults); self-propelling or being pushed by an assistant, indoor and outdoor use. |



Sunrise Medical GmbH & Co. KG

Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.de
www.SunriseMedical.eu

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.

