



JAY®

Edition: 11.24
Valid until: 11.25

Balance

JAY Balance International

Invoice to:
Name:
Street:
Town: Postal code:
E-Mail:
Tel: Fax:

Deliver to:
Name:
Street:
Town: Postal code:
E-Mail:
Tel: Fax:

Order date:
Purchase order:
Number required:



Marked for:
Sunrise order no:

- ☐ Order
☐ Calculation

Jay® Balance Standard Features

Fluid pad (Flow Fluid or Cryo Fluid) or Air (single or dual valve) insert to optimize immersion and accommodate positioning needs.
Positioning accessories available. Innovative dual-cover system: water-resistant inner cover, easy to clean, combined with 3 types of machine-washable outer covers: stretch, microclimatic and incontinent.
Available in two heights: 75 mm
2 Year guarantee; 6 months on cover.

Max. user weight:

150 kg (cushion width 34 to 50 cm)

225 kg (cushion width 56 to 60 cm)

USER WEIGHT kg

Mandatory information to ensure correct configuration

Reimbursement code:

Application:

Client at high risk of skin breakdown with or without high postural needs.

Symmetrical to significant asymmetrical posture requiring laterally, forward/rearward stability and/or positioning.

May be unable to reposition or perform an independent pressure relief.

To order a JAY Balance cushion with 75 mm height, follow the steps below to build your part number and then choose your accessories

fixed	Width	Depth	Insert	Outer Cover	Pos. Kit	fixed
JBALC	40	42	FL	M	P	EU
JAY Balance Cushion	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	

Example: JBALC4042FLMPEU

CUSHION WIDTH AND DEPTH - STEP 1 + STEP 2

Standard Cushion (width: 34 - 50 cm; max user weight 150 kg)

Bariatric Cushion (width: 56 - 60 cm; max user weight 225 kg)

Fluid pad (Flow Fluid)
or Air (single or dual valve)

Base Price	831 €
Base Price	879 €


Fluid pad (Cryo Fluid)

Base Price	955 €
Base Price	1.003 €

Width		34	36	38	40	42	44	46	48	50	56	60
Depth	34	○*										
	36	○*	○	○								
	38	○*	○	○	○	○	○					
	40	○*	○	○	○	○	○					
	42	○*	○	○	○	○	○	○	○	○		
	44		○	○	○	○	○	○	○	○		
	46		○	○	○	○	○	○	○	○	○	○
	48					○	○	○	○	○	○	○
	50						○	○	○	○	○	○
PLA Size	56										○	○
	60										○	○
		A			B					C		

* Not available with Cryo Fluid

Jay Balance Measuring Guide



Width X cm - Measurement made from outside edge to outside edge of cushion (including cover)

Depth Y cm - Measurement made from front of seat notch to front topside (or bottom of contracture cut) of cushion (including cover)

WIDTH x DEPTH

STEP 1 STEP 2



For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2025

JAY BALANCE CUSHION

JAY Balance Lead Time - for the below listet cushions the lead time is 5 days; for all others the lead time is 15 days

Seat Width	Seat Depth	Balance Part No.	Seat Width	Seat Depth	Balance Part No.
34	34	JBALC3434FLMPEU	44	46	JBALC4446ADMPEU
	36	JBALC3640FLMPEU		46	JBALC4446ASMPEU
	40	JBALC3440FLMPEU		46	JBALC4446FLSPEU
	36	JBALC3636FLMPEU		46	JBALC4446FLIEU
	40	JBALC3640ASMEU		46	JBALC4446FLIEU
	40	JBALC3640FLMPEU		46	JBALC4446FLMPEU
	40	JBALC3640FLMPEU		46	JBALC4446FLMPEU
	42	JBALC3642FLMPEU		48	JBALC4448FLMPEU
	42	JBALC3642FLMPEU		48	JBALC4448FLMPEU
	44	JBALC3644FLMPEU		50	JBALC4450FLMPEU
	46	JBALC3646FLMPEU		50	JBALC4450FLMPEU
	46	JBALC3646FLMPEU		50	JBALC4450FLMPEU
	40	JBALC4040ADMEU		42	JBALC4642FLIEU
	40	JBALC4040ASMEU		42	JBALC4642FLMPEU
	40	JBALC4040FLIEU		42	JBALC4642FLMPEU
	42	JBALC4042ADMEU		44	JBALC4644FLMPEU
	42	JBALC4042ASMEU		44	JBALC4644FLMPEU
	42	JBALC4042FLMPEU		46	JBALC4646ADMPEU
	44	JBALC4044ASMEU		46	JBALC4646FLSPEU
	44	JBALC4044FLMPEU		46	JBALC4646ASMPEU
	44	JBALC4044FLMPEU		46	JBALC4646ADMPEU
	46	JBALC4046ADMEU		46	JBALC4646ASMEU
	46	JBALC4046FLIEU		46	JBALC4646FLIEU
	46	JBALC4046FLMPEU		46	JBALC4646FLMPEU
	46	JBALC4046FLMPEU		46	JBALC4646FLMPEU
	40	JBALC4240FLIEU		48	JBALC4648FLIEU
	40	JBALC4240FLMPEU		48	JBALC4648FLMPEU
	40	JBALC4240FLMPEU		48	JBALC4648FLMPEU
	42	JBALC4242ADMEU		50	JBALC4650FLIEU
	42	JBALC4242ADMPEU		50	JBALC4650FLMPEU
	42	JBALC4242ASMPEU		50	JBALC4650FLMPEU
	42	JBALC4242FLIEU		50	JBALC4650FLMPEU
	42	JBALC4242FLIEU		50	JBALC4650FLMPEU
	42	JBALC4242FLMPEU		50	JBALC4650FLMPEU
	42	JBALC4242FLMPEU		50	JBALC4650FLMPEU
	42	JBALC4242FLSPEU		50	JBALC4650FLMPEU
	44	JBALC4244ADMEU		42	JBALC4842FLMPEU
	44	JBALC4244FLIEU		44	JBALC4844FLMPEU
	44	JBALC4244FLMPEU		44	JBALC4844FLMPEU
	44	JBALC4244FLMPEU		46	JBALC4846FLIEU
	46	JBALC4246ADMPEU		46	JBALC4846FLMPEU
	46	JBALC4246FLIEU		46	JBALC4846FLMPEU
	46	JBALC4246FLMPEU		48	JBALC4848FLMPEU
	46	JBALC4246FLMPEU		48	JBALC4848FLMPEU
	46	JBALC4246FLMPEU		50	JBALC4850FLIEU
	48	JBALC4248FLIEU		50	JBALC4850FLMPEU
	48	JBALC4248FLMPEU		50	JBALC4850FLMPEU
	48	JBALC4248FLMPEU		50	JBALC4850FLMPEU
	48	JBALC4248FLSPEU		50	JBALC4850FLMPEU
40	40	JBALC4440FLMPEU	46	42	JBALC5042FLMPEU
	40	JBALC4440FLMPEU		42	JBALC5042FLMPEU
	40	JBALC4440FLMPEU		44	JBALC5044FLMPEU
	42	JBALC4442FLIEU		44	JBALC5044FLMPEU
	42	JBALC4442FLMPEU		46	JBALC5046FLMPEU
	42	JBALC4442FLMPEU		46	JBALC5046FLMPEU
	42	JBALC4442FLMPEU		48	JBALC5048FLMPEU
	42	JBALC4442FLMPEU		48	JBALC5048FLMPEU
	44	JBALC4444ADMEU		50	JBALC5048FLMPEU
	44	JBALC4444ADMPEU		50	JBALC5050ASMPEU
	44	JBALC4444ASMEU		50	JBALC5050FLMPEU
	44	JBALC4444ASMEU		50	JBALC5050FLMPEU
	44	JBALC4444FLIEU		50	JBALC5050FLMPEU
	44	JBALC4444FLMPEU		50	JBALC5050FLMPEU
	44	JBALC4444FLMPEU		50	JBALC5050FLMPEU
	44	JBALC4444FLSPEU		50	JBALC5050FLMPEU
46	46	JBALC4446ADMEU	56	46	JBALC5646FLMPEU
	46	JBALC4446ADMEU		56	JBALC5656FLMPEU



For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2025

JAY BALANCE CUSHION

PELVIC LOADING AREA (PLA) INSERT* - STEP 3

	Select	Part #
FLUID - Fluid Pad Option	☐	FL
CRYO FLUID - Fluid Pad Option	☐	CFL
AIR Insert Option - Single Valve**	☐	AS
AIR Insert Pad Option - Dual Valve**	☐	AD



Fluid - FL



Cryo Fluid - CFL

INSERT

STEP 3

(FL, CFL AS or AD)



Air single valve - AS



Air dual valve - AD

* The insert pad sizes supplied are based on PLA size and automatically determined when you select the cushion (see matrix)

**Air Bladder Options shipped with instructions, air pump, patch kit and nozzle fittings

OUTER COVER OPTION - STEP 4

	Select	Part #
Microclimatic	☐	M
Stretch	☐	S
Incontinent	☐	I

OUTER COVER

STEP 4

(M, S or I)

All cushions are delivered with an water resistant inner cover.

POSITIONING KIT * - STEP 5

If you wish your Jay Balance to have positioning components, please write "P" in the Step 5 box. Otherwise, leave the box empty to get a standard cushion (without positioning components). When ordering a positioning version, you will get the following positioning kit:

Positioning Kit (containing each 1 of the components mentioned below)

Lateral Thigh 2.5 cm Left
Lateral Thigh 2.5 cm Right
Medial Thigh 2.5 cm (1") (each)
Lateral Pelvic 2.5 cm (2") Left
Lateral Pelvic 2.5 cm (2") Right

Part #
P

Charge
68 €

Lateral Pelvis LH

Lateral Pelvis RH



Lateral Thigh LH - 2,5 cm

Medial Thigh - 2,5 cm

Lateral Thigh RH - 2,5 cm

POS. KIT

STEP 5

(P or blank)

* all positioning components are delivered with a stretchable water resistant inner cover to accommodate these components onto the Balance cushion.

JAY Balance Smart Part Number

When ordering a JAY Balance cushion with a height of 75 mm, please follow the next Notes

Width	Depth	Insert	Outer Cover	Pos. Kit		
JBALC	40	42	FL	M	P	EU
fixed	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	fixed

The part # above (JBALC4042FLMPEU) describes a Jay Balance cushion with 40 width and 42 depth, fluid insert, microclimatic outer cover and positioning components.

Additional explanations:

When selecting the PLA insert and the positioning option for your cushion, the water resistant inner cover will be automatically selected (here in this example it will be an inner cover which accommodates fluid insert and positioning components).

If any modifications of your cushion is needed please contact your Sunrise Medical authorized Customer Service.

JAY BALANCE COSTS

	Unit price	Total price
Cushion Price		
Option Price		
Total Cost		

Sunrise Medical Customer service
Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg
Tel. +497253 980-400
Fax. +497253 980-222

For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2025

Spare Parts	Charge
Foam Base	Charge

Cushion base (ww = width; dd = depth)		
☐	BJBALCwwddEU Base assembly for 75 mm height cushion	252 €
Covers		Charge
Spare Outer Cover (ww = width; dd = depth)		
☐	CJBALCwwddSEU Stretch outer cover for standard fluid and air cushion	134 €
☐	CJBALCwwddMEU Microclimatic outer cover for standard fluid and air cushion	134 €
☐	CJBALCwwddIEU Incontinent outer cover for standard and air cushion	134 €
☐	CCJBALCwwddSEU Stretch outer cover for Cryo fluid cushion	134 €
☐	CCJBALCwwddMEU Microclimatic outer cover for Cryo fluid cushion	134 €
☐	CCJBALCwwddIEU Incontinent outer cover for Cryo fluid	134 €
Spare Inner Cover for 75 mm height cushion (ww = width; dd = depth)		
☐	CIFLJBALCwwddEU Inner cover for standard and Cryo fluid cushion	130 €
☐	CIFLPJBALCwwddEU Inner cover for standard and Cryo Fluid and positioning cushion	130 €
☐	CIAJBALCwwddEU Inner cover for air cushion	130 €
☐	CIAPJBALCwwddEU Inner cover for air and positioning cushion	130 €
Note: Please note that inner and outer spare covers fit to both the new 75 mm high Balance cushion and the former 85 mm high Balance cushion.		
Solid Seat (ww = width; dd = depth)		
☐	JSIwwdd Solid Seat insert (wood plate)	93 €
Fluid & Air inserts *		Charge
PLA A		
☐	FLJBALCAEU Fluid insert for PLA size A	478 €
☐	CFLJBALAEU Cryo Fluid insert for PLA size A	603 €
☐	J3CSAASEU Air insert - single valve - for PLA size A	478 €
☐	J3CSAADEU Air insert - dual valve - for PLA size A	507 €
PLA B		
☐	FLJBALCBEU Fluid insert for PLA size B	478 €
☐	CFLJBALBEU Cryo Fluid insert for PLA size B	603 €
☐	J3CSBASEU Air insert - single valve - for PLA size B	478 €
☐	J3CSBADEU Air insert - dual valve - for PLA size B	507 €
PLA C		
☐	FLJBALCCEU Fluid insert for PLA size C	478 €
☐	CFLJBALCEU Cryo Fluid insert for PLA size C	603 €
☐	J3CSCASEU Air insert - single valve - for PLA size C	478 €
☐	J3CSCADEU Air insert - dual valve - for PLA size C	507 €
* Note: A Balance cushion with Air insert has a different inner cover as with fluid insert. In case you want to change the insert for an existing Balance cushion, you need to order a suitable inner cover.		
 		For more information visit our webpage www.sunrisemedical.com © Sunrise Medical 2025

Spare Parts			
JAY Balance Positioning kit *			Charge
	for depth 34-40 cm	for depth 42-60 cm	
Lateral Thigh Left 2.5cm Left	☐ JBALCLTLH-S	☐ JBALCLTLH	29 €
Lateral Thigh Left 2.5cm Right	☐ JBALCLTRH-S	☐ JBALCLTRH	29 €
Medial Thigh 2.5cm (each)	☐ JBALCMT-S	☐ JBALCMT	29 €
Lateral Pelvic 2.5 cm Left	☐ JBALCLPLH-S	☐ JBALCLPLH	35 €
Lateral Pelvic 2.5 cm Right	☐ JBALCLPRH-S	☐ JBALCLPRH	35 €
Positioning kit	☐ JBALCPKIT-S	☐ JBALCPKIT	139 €
(Positioning kit is containing 1 of each of the components above)			
* Note: A Balance Standard cushion has a different inner cover as the Positioning version. In case you want to change from Standard to Positioning (and versa vice) you need to order a suitable inner cover			
JAY Balance Obliquity Pads			Charge
Order 1 or 2 obliquity pads per side, depending on the level of obliquity that needs to be corrected. The pad has hook/loop on it and can be stacked (maximum 2).			
Obliquity pads can be stacked the following way: for PLA size A: order the obliquity pad size A - one side (RH or LH) - and additionally order the size B (RH or LH). For PLA size B, order the obliquity pad size B - one side (RH or LH) - and additionally order the size C (RH or LH). For PLA size C, order the obliquity pad size C - one side (RH or LH) - and additionally order the size D (RH or LH). Please notice that size D can only be ordered in combination with size C otherwise it will not fit in the PLA.			
☐ JBALCOP-ALH	Obliquity pad (left) for PLA A		52 €
☐ JBALCOP-BLH	Obliquity pad (left) for PLA B		52 €
☐ JBALCOP-CLH	Obliquity pad (left) for PLA C		52 €
☐ JBALCOP-DLH	Additional obliquity pad (left) that can be stacked over obliquity pad for PLA C		52 €
☐ JBALCOP-ARH	Obliquity pad (right) for PLA A		52 €
☐ JBALCOP-BRH	Obliquity pad (right) for PLA B		52 €
☐ JBALCOP-CRH	Obliquity pad (right) for PLA C		52 €
☐ JBALCOP-DRH	Additional obliquity pad (right) that can be stacked over obliquity pad for PLA C		52 €
Note: Cushions or inner covers purchased before February 2017 do not contain appropriate velcro attachment for the obliquity pads.			
 			
Sunrise Medical Customer service Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg Tel. +497253 980-400 Fax. +497253 980-222			For more information visit our webpage www.sunrisemedical.com © Sunrise Medical 2025