

# Cirrus G5

Tilt-in-space wheelchairs designed for pressure relief and positioning

## ADJUSTMENTS - MADE EASY!



Easy to operate sliding seat for safe and easy transfer in and out of the wheelchair. By sliding the seat backwards, the user is positioned correctly in the wheelchair, reducing the need for heavy lifting.



Individual adjustments and re-setting of the wheelchair are simplified through the clear measuring scales.



Easy to use and identify handles and grips, making adjusting the chair simple.

Breezy's range of tilting wheelchairs are a great option for more dependent wheelchair users with high seating comfort and/or changing needs. Designed with both users and caregivers in mind, they offer full adjustability and excellent positioning options, allowing for quick and easy personalisation to suit individual need. Due to the tilt function, the seating position can easily be adapted to activity and change in pressure distribution.

# Breezy® Cirrus G5

## Cirrus'ly Great Value.....Seriously Great Features

The Breezy Cirrus G5 Tilt-in-space wheelchair is designed with both client and carer in mind to provide ultimate comfort and flexibility.

The wheelchair features back recline and tilt which provides an easy and comfortable way to change position throughout the day whilst assisting in pressure management.

The heritage of the Cirrus dates back to 1987. Now, a new age has been born with the Cirrus G5, which offers a very competitively priced inclusive specification with a multitude of options to allow you to meet higher postural seating needs - all available on the standard prescription form.



16", 22" and 24" rear wheel options available

## Technical Specification

<b>Seat Width:</b>	420, 480, 540mm
<b>Seat Depth:</b>	390 - 550mm
<b>Rear Wheel:</b>	16, 22 and 24"
<b>Seat Height:</b>	380 - 550mm
<b>Backrest Height:</b>	550 and 650mm
<b>Overall Width:</b>	Seat width + 220mm (+ 240 with drum brake)
<b>Seat Tilt:</b>	0 - 33 degrees
<b>Back Recline:</b>	0 - 35 degrees
<b>Max. User Weight:</b>	135 kg (21 stone)



With a full 33° tilt and 35° backrest adjustment, the Breezy Cirrus G5 provides an easy and comfortable way to change position throughout the day - all made possible whilst the user is sitting in the wheelchair without the need for any adjustment tools!

## Standard Features

The standard specification of the Breezy Cirrus G5 includes a multitude of features and options.

If you require alternative selections, or need to meet higher clinical needs, simply select the appropriate options on the prescription form.



## On chair adjustment available on the Cirrus G5



Option of flip back or single post height adjustable armrest. 3cm seat width adjustment available with single post



Measuring markers for ease of adjustment



Built in adjustment for legrests



Buckled locked tension adjustable straps for positioning



Contrast colour handles for ease of use

Ease of adjustment with the Cirrus G5 allows you to achieve greater comfort and meet postural needs with the client in the chair.



### Armrests options available

- Single post tool height adjustable, removable with seat width adjustment.
- Single post flip back height adjustable, removable which follows backrest. There are also options of armpad covers and hemi arm pads if required.



### Seating options available

Option of seat plate with pelvic well and a variety of different cushions, including the standard Nylon with PU coating and the Sedeo Pro seat cushions and our extensive JAY range. The prescription form also offers a variety of different seat accessories such as two and four point positioning harnesses and pommels. Most importantly, you have the ability to have a sliding or fixed seat board depending on your client's needs and their environment.



### Backrest options available

The standard backrest is tension adjustable and includes swing-away lateral supports. You are able to choose JAY backs or fit third party backrests for more complex seating needs. Our new Sedeo Lite backrest is also compatible.



### Legrest options available

You can achieve an angle of between 0-90 degree with our ELRs. There is also the option to have tool adjustable legrests that can be altered between 70, 80 and 90 degrees. As well as many other accessories including foot boxes, cast supports, amputee pads, knee pads and calf supports.



### Other options available

There are alternative headrest options to suit your requirements, along with a variety of sizing options to choose from and various options of push handles to suit the carer's needs.