



JAY®

Balance

JAY Balance_OF_EU Template_2022-02-22

Invoice to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Deliver to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:



Order date:
 Purchase order:
 Number required:

Marked for:
 Sunrise order no:

- Order
- Calculation

JAY® Balance Standard Features

Fluid pad (Flow Fluid or Cryo Fluid) or Air (single or dual valve) insert to optimize immersion and accommodate positioning needs.
 Positioning accessories available. Innovative dual-cover system: water-resistant inner cover, easy to clean, combined with 3 types of machine-washable outer covers: stretch, microclimatic and incontinent.
 Available in two heights: 75 mm and 100 mm
 2 Year guarantee; 6 months on cover.

Max. user weight:
 150 kg (cushion width 34 to 50 cm)
 225 kg (cushion width 56 to 60 cm)

USER WEIGHT kg
 Mandatory information to ensure correct configuration

Reimbursement code:

Application:

Client at high risk of skin breakdown with or without high postural needs.
 Symmetrical to significant asymmetrical posture requiring laterally, forward/rearward stability and/or positioning.
 May be unable to reposition or perform an independent pressure relief.

To order a JAY Balance cushion with 75 mm height, follow the steps below to build your part number and then choose your accessories

fixed	Width	Depth	Insert	Outer Cover	Pos. Kit	fixed
JBALC	40	42	FL	M	P	EU
JAY Balance Cushion	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	

Example: JBALC4042FLMPEU

To order a JAY Balance cushion with 100 mm height, follow the steps below to build your part number and then choose your accessories

fixed	Height	Width	Depth	Insert	Outer Cover	Pos. Kit	fixed
JBALC	DC	40	42	FL	M	P	EU
JAY Balance Cushion	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	STEP 6	
Deep Contour							

Example: JBALCDC4042FLMPEU

CUSHION WIDTH AND DEPTH - STEP 1 + STEP 2

Standard Cushion (width: 34 - 50 cm; max user weight 150 kg)
 Bariatric Cushion (width: 56 - 60 cm; max user weight 225 kg)


Width

cm	34	36	38	40	42	44	46	48	50	56	60
34	○*										
36	○*	○	○								
38	○*	○	○	○	○	○					
40	○*	○	○	○	○	○					
42	○*	○	○	○	○	○	○	○	○		
44		○	○	○	○	○	○	○	○		
46		○	○	○	○	○	○	○	○	○	○
48					○	○	○	○	○	○	○
50						○	○	○	○	○	○
56										○	○
60										○	○
PLA Size	A			B				C			

Depth

* Not available with Cryo Fluid

JAY Balance Measuring Guide



Width X cm - Measurement made from outside edge to outside edge of cushion (including cover)

Depth Y cm - Measurement made from front of seat notch to front topside (or bottom of contracture cut) of cushion (including cover)

WIDTH x DEPTH
 STEP 1 STEP 2



For more information visit our webpage
www.sunrisemedical.com
 © Sunrise Medical 2022

JAY BALANCE CUSHION

JAY Balance Lead Time - for the below listet cushions the lead time is 5 days; for all others the lead time is 15 days

Seat Width	Seat Depth	Balance Part No.		
34	34	JBALC3434FLMPEU		
	36	JBALC3640FLMEU		
	40	JBALC3440FLMEU		
36	36	JBALC3636FLMEU		
	40	JBALC3640ASMEU		
		JBALC3640FLMEU		
		JBALC3640FLMPEU		
	42	JBALC3642FLMEU		
	44	JBALC3644FLMPEU		
46	JBALC3646FLMEU JBALC3646FLMPEU			
40	40	JBALC4040ADMEU		
		JBALC4040ASMEU		
		JBALC4040FLIEU		
	42	JBALC4042ADMEU		
		JBALC4042ASMEU		
		JBALC4042FLMPEU		
	44	JBALC4044ASMEU		
		JBALC4044FLMEU JBALC4044FLMPEU		
46	JBALC4046ADMEU			
	JBALC4046FLIEU			
	JBALC4046FLMEU JBALC4046FLMPEU			
42	40	JBALC4240FLIEU		
		JBALC4240FLMEU		
		JBALC4240FLMPEU		
	42	JBALC4242ADMEU JBALC4242ADMPEU JBALC4242ASMPEU JBALC4242FLIEU JBALC4242FLIPEU JBALC4242FLMEU JBALC4242FLMPEU JBALC4242FLSPEU		
			44	JBALC4244ADMEU
				JBALC4244FLIEU
				JBALC4244FLMEU JBALC4244FLMPEU
			46	JBALC4246ASMEU
				JBALC4246ADMEU
	JBALC4246FLIEU JBALC4246FLIPEU JBALC4246FLMEU JBALC4246FLMPEU			
	48	JBALC4248FLIEU		
		JBALC4248FLIPEU		
		JBALC4248FLMEU JBALC4248FLMPEU JBALC4248FLSEU		
	44	40	JBALC4440FLMEU	
JBALC4440FLIEU				
JBALC4440FLMPEU				
42		JBALC4442FLIPEU		
		JBALC4442FLIEU		
		JBALC4442FLMEU JBALC4442FLMPEU		
44		JBALC4444ADMEU		
		JBALC4444ADMPEU		
		JBALC4444ASMPEU		
		JBALC4444ASMEU		
46	JBALC4444FLIPEU			
	JBALC4444FLIEU			
	JBALC4444FLMEU JBALC4444FLMPEU JBALC4444FLSPEU			
46	JBALC4446ADMEU			

Seat Width	Seat Depth	Balance Part No.	
44	46	JBALC4446ADMPEU	
		JBALC4446ASMPEU	
		JBALC4446FLSPEU	
	48	JBALC4446FLIPEU	
		JBALC4446FLIEU	
		JBALC4446FLMEU JBALC4446FLMPEU	
50	JBALC4448FLMEU		
	JBALC4448FLMPEU		
	JBALC4450FLMEU		
46	42	JBALC4450FLIPEU	
		JBALC4450FLIEU	
		JBALC4450ASMPEU JBALC4450FLMPEU	
	44	JBALC4642FLIEU	
		JBALC4642FLMEU	
		JBALC4642FLMPEU	
	46	46	JBALC4644FLMEU
			JBALC4644FLMPEU
			JBALC4646ADMEU
			JBALC4646FLSPEU
48	46	JBALC4646ASMPEU	
		JBALC4646ADMPEU	
		JBALC4646ASMEU	
50	46	JBALC4646FLIEU	
		JBALC4646FLIPEU	
		JBALC4646FLMEU JBALC4646FLMPEU	
48	48	JBALC4648FLIEU	
		JBALC4648FLMEU	
		JBALC4648FLMPEU	
50	42	JBALC4650FLIEU	
		JBALC4650FLMEU	
		JBALC4650FLIPEU	
	44	JBALC4650FLMPEU	
		JBALC4650FLSEU	
		JBALC4842FLMPEU	
48	44	JBALC4844FLMPEU	
		JBALC4844FLMEU	
	46	JBALC4846FLIEU	
		JBALC4846FLMEU	
		JBALC4846FLMPEU	
50	48	JBALC4848FLMEU	
		JBALC4848FLMPEU	
		JBALC4850FLIEU	
56	46	JBALC4850FLMEU	
		JBALC4850FLMPEU	
		JBALC4850FLSPEU	
	50	42	JBALC5042FLMEU
			JBALC5042FLMPEU
			JBALC5044FLMEU
56	44	JBALC5044FLMPEU	
		JBALC5046FLIEU	
		JBALC5046FLMEU	
50	46	JBALC5046FLMPEU	
		JBALC5048FLMEU	
		JBALC5048FLMPEU	
56	50	JBALC5050ASMPEU	
		JBALC5050FLIEU	
		JBALC5050FLMEU JBALC5050FLMPEU	
56	46	JBALC5646FLMEU	
		JBALC5650FLMEU	
	50	JBALC5650ASMEU	
		JBALC5650FLMPEU	
56	JBALC5656FLMPEU		
	JBALC5656FLSPEU		




For more information
visit our webpage

© Sunrise Medical 2022

JAY BALANCE CUSHION

PELVIC LOADING AREA (PLA) INSERT* - STEP 3

	Select	Part #	Charge
FLUID - Fluid Pad Option	<input type="checkbox"/>	FL	-
CRYO FLUID - Fluid Pad Option	<input type="checkbox"/>	CFL	-
AIR Insert Option - Single Valve**	<input type="checkbox"/>	AS	add XX €
AIR Insert Pad Option - Dual Valve**	<input type="checkbox"/>	AD	add XX €



INSERT

STEP 3

(FL,CFL AS or AD)

* The insert pad sizes supplied are based on PLA size and automatically determined when you select the cushion (see matrix)

**Air Bladder Options shipped with instructions, air pump, patch kit and nozzle fittings

OUTER COVER OPTION - STEP 4

	Select	Part #	Charge
Microclimatic	<input type="checkbox"/>	M	-
Stretch	<input type="checkbox"/>	S	-
Incontinent	<input type="checkbox"/>	I	-

OUTER COVER

STEP 4

(M, S or I)

All cushions are delivered with an water resistant inner cover.

POSITIONING KIT * - STEP 5

If you wish your JAY Balance to have positioning components, please write "P" in the Step 5 box. Otherwise, leave the box empty to get a standard cushion (without positioning components). When ordering a positioning version, you will get the following positioning kit:

Positioning Kit (containing each 1 of the components mentioned below)	Part #	Charge
	P	€
Lateral Thigh 2.5 cm Left		
Lateral Thigh 2.5 cm Right		
Medial Thigh 2.5 cm (1") (each)		
Lateral Pelvic 2.5 cm (2") Left		
Lateral Pelvic 2.5 cm (2") Right		



POS. KIT

STEP 5

(P or blank)

* all positioning components are delivered with a stretchable water resistant inner cover to accommodate these components onto the Balance cushion.

JAY Balance Smart Part Number

When ordering a JAY Balance cushion with a height of 75 mm, please follow the next steps

JBALC	Width	Depth	Insert	Outer Cover	Pos. Kit	EU
fixed	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	fixed

Notes

The part # above (JBALC4042FLMPEU) describes a Jay Balance cushion with 40 width and 42 depth, fluid insert, microclimatic outer cover and positioning components.

When ordering a JAY Balance cushion with a height of 100 mm, please follow the next steps

JBALC	Height	Width	Depth	Insert	Outer Cover	Pos. Kit	EU
fixed	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	STEP 6	fixed

Additional explanations:

When selecting the PLA insert and the positioning option for your cushion, the water resistant inner cover will be automatically selected (here in this example it will be an inner cover which accommodates fluid insert and positioning components).

If any modifications of your cushion is needed please contact your Sunrise Medical authorized Customer Service.

JAY BALANCE COSTS	Unit price	Total price
Cushion Price		
Option Price		
Total Cost		
	For more information visit our webpage www.sunrisemedical.eu © Sunrise Medical 2022	

Spare Parts	Charge
-------------	--------

Foam Base	Charge
Cushion base (ww = width; dd = depth)	
<input type="checkbox"/> BJBALCwwddEU Base assembly for 75 mm height cushion	€
<input type="checkbox"/> BJBALCDCwwddEU Base assembly for 100 mm height cushion	€
<input type="checkbox"/> MJBALCwwddEU Base assembly for 85 mm height cushion (via JAY Your Way)	€

Covers	Charge
--------	--------

Spare Outer Cover for 75 mm height cushion (ww = width; dd = depth)	
<input type="checkbox"/> CJBALCwwddSEU Stretch outer cover for standard fluid and air cushion	€
<input type="checkbox"/> CJBALCwwddMEU Microclimatic outer cover for standard fluid and air cushion	€
<input type="checkbox"/> CJBALCwwddIEU Incontinent outer cover for standard and air cushion	€
<input type="checkbox"/> CCJBALCwwddSEU Stretch outer cover for Cryo fluid cushion	€
<input type="checkbox"/> CCJBALCwwddMEU Microclimatic outer cover for Cryo fluid cushion	€
<input type="checkbox"/> CCJBALCwwddIEU Incontinent outer cover for Cryo fluid	€
Spare Outer Cover for 100 mm height cushion (ww = width; dd = depth)	
<input type="checkbox"/> CJBALCDCwwddSEU Stretch outer cover for standard fluid and air cushion	€
<input type="checkbox"/> CJBALCDCwwddMEU Microclimatic outer cover for standard fluid and air cushion	€
<input type="checkbox"/> CJBALCDCwwddIEU Incontinent outer cover for standard and air cushion	€
<input type="checkbox"/> CCJBALCDCwwddSEU Stretch outer cover for Cryo fluid cushion	€
<input type="checkbox"/> CCJBALCDCwwddMEU Microclimatic outer cover for Cryo fluid cushion	€
<input type="checkbox"/> CCJBALCDCwwddIEU Incontinent outer cover for Cryo fluid	€
Spare Inner Cover for 75 mm height cushion (ww = width; dd = depth)	
<input type="checkbox"/> CIFLJBALCwwddEU Inner cover for standard and Cryo fluid cushion	€
<input type="checkbox"/> CIFLPJBALCwwddEU Inner cover for standard and Cryo Fluid and positioning cushion	€
<input type="checkbox"/> CIAJBALCwwddEU Inner cover for air cushion	€
<input type="checkbox"/> CIAPJBALCwwddEU Inner cover for air and positioning cushion	€
Spare Inner Cover for 100 mm height cushion (ww = width; dd = depth)	
<input type="checkbox"/> CIFLJBALCDCwwddEU Inner cover for standard and Cryo fluid cushion	€
<input type="checkbox"/> CIFLPJBALCDCwwddEU Inner cover for standard and Cryo Fluid and positioning cushion	€
<input type="checkbox"/> CIAJBALCDCwwddEU Inner cover for air cushion	€
<input type="checkbox"/> CIAPJBALCDCwwddEU Inner cover for air and positioning cushion	€
Note: Please note that inner and outer spare covers fit to both the new 75 mm high Balance cushion and the former 85 mm high Balance cushion.	

Solid Seat	Charge
<input type="checkbox"/> JSIwwdd Solid Seat insert (wood plate)	€

Spare Parts

Fluid & Air inserts

Charge

PLA A

<input type="checkbox"/>	FLJBALCAEU	Fluid insert for PLA size A	€
<input type="checkbox"/>	FLJBALCDCAEU	Deep Contour Fluid insert for PLA size A	€
<input type="checkbox"/>	CFLJBALAEU	Cryo Fluid insert for PLA size A	€
<input type="checkbox"/>	CFLJBALDCAEU	Deep Contour Cryo Fluid insert for PLA size A	€
<input type="checkbox"/>	J3CSAASEU	Air insert - single valve - for PLA size A	€
<input type="checkbox"/>	J3CSAADEU	Air insert - dual valve - for PLA size A	€

PLA B

<input type="checkbox"/>	FLJBALCBEU	Fluid insert for PLA size B	€
<input type="checkbox"/>	FLJBALCDCBEU	Deep Contour Fluid insert for PLA size B	€
<input type="checkbox"/>	CFLJBALBEU	Cryo Fluid insert for PLA size B	€
<input type="checkbox"/>	CFLJBALDCBEU	Deep Contour Cryo Fluid insert for PLA size B	€
<input type="checkbox"/>	J3CSBASEU	Air insert - single valve - for PLA size B	€
<input type="checkbox"/>	J3CSBADEU	Air insert - dual valve - for PLA size B	€

PLA C

<input type="checkbox"/>	FLJBALCCEU	Fluid insert for PLA size C	€
<input type="checkbox"/>	CFLJBALDCCEU	Deep Contour Fluid insert for PLA size C	€
<input type="checkbox"/>	CFLJBALCEU	Cryo Fluid insert for PLA size C	€
<input type="checkbox"/>	FLJBALCDCCEU	Deep Contour Cryo Fluid insert for PLA size C	€
<input type="checkbox"/>	J3CSCASEU	Air insert - single valve - for PLA size C	€
<input type="checkbox"/>	J3CSCADEU	Air insert - dual valve - for PLA size C	€

Note: A Balance cushion with Air insert has a different inner cover as with fluid insert. In case you want to change the insert for an existing Balance cushion, you need to order a suitable inner cover.

JAY Balance Positioning kit

Charge

	for depth 34-40 cm	for depth 42-60 cm	
Lateral Thigh Left 2.5cm Left	<input type="checkbox"/> JBALCLTLH-S	<input type="checkbox"/> JBALCLTLH	€
Lateral Thigh Left 2.5cm Right	<input type="checkbox"/> JBALCLTRH-S	<input type="checkbox"/> JBALCLTRH	€
Medial Thigh 2.5cm (each)	<input type="checkbox"/> JBALCMT-S	<input type="checkbox"/> JBALCMT	€
Lateral Pelvic 2.5 cm Left	<input type="checkbox"/> JBALCLPLH-S	<input type="checkbox"/> JBALCLPLH	€
Lateral Pelvic 2.5 cm Right	<input type="checkbox"/> JBALCLPRH-S	<input type="checkbox"/> JBALCLPRH	€
Positioning kit	<input type="checkbox"/> JBALCPKIT-S	<input type="checkbox"/> JBALCPKIT	€

(Positioning kit is containing 1 of each of the components above)

Note: A Balance cushion with Positioning options has a different inner cover as the Positioning version. In case you want to change from Standard to Positioning (and versa vice) you need to order a suitable inner cover

JAY Balance Obliquity Pads

Charge

Order 1 or 2 obliquity pads per side, depending on the level of obliquity that needs to be corrected. The pad has hook/loop on it and can be stacked (maximum 2).

Obliquity pads can be stacked the following way: for PLA size A: order the obliquity pad size A - one side (RH or LH) - and additionally order the size B (RH or LH). For PLA size B, order the obliquity pad size B - one side (RH or LH) - and additionally order the size C (RH or LH). For PLA size C, order the obliquity pad size C - one side (RH or LH) - and additionally order the size D (RH or LH). Please notice that size D can only be ordered in combination with size C otherwise it will not fit in the PLA.

<input type="checkbox"/>	JBALCOP-ALH	Obliquity pad (left) for PLA A	€
<input type="checkbox"/>	JBALCOP-BLH	Obliquity pad (left) for PLA B	€
<input type="checkbox"/>	JBALCOP-CLH	Obliquity pad (left) for PLA C	€
<input type="checkbox"/>	JBALCOP-DLH	Additional obliquity pad (left) that can be stacked over obliquity pad for PLA C	€
<input type="checkbox"/>	JBALCOP-ARH	Obliquity pad (right) for PLA A	€
<input type="checkbox"/>	JBALCOP-BRH	Obliquity pad (right) for PLA B	€
<input type="checkbox"/>	JBALCOP-CRH	Obliquity pad (right) for PLA C	€
<input type="checkbox"/>	JBALCOP-DRH	Additional obliquity pad (right) that can be stacked over obliquity pad for PLA C	€

Note: Cushions or inner covers purchased before February 2017 do not contain appropriate velcro attachment for the obliquity pads.

For more information
visit our webpage
www.sunrisemedical.eu
© Sunrise Medical 2022