

Q500® M/F/R

with SEDEO® PRO



07/2024

ORDER FORM


Errors and omissions excepted.

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me®	User Weight Capacity: 160 kg	
Purchase Order No.:			User Weight: kg	
Art. No.: Q500 M: XXX, Q500 F: XXX, Q500 R: XXX,			Reimbursement Code: xxx	
Order Date:	Day	Month	Year	

Base Price € **XX.XXX,-**

Client Measurements

a) Top Of Head To Seat	mm
b) Top Of Shoulder To Seat	mm
c) Underarm To Seat	mm
d) Elbow To Seat	mm
e) Top Of Pelvis To Seat	mm
f) Shoulder Width	mm
g) Trunk Width	mm
i) Hip Width	mm
j) Actual Sitting Depth	mm
m) Foot Length	mm
n) Lower Leg Length	mm
o) Actual Seat Width	mm



Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	SEDEO PRO
Typical WC Weight	Without Batteries 140 kg	Seat Size Range	Width 380 - 560 mm, Depth 400 - 560 mm
	With Batteries 165 kg	Seat To Floor Height Range	415 mm - 490 mm

MWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" Drive Wheels (Black)
Base Size	Width Min. 620 mm	Tires / Castors	7" Front / 7" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	50, 60 And 80 Ah Batteries
Drive Motors	2-Pole 6 kph, 4-Pole 6/8/10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

FWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 610 mm	Tires / Castors	9" Rear (Black), Option 10"
Drive Control	VR2, Option R-net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	50, 60 And 80 Ah Batteries
Drive Motors	2-Pole 6 kph, 4-Pole 6/8/10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	65 mm	Frame Colour	Black

RWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 610 mm	Tires / Castors	9" Front (Black), Option 10"
Drive Control	VR2, Option R-net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	50, 60 And 80 Ah Batteries
Drive Motors	2-Pole 6 kph, 4-Pole 6/8/10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	65 mm	Frame Colour	Black

Base

Q500 Drive Bases

		€
QOP000023	Q500 M (MWD Base Chair)	–
QOP000024	Q500 F (FWD Base Chair)	–
QOP000029	Q500 R (RWD Base Chair)	–

Suspension Setting (MWD Only)

		€
QOP150001	Standard Suspension	–
QOP150002	Firm Suspension Recommended with QOP130009 Seat elevator/tilt module, Recommended with Powered tilt (QOP130004/QOP130012/ QOP130007 for clients over 75 kg user weight	XXX,-

Motor Packages

		€
MWD Base		
QOP180011	2-Pole Motor, 6 kph / 4 mph, Standard	XXX,-
MWD/FWD/RWD		
QOP180012	2-Pole Motor, 6 kph / 4 mph, Performance Only With Drive Wheel 14"	XXX,-
QOP180009	4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed)	XXX,-
QOP180010	4-Pole Motor, 8 kph / 5 mph	XXX,-
QOP180007	4-Pole Motor, 10 kph / 6.5 mph If FWD Selected: Highly Recommended With Traction Control System	XXX,-
QOP180008	4-Pole Motor, 12.5 kph / 8 mph If FWD Selected: Only In Combination With Traction Control System	XXX,-

Tracking

		€
QOP110009	Traction Control System For R-Net (GYRO Technology)	XXX,-
QOP110014	Omit Traction Control System For R-Net	XXX,-

Drive Wheels 

		€
MWD/FWD/RWD		
QOP070065	Drive Wheel, 13", Pneumatic, Knobbly, Black	XXX,-
QOP070069	Drive Wheel, 13", Solid, Knobbly, Black	XXX,-
QOP070034	Drive wheel, 14", Pneumatic, Knobbly, Black	XXX,-
QOP070033	Drive wheel, 14", Solid, Knobbly, Black	XXX,-
FWD/RWD Base		
QOP070032	Drive Wheel, 14", Pneumatic, V-Type, Black	XXX,-
QOP070031	Drive Wheel, 14", Solid, V-Type, Black	XXX,-
QOP070091	Drive Wheel, 14", Pneumatic, Trelleborg, Grey Only With QOP080065 Castor Wheel, 9", Solid, Grey; Only with QOP180007 4-Pole Motor 10 kph; Only for Norway, Only For NL Supplied Products	XXX,-

Castor Wheels 

		€
MWD Base		
QOP080213	Castor Wheel, 7", Solid, Black	XXX,-
FWD		
QOP080210	Castor Wheel, 9", Solid, Black	XXX,-
QOP080211	Castor Wheel, 9", Pneumatic, Black	XXX,-
QOP080065	Castor Wheel, 9", Solid, Grey Only with QOP070091 Drive Wheel 14" Pneumatic Trelleborg Grey; Only for Norway, Only For NL Supplied Products	
RWD Base		
QOP080210	Castor Wheel, 9", Solid, Black	XXX,-
QOP080211	Castor Wheel, 9", Pneumatic, Black	XXX,-
QOP080058	Castor Wheel, 10", Solid, Black Only With Drive Wheel 14"	XXX,-
QOP080212	Castor Wheel, 10", Pneumatic, Black Only With Drive Wheel 14"	XXX,-

Mudguards (FWD/RWD Base Only)  €

QOP010007	Mudguard, Drive Wheel	XXX,-
QOP010008	Mudguard, Castor Wheel, 9"	XXX,-

Kerb Climber (RWD Only)  €

QOP090038	Kerb Climber	XXX,-
-----------	--------------	-------

Colour Accents €

QOP140005	Apple Green (RAL 110 70 77)	XXX,-
QOP140004	Black (RAL 9005)	XXX,-
QOP140003	Blue (RAL 260 40 45)	XXX,-
QOP140006	Orange (RAL 060 60 80)	XXX,-
QOP140002	Red (RAL 040 40 67)	XXX,-
QOP140007	Silver (RAL 9006)	XXX,-
QOP140001	White (RAL 9010)	XXX,-

LED Lights & Indicators  €

QOP110201	LED Light And Indicator QOP110030 Lighting Module VR2 Or QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	XXX,-
QOP110202	LED Fender Light QOP110030 Lighting Module VR2 or QOP110032 Actuator Control & Lights Module required – Not Road Approved And Without Indicators; On Q500 M Only With QOP040123 Armrest, Reclining, Flip Back or QOP040070 Armrest, Reclining, Multi Adjustable, Flip Back or QOP040073 Armrest, Reclining, Multi Adjustable, Flip Back, With Lock Function; Not With QMP130012 Posterior Tilt, Powered, -5° - 25°; On Q500 F and R Only With QOP010007 Mudguard, Drive Wheel; Not With SH415	XXX,-
QOP110209	Omit Light And Indicator	XXX,-

Seat Heights Measured From Front Of Seat Pan To Floor. Obtain Quote To Verify Seat To Floor Height €

SH415	415 mm Only With 60 Ah Batteries, Not With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	-
SH425	425 mm Only With 60 Ah Batteries, Only With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	-
SH44	440 mm Not With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	-
SH45	450 mm Only With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	-
SH465	465 mm	-
SH49	490 mm	-

Batteries & Chargers  €**Batteries**

QOP160003	Battery, 50 Ah, Gel, Maintenance-Free	XXX,-
QOP160001	Battery, 60 Ah, Gel, Maintenance-Free	XXX,-
QOP160005	Battery, 80 Ah, Gel, Maintenance-Free	XXX,-
QOP160033	Battery, 60 Ah (C20), Gel, Maintenance-Free, Set, Type Exide Only for NL Helmond supplied products	XXX,-
QOP160034	Battery, 70 Ah (C20), Gel, Maintenance-Free, Set, Type Exide Only for NL Helmond supplied products	XXX,-
QOP160009	Omit Batteries	XXX,-

Chargers

QOP160101	Charger, Automatic, 220 V / 8 A	XXX,-
QOP160100	Charger, Automatic, 220 V / 8 A, UK	XXX,-
QOP160102	Charger, Classic, 220 V / 8 A, IP54	XXX,-
QOP160103	Charger, Automatic, 220 V / 10 A, UK	XXX,-
QOP160104	Charger, Automatic, 220 V / 10 A	XXX,-
QOP160029	Charger, Automatic, Vitron Impulse25 smart (IXP4) 220 V / 8 A Only for Norway	XXX,-
QOP160109	Omit Charger	XXX,-
QOP160035	Charger, Automatic, IVRA, 220V / 8A Only for NL Helmond supplied products	XXX,-

Powered Seating

Power Seating Options

€

QOP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	–
QOP130004	Posterior Tilt, Powered, 0° - 30° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	–
QOP130012	Posterior Tilt, Powered, -5° - 25° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required, Not With Seat Height 415 Or 440 mm	–
QOP130007	Posterior Tilt, Powered, 0° - 50° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	–
QOP130005	Posterior Tilt, Powered, -5° - 45° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required, Not With Seat Height 415 Or 440 mm	–
QOP130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	–
QOP130006	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered Not With QOP110036 Remote Joystick, VR2; In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	–

C-Me Function C-Me Function Require QMP130009 / QMP130006. Not available with QMP030084 /QMP030085 Backrest, Manual Recline

€

QOP110034	C-Me Function, Enabled 4.5 kph Maximum Speed When Lifted (Tilt Angle <9°), Not Available With 50° Lift/Tilt In Combination With 90°-150° Power Recline	–
QOP110035	C-Me Function, Disabled 1.5 kph Maximum Speed When Lifted	XXX,-

Backrest Frame Type

€

QOP030083	Backrest, Tool Angle Adjustable, 90° - 120°	–
QOP030084	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right Not On Same Joystick Side	XXX,-
QOP030085	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left Not On Same Joystick Side	XXX,-
QOP030086	Backrest, Powered Recline, 90° - 130° Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	XXX,-
QOP030087	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	XXX,-
QOP030229	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear with Relax function Not With QOP110036 Remote Joystick, VR2. In Case Of QOP110004 Power Module, R-Net, 120 A: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required; Only With QMP130006 Seat Elevator (300mm) Tilt Module (Tilt 0° - 50°), Powered	XXX,-

Backrest Frame Option

€

QOP030098	Backrest, Flip Down, Quick Release	XXX,-
-----------	------------------------------------	-------

Back Cushion Type

€

SEDEO PRO (See Section A)

QOP030104	SEDEO PRO, Back Cushion	–
-----------	-------------------------	---

JAY J3 (See Section B)

QOP030054	JAY J3, Back Cushion	XXX,-
-----------	----------------------	-------

Interfaces

QOP030106	SEDEO PRO, Interface Only	–
QOP030107	SEDEO PRO, Interface Only For Back Cushion, Flexible	–
QOP030055	JAY J3, Interface Only	–
QOP030209	Omit Back Cushion	–

A – SEDEO PRO Back Size

€

Back Cushion Width

BW38	380 mm, For Effective Back Width 380 mm	–
BW42	420 mm, For Effective Back Width 380 - 420 mm	–
BW46	460 mm, For Effective Back Width 420 - 460 mm	–
BW50	500 mm, For Effective Back Width 460 - 500 mm	–
BW54	540 mm, For Effective Back Width 500 - 540 mm	–
BW56	560 mm, For Effective Back Width 520 - 560 mm	–

Back Cushion Height

BH44	440 mm, For Effective Back Height 440 - 500 mm Only With SEP030116 SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps & QOP030117 SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adj Straps	–
BH50	500 mm, For Effective Back Height 500 - 560 mm	–
BH56	560 mm, For Effective Back Height 560 - 620 mm	–

A – SEDEO PRO Back Cushion

€

QOP030109	SEDEO PRO, Back Cushion, Profiled	XXX,-
QOP030113	SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	XXX,-
QOP030114	SEDEO PRO, Back Cushion, Deep Contour, Passive	XXX,-
QOP030115	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Passive, Adjustable Straps Only Available For NL	XXX,-
QOP030116	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	XXX,-
QOP030117	SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	XXX,-
	SEDEO PRO, Back Cushion, Adaptive (Choose B4Me Requirements Separately)	XXX,-

A – SEDEO PRO Back Cushion Cover

€

QOP030059	Back Cushion Cover, Incontinence Protection	–
QOP030108	Back Cushion Cover, Comfort	XXX,-
QOP030230	Back Cushion Cover, Comfort, Double zipper only for QOP030116 Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps and QOP030117 Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps. Only For Norway, Only For NL Supplied Products	XXX,-

A – SEDEO Lateral Thoracic Supports Only With QOP030109 SEDEO PRO, Back Cushion, Profiled & QOP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps

€

QOP030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	XXX,-
QOP030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	XXX,-

B – JAY J3 Back Size

€

JAY030111	JAY J3 Back Size, 410 mm	XXX,-
JAY030112	JAY J3 Back Size, 460 mm	XXX,-
JAY030113	JAY J3 Back Size, 510 mm	XXX,-

B – JAY J3 Back, Mid Deep Contour

€

JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	XXX,-
-----------	---	-------

B – JAY J3 Back, Deep Contour

€

JAY030021	JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm	XXX,-
-----------	---	-------

Headrest Pad

€

QOP030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm Only Use With QOP030433 or QMP030434 SEDEO Headrest Mount	XXX,-
QOP030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm Only Use With QOP030433 or QMP030434 SEDEO Headrest Mount	XXX,-
QOP030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm Only Use With QOP030433 or QMP030434 SEDEO Headrest Mount	XXX,-
WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm Only Use With WHT030408 / WHT030409 / WHT030411 Whitmyer Headrest Mount	XXX,-
WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm Only Use With WHT030408 / WHT030409 / WHT030411 Whitmyer Headrest Mount	XXX,-

Headrest Mount Choose A Whitmyer Headrest Mount And Headrest Receiver In Case You Would Like To Add A Head Control In The Special Control Section

QOP030433	SEDEO, Headrest Mount	XXX,-
QOP030434	SEDEO, Headrest Mount, With Quick Release	XXX,-
WHT030408	Whitmyer, Headrest Mount, LINX ²	XXX,-
WHT030409	Whitmyer, Headrest Mount, AXYS	XXX,-
WHT030411	Whitmyer, Headrest Mount, COBRA XTRA Flip Back	XXX,-

Headrest Receiver

QOP030119	Headrest Receiver, Back Cushion Mount Only With QOP030109 SEDEO PRO, Back Cushion, Profiled & QOP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	XXX,-
QOP030120	Headrest Receiver, Backrest Frame Mount	XXX,-

Push Handles

QOP030204	Push Handles, Height-Adjustable	XXX,-
-----------	---------------------------------	-------

Seat Rail

QOP020085	Seat Rail, Short, For Seat Depth 400 - 560 mm	-
QOP020069	Seat Rail, Long, For Seat Depth 460 - 560 mm Not With Seat Width 380 - 420 mm	XXX,-

Seat Cushion Type

SEDEO (See Section C)

QOP020027	SEDEO, Seat Cushion	-
-----------	---------------------	---

JAY (See Section D)

QOP020515	JAY, Seat Cushion	XXX,-
-----------	-------------------	-------

Omit (Select Seat Width And Depth in Section C)

QOP020110	Omit Seat Cushion	XXX,-
-----------	-------------------	-------

C – SEDEO Seat Size

Seat Cushion Width

SW38	380 mm, For Effective Seat Width 380 mm	-
SW42	420 mm, For Effective Seat Width 380 - 420 mm	-
SW46	460 mm, For Effective Seat Width 420 - 460 mm	-
SW50	500 mm, For Effective Seat Width 460 - 500 mm	-
SW54	540 mm, For Effective Seat Width 500 - 540 mm	-
SW56	560 mm, For Effective Seat Width 520 - 560 mm	-

Seat Cushion Depth

SD40	400 mm, Comes With Cushion Size 460 mm Not with QOP090128 or QOP090129 Positioning belt Short or Long	-
SD42	420 mm, Comes With Cushion Size 460 mm	-
SD44	440 mm, Comes With Cushion Size 460 mm	-
SD46	460 mm, Comes With Cushion Size 500 mm	-
SD48	480 mm, Comes With Cushion Size 500 mm	-
SD50	500 mm, Comes With Cushion Size 500 mm	-
SD52	520 mm, Comes With Cushion Size 560 mm	-
SD54	540 mm, Comes With Cushion Size 560 mm	-
SD56	560 mm, Comes With Cushion Size 560 mm	-

C – SEDEO Seat Cushions  €

QOP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	XXX
QOP020074	SEDEO, Seat Cushion, Memory, Height 70 mm Only With QOP020518 Seat Cushion Cover, Incontinence Protection, Only For NL	XXX
QOP020075	SEDEO, Seat Cushion, Support, Height 90 mm	XXX
QOP020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm	XXX
QOP020077	SEDEO, Seat Cushion, Form, Height 110 mm Only With QOP020518 Seat Cushion Cover, Incontinence Protection, Only For NL	XXX
QOP020078	SEDEO, Seat Cushion, Stabilizer, Height 90 mm	XXX
	SEDEO, Seat Cushion, Adaptive (Choose B4Me Requirements Separately)	XXX

C – SEDEO Seat Cushion Cover €

QOP020518	Seat Cushion Cover, Incontinence Protection	XXX
QOP020080	Seat Cushion Cover, Comfort	XXX

C – SEDEO Seat Boards €

QOP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	XXX
QOP020083	SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear	XXX
QOP020084	SEDEO, Seat Pan, 75 mm, Depth-Adjustable, Cushions Guide At Rear Only For NL	XXX

D – JAY Seat Size €**Seat Cushion Width**

SW40	400 mm	–
SW42	420 mm	–
SW44	440 mm	–
SW46	460 mm	–
SW48	480 mm	–
SW50	500 mm, Comes With Cushion Size 500 mm	–
SW52	520 mm, Comes With Cushion Size 500 mm	–
SW54	540 mm, Comes With Cushion Size 560 mm	–
SW56	560 mm, Comes With Cushion Size 560 mm	–

Seat Cushion Depth

SD40	400 mm	–
SD42	420 mm	–
SD44	440 mm	–
SD46	460 mm	–
SD48	480 mm	–
SD50	500 mm	–
SD52	520 mm	–
SD54	540 mm Only With Seat Width 560 mm	–
SD56	560 mm Only With Seat Width 560 mm	–

D – JAY Seat Cushions 

€

JAY020052	JAY LowZone Gel Switch with Wipedown Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020006	JAY Easy Visco With Flat Base	XXX,-
JAY020053	JAY LowZone+ Flat with Wipedown Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020063	JAY LowZone+ Flat with Microclimate Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020054	JAY LowZone+ Contour with Wipedown Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020064	JAY LowZone+ Contour with Microclimate Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020055	JAY ClimaZone Flat with Wipedown Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020065	JAY ClimaZone Flat with Coolover Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020056	JAY ClimaZone Contour with Wipedown Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020066	JAY ClimaZone Contour with Coolover Cover, Height 80 mm UK & Ireland only; Only for Q500 M Base	XXX,-
JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	XXX,-
JAY020013	JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm	XXX,-
JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	XXX,-
JAY020017	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	XXX,-

Knee & Thigh Supports

€

QOP020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger Receivers Only, QOP090122 Transfer Bar	XXX,-
QOP020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger Receivers Only, QOP090122 Transfer Bar	XXX,-
QOP020113	Thigh Support, Lateral, Small, 120 x 120 mm, Removable, Right Not With QOP040035 Side Guard, Wedge Shaped, Padded	XXX,-
QOP020114	Thigh Support, Lateral, Small, 120 x 120 mm, Removable, Left Not With QOP040035 Side Guard, Wedge Shaped, Padded	XXX,-
QOP020121	Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Right Not With QOP040035 Side Guard, Wedge Shaped, Padded	XXX,-
QOP020122	Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Left Not With QOP040035 Side Guard, Wedge Shaped, Padded	XXX,-

Armrest Types 

€

QOP040019	Armrest, Single Post, Flip Back In case of MWD base, not with Seat height 415/seat height 425 & Tilt, Powered, 0° - 30°/ Tilt, Powered, 0° - 50°	-
QOP040021	Armrest, Single Post, Removable In case of MWD base, not with Seat height 415/seat height 425 & Tilt, Powered, 0° - 30°/ Tilt, Powered, 0° - 50°	XXX,-
QOP040022	Armrest, Single Post, SEDEO LITE Arm Pad, Integral Foam, 380 mm Right/Left Included; Side Guard Included; Only With QOP110082 Joystick Mounting, Parallel, Swing Away; Not With QOP090040 Tray Table, Width-Adjustable, Swing Away & QOP090046 Tray Table, Half-Side, Swing Away; Not With Special Controls Mounted On Armrest; Not With QOP090120 Bag For Armrest	XXX,-
QOP040123	Armrest, Reclining, Flip Back	XXX,-
QOP040070	Armrest, Reclining, Multi Adjustable, Flip Back	XXX,-
QOP040073	Armrest, Reclining, Multi Adjustable, Flip Back, With Lock Function Only For Norway, Only For NL Supplied Products	XXX,-

Arm Pads 

€

QOP040023	Arm Pad, Integral Foam, 300 mm, Right	XXX,-
QOP040024	Arm Pad, Integral Foam, 300 mm, Left	XXX,-
QOP040025	Arm Pad, Integral Foam, 400 mm, Right	-
QOP040026	Arm Pad, Integral Foam, 400 mm, Left	-
QOP040027	Arm Pad, Integral Foam, 450 mm, Right With QOP040019 Flip Back Armrest, Seat Width +20 mm	XXX,-
QOP040028	Arm Pad, Integral Foam, 450 mm, Left With QOP040019 Flip Back Armrest, Seat Width +20 mm	XXX,-

Comfort Arm Pads		
QOP040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	XXX,-
QOP040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	XXX,-
QOP040074	Arm Pad, Comfort, Padded, Stretch Rubber, 300 mm, Right Only For Norway, Only For NL Supplied Products	XXX,-
QOP040075	Arm Pad, Comfort, Padded, Stretch Rubber, 300 mm, Left Only For Norway, Only For NL Supplied Products	XXX,-
QOP040076	Arm Pad, Comfort, Padded, Stretch Rubber, 400 mm, Right Only For Norway, Only For NL Supplied Products	XXX,-
QOP040077	Arm Pad, Comfort, Padded, Stretch Rubber, 400 mm, Left Only For Norway, Only For NL Supplied Products	XXX,-
QOP040078	Arm Pad, Comfort, Padded, Removable Cover, 300 mm, Right Only For Norway, Only For NL Supplied Products	XXX,-
QOP040079	Arm Pad, Comfort, Padded, Removable Cover, 300 mm, Left Only For Norway, Only For NL Supplied Products	XXX,-
QOP040080	Arm Pad, Comfort, Padded, Removable Cover, 400 mm, Right Only For Norway, Only For NL Supplied Products	XXX,-
QOP040081	Arm Pad, Comfort, Padded, Removable Cover, 400 mm, Left Only For Norway, Only For NL Supplied Products	XXX,-

Hemiplegia Arm Pads		
QOP040031	Arm Pad, Hemiplegia, Integral Foam, Right Not With QOP110052 Joystick Mount, Right, QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away, QOP090110 Rear View Mirror, Right; not with QOP110362 Hand Control Mounting Location Right; not with QOP110119/QOP110123 OMNI2 Module Mount Right	XXX,-
QOP040032	Arm Pad, Hemiplegia, Integral Foam, Left Not With QOP110053 Joystick Mount, Right, QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away, QOP090118 Rear View Mirror, Left; not with QOP110361 Hand Control Mounting Location Left; not with QOP110120/QOP110124 OMNI2 Module Mount Left	XXX,-

Side Guard Pads For Armrest €

Side Guard Pads For Armrest, Single Post, Flip Back / Removable		
QOP040033	Side Guard Not With QOP040123 Armrest, Reclining, Flip Back	XXX,-
QOP040034	Side Guard, Integral Foam Seat Width -20 mm, Not With QOP040123 Armrest, Reclining, Flip Back	XXX,-
QOP040035	Side Guard, Wedge Shaped, Padded Seat Width -80 mm, Only With QOP040021 Armrest, Single Post, Removable	XXX,-

Side Guard Pads For Armrest, Reclining, Flip Back		
QOP040036	Side Guard, Mounted On Seat Rail Bracket Only With QOP040123 Armrest, Reclining, Flip Back	XXX,-

Centre Mount Legrest €

QOP050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard In Combination With Lift Or Tilt, Choose QOP050087 Centre Calf Support, Integral Foam Or QOP050092 Calf Pad, Integral Foam, Pair; Limited Knee Angle Possibilities With RWD/HWD	XXX,-
-----------	---	-------

Pad Options		
QOP050087	Centre Calf Support, Integral Foam	XXX,-
QOP050092	Calf Pad, Integral Foam, Pair	XXX,-

Fixed Centre Mount Legrest €

QOP050140	Fixed Centre Mount Legrest Not with RWD; In Combination With Lift Or Tilt Choose QPP050092 Calf Pad, Intergral Foam, Pair; Only For UK Supplied Products	XXX,-
-----------	---	-------

Footplates		
QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	XXX,-

Pad Options		
QOP050092	Calf Pad, Integral Foam, Pair	XXX,-

Power Centre Mount Legrest 		€
QOP050076	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension In Combination With Lift Or Tilt, Choose QOP050092 Calf Pad, Integral Foam, Pair; Not with RWD	XXX,-
QOP050071	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor Not With VR2: QOP110036, QOP110037, QOP110040; In Case Of R-Net: QOP110032 required; In Combination With Lift Or Tilt, Choose QOP050092 Calf Pad, Integral Foam, Pair; Not with RWD	XXX,-
QOP050182	Centre Mount Legrest, Powered, 90° - 180°, 135 mm Extension, Short, LLL 300mm In Combination With Lift Or Tilt, Choose QOP050092 Calf Pad, Integral Foam, Pair; Only For Norway, Only For NL Supplied Products; Not with RWD	
Footplates		
QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	XXX,-
QOP050073	Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	XXX,-
Pad Options		
QOP050092	Calf Pad, Integral Foam, Pair	XXX,-
SEDEO PRO Swing Away Legrest 		€
QOP050015	SEDEO PRO, Legrest, 70°, Swing Away	XXX,-
QOP050016	SEDEO PRO, Legrest, Manual Elevating, Swing Away	XXX,-
QOP050017	SEDEO PRO, Legrest, Power Elevating, Swing Away	XXX,-
QOP050183	SEDEO PRO, Legrest, Power Elevating, Swing Away, Short, LLL 300 mm	XXX,-
QOP050000	Omit Legrest	XXX,-
Hanger Receivers		
QOP050170	SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thickness 70 mm And Higher	XXX,-
QOP050171	SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thickness 50 mm And Lower	XXX,-
QOP050080	SEDEO PRO, Hanger Receivers Only, High	XXX,-
QOP050172	SEDEO PRO, Hanger Receivers Only, Low	XXX,-
Footplates		
QOP050088	Footplate, For Seat Width 385 - 400 mm	XXX,-
QOP050089	Footplate, For Seat Width 420 - 460 mm	XXX,-
QOP050090	Footplate, For Seat Width 480 - 520 mm	XXX,-
QOP050091	Footplate, For Seat Width 540 - 560 mm	XXX,-
QOP050081	Footplates, Coupled Choose Also A Footplate; Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away	XXX,-
Support Options 		
QOP050082	Calf Support, Integral Foam, Depth-Adjustable, Pair Not On RWD With QOP050015 SEDEO PRO, Legrest, 70°, Swing Away	XXX,-
QOP050083	Calf Support, Padded, Depth-Adjustable, Pair Not On RWD With QOP050015 SEDEO PRO, Legrest, 70°, Swing Away	XXX,-
QOP050084	Heel Loop	XXX,-
QOP050085	Calf Strap Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away	-
QOP050200	Knee Pad, Integral Foam	XXX,-
Amputee Support		
QOP050142	Amputee Support, Integral Foam, Right	XXX,-
QOP050143	Amputee Support, Integral Foam, Left	XXX,-
QOP050013	Amputee Support, Padded, Right Extra On A Set Of Swing Away Legrests, Attached Lose	XXX,-
QOP050014	Amputee Support, Padded, Left Extra On A Set Of Swing Away Legrests, Attached Lose	XXX,-

Positioning Supports & Accessories		€
QOP090128	Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm	–
QOP090129	Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm	XXX,-
QOP090110	Rear View Mirror, Right Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away	XXX,-
QOP090118	Rear View Mirror, Left Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away	XXX,-
QOP090009	Crutch Holder Not With QOP030086 Powered Recline, 90° - 130°, QOP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear	XXX,-
QOP090119	Brackets For Bag On Backrest Frame Not With QOP030086 Powered Recline, 90° - 130°, QOP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear	XXX,-
QOP090120	Bag For Armrest Not With QOP040123 Armrest, Reclining, Flip Back	XXX,-
QOP090122	Transfer Bar Not With QOP020342/QMP020341 Knee Support, Lateral, Swing Away, Right/Left & QOP050015-17 SEDEO PRO, Legrest, Swing Away & QOP050080 SEDEO PRO, Hanger Receivers Only	XXX,-
QOP030605	Abduction Pommel Not With (Power) Centre Mount Legrests	XXX,-
QOP090205	Dahl Lockplate Mounted Only For Norway, Only For NL Supplied Products	
Tray Table (Mounted Opposite Joystick)		
QOP090040	Tray Table, Width-Adjustable, Swing Away	XXX,-
QOP090046	Tray Table, Half-Side, Swing Away Not with QOP090126 raised Edges, Acrylic For Tray Table	XXX,-
Tray Table Accessories		
QOP090126	Raised Edges, Acrylic For Tray Table	XXX,-

Joystick / Electronics

Remote Joystick – VR2 		€
QOP110036	Remote Joystick, VR2	–
QOP110037	Remote Joystick, VR2, 2 Actuators	XXX,-
QOP110040	Remote Joystick, VR2, 2 Actuators And Lights	XXX,-

Joystick Mounting 		€
QOP110050	Joystick Mounting, Fixed	–
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	XXX,-
QOP110082	Joystick Mounting, Parallel, Swing Away Only With QOP040022 Armrest, Single Post, SEDEO LITE	XXX,-

Joystick Mounting Location		€
QOP110052	Joystick Mounting Location, Right	–
QOP110053	Joystick Mounting Location, Left	–

Joystick Knobs (Optional) 		€
QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	XXX,-
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	XXX,-
QOP090508	Joystick Knob, Ball, Ø 35 mm	XXX,-
QOP090045	Joystick Knob, Allround, Ball, Ø 48 mm	XXX,-
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	XXX,-
QOP090512	Joystick Knob, Stick, 80 x 15 mm	XXX,-
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	XXX,-

Electronic Programmers – VR2 €

QOP090065	Mobility Programmer With Dongle, VR2	XXX,-
-----------	--------------------------------------	-------

Power Module – VR2 €

QOP110002	Power Module, VR2, 90 A Power Legrests Not Available In Combination With C-Me	XXX,-
-----------	--	-------

Lighting Module – VR2 €

QOP110030	Lighting Module, VR2 If LED Light And Indicator Is Chosen, Then This Is Required	XXX,-
-----------	---	-------

Attendant Controls €

QOP110021	Attendant Control, Dual	XXX,-
QOP110020	Joystick Mounted As Attendant Control Choose Also Remote Joystick	XXX,-

Attendant Mounting Location

QOP110025	Attendant Mounting Location, Right	–
QOP110026	Attendant Mounting Location, Left	–

Remote Joystick – R-Net €

QOP110015	Remote Joystick, R-Net, LED, Actuators	–
QOP110016	Remote Joystick, R-Net, LED, Actuators And Lights	–
QOP110017	Remote Joystick, R-Net, Colour LCD	–
QOP110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	XXX,-
QOP110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	XXX,-
QOP110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	XXX,-
QOP110019	Omit Remote Joystick, R-Net	XXX,-

Joystick Mounting €

QOP110050	Joystick Mounting, Fixed	–
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	XXX,-
QOP110082	Joystick Mounting, Parallel, Swing Away Only With QOL040022 Armrest, Single Post, SEDEO LITE	XXX,-

Joystick Mounting Location €

QOP110052	Joystick Mounting Location, Right	–
QOP110053	Joystick Mounting Location, Left	–

Display Language Selection €

QOP119000	German	QOP119004	Portuguese	QOP119009	Swedish	QOP119013	Slovak
QOP119001	English	QOP119005	Italian	QOP119010	Greek	QOP119014	Polish
QOP119002	French	QOP119006	Dutch	QOP119011	Finnish	QOP119015	Danish
QOP119003	Spanish	QOP119008	Norwegian	QOP119012	Czech		

Joystick Knobs (Optional) €

QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	XXX,-
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	XXX,-
QOP090508	Joystick Knob, Ball, Ø 35 mm	XXX,-
QOP090045	Joystick Knob, Allround, Ball Ø 48 mm	XXX,-
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	XXX,-
QOP090512	Joystick Knob, Stick, 80 x 15 mm	XXX,-
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	XXX,-

Electronic Programmers – R-Net €

QOP090066	PC Programmer With Dongle, R-Net	XXX,-
-----------	----------------------------------	-------

Power Modules – R-Net €

QOP110003	Power Module, R-Net, 90 A	XXX,-
QOP110004	Power Module, R-Net, 120 A	XXX,-

Seating / Lighting Module Control – R-Net €

QOP110032	Actuator Control & Lights Module, R-Net, ISM-6L	XXX,-
-----------	---	-------

Attendant Controls €

QOP110021	Attendant Control, Dual	XXX,-
QOP110020	Joystick Mounted As Attendant Control Choose Also Remote Joystick	XXX,-

Attendant Mounting Location

QOP110025	Attendant Mounting Location, Right	–
QOP110026	Attendant Mounting Location, Left	–
QOP110026	Attendant Mounting Location, Left	–

Table Top Controls €

QOP110143	Table Top Control, LCD, Centre Tray Mounted	XXX,-
-----------	---	-------

Tray Table Mounting Location

QOP110080	Tray Table Mounting Location, Left	–
QOP110081	Tray Table Mounting Location, Right	–

#1 Ctrl+5 Switch Box Options €

QOP110043	#1 Ctrl+5 Switch Box With Free Assignable Functions If combined with attendant controll, Ctr/+5 will be mounted to the armrest	XXX,-
-----------	---	-------

#1 Ctrl+5 Switch Box Mounting Location

QOP110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick Not Possible With QOP110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)	–
QOP110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	–

#1 Ctrl+5 Switch Type

QOP110175	#1 Ctrl+5 Switch Type With Push Buttons	–
QOP110176	#1 Ctrl+5 Switch Type With Toggle Switches	–

Direct Access Only For Norway, Only For NL Supplied Products

QOP110174	Direct Access, Behind Joystick	XXX,-
QOP110188	Direct Access, Opposite Joystick	XXX,-

Special Controls**OMNI2 Input Control Module – R-Net** Choose OMNI2 Mount And Input Device in Sections Below And Choose Switches in Section E €

QOP110116	OMNI2 Input Control Module, R-Net, Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker	XXX,-
QOP110117	OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker	XXX,-
QOP110118	OMNI2 Input Control Module, R-Net, Scanning Mode, Traction Control System (GYRO Technology), Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker	XXX,-
QOP110008	Programming Service Not Discountable	XXX,-

OMNI2 Module Mount If Used With Standard Joystick Selection, It Must Be Mounted Opposite Side €

QOP110119	OMNI2 Module, Fixed Mount, Right	–
QOP110120	OMNI2 Module, Fixed Mount, Left	–
QOP110123	OMNI2 Module, Swing Away Mount, Right Not With QOP090110 Rear View Mirror Right	XXX,-
QOP110124	OMNI2 Module, Swing Away Mount, Left Not With QOP090118 Rear View Mirror Left	XXX,-

Hand Control For Switch Selection See Section E

€

Proportional Input Devices

QOP110301	Hand Control, Switch-It MicroGuide	XXX,-
QOP110302	Hand Control, Switch-It MicroPilot	XXX,-
QOP110304	Hand Control, Switch-It VersaGuide	XXX,-
QOP110305	Hand Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP110306	Hand Control, Allround Joystick	XXX,-
QOP110307	Hand Control, Allround Joystick Light Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm	XXX,-
QOP110308	Hand Control, Allround Joystick Heavy Duty	XXX,-
QOP110309	Hand Control, Allround Joystick, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	XXX,-
QOP110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	XXX,-
QOP110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	XXX,-

Proportional Input Devices, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker

QOP110044	Hand Control, Switch-It MicroGuide	XXX,-
QOP110045	Hand Control, Switch-It MicroPilot	XXX,-
QOP110047	Hand Control, Switch-It VersaGuide	XXX,-
QOP110048	Hand Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP110049	Hand Control, Allround Joystick	XXX,-
QOP110059	Hand Control, Allround Joystick Light Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm	XXX,-
QOP110060	Hand Control, Allround Joystick Heavy Duty	XXX,-

Hand Control Mounting Platforms

€

QOP110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray Only For Micro Series Joysticks	XXX,-
-----------	--	-------

Hand Control Mounting 

QOP110350	Hand Control Mounting, Link-It To The Armrest, Fixed	XXX,-
QOP110353	Hand Control Mounting, Link-It To The Armrest, Swing Away	XXX,-

Hand Control Mounting Location

QOP110362	Hand Control Mounting Location, Right	XXX,-
QOP110361	Hand Control Mounting Location, Left	XXX,-

Head Control For Switch Selection See Section E. Choose A Whitmyer Headrest Mount On Previous Pages 

€

Proportional Input Devices

QOP110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	XXX,-
QOP110371	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	XXX,-
QOP110208	Head Control, Switch-It Vigo Head Drive Only with headrest pad WHT030402 or WHT030403	XXX,-

Proportional Input Devices, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker

QOP110061	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	XXX,-
QOP110062	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	XXX,-
QOP110210	Head Control, Switch-It Vigo Head Drive Only with headrest pad WHT030402 or WHT030403	XXX,-

Chin Control 

€

Proportional Input Devices

QOP110470	Chin Control, Switch-It MicroGuide	XXX,-
QOP110471	Chin Control, Switch-It MicroPilot	XXX,-
QOP110472	Chin Control, Switch-It VersaGuide	XXX,-
QOP110473	Chin Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP110474	Chin Control, Allround Joystick	XXX,-
QOP110475	Chin Control, Allround Joystick Light Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP110476	Chin Control, Allround Joystick, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	XXX,-
QOP110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 Module Not Required, Obligated to Have Master R-Net Module; Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-

Proportional Input Devices, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker

QOP110070	Chin Control, Switch-It MicroGuide	XXX,-
QOP110071	Chin Control, Switch-It MicroPilot	XXX,-
QOP110072	Chin Control, Switch-It VersaGuide	XXX,-
QOP110073	Chin Control, Switch-It VersaGuide EZ Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP110074	Chin Control, Allround Joystick	XXX,-
QOP110075	Chin Control, Allround Joystick Light Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide	XXX,-

Chin Control Mounting 

QOP110480	Chin Control Mounting, Link-It To Chin Harness, Small If You Need Switches Mounted On The Harness, You Can Only Use Allround Joystick	XXX,-
QOP110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	XXX,-
QOP110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	XXX,-
QOP110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	XXX,-
QOP110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	XXX,-
QOP110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	XXX,-
QOP110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	XXX,-

Link-It Extension Kit

QOP110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	XXX,-
QOP110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	XXX,-

Sip 'N Puff For Switch Selection See Section E

€

QOP110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest Headrest To Be Selected	XXX,-
Sip 'N Puff, Price Inclusive OMNI2 Module And Traction Control System (GYRO Technology) OMNI2 Module Incl. Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker		
QOP110076	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest Headrest To Be Selected	XXX,-

E – Switches Select 1 Switch And Mounting Position Accordingly. For A Second Switch Please Use #2 Switch Selection Part.
 Single Function For Special Input Devices 

€

#1 Switch Selection

QOP110611	#1 Single Switch, Ø 64 mm, Black	XXX,-
QOP110612	#1 Single Switch, Ø 64 mm, Red	XXX,-
QOP110613	#1 Single Switch, Ø 64 mm, Green	XXX,-
QOP110523	#1 Egg Switch – Black	XXX,-
QOP110524	#1 MicroLite Switch – Black	XXX,-
QOP110651	#1 Single Switch, Ø 30 mm, Black	XXX,-
QOP110654	#1 Single Switch, Ø 30 mm, Blue	XXX,-
QOP110653	#1 Single Switch, Ø 30 mm, Green	XXX,-
QOP110652	#1 Single Switch, Ø 30 mm, Red	XXX,-
QOP110655	#1 Single Switch, Ø 30 mm, Yellow	XXX,-
QOP110534	#1 Twister Switch – Black	XXX,-
QOP110535	#1 Twister Switch – Green	XXX,-
QOP110536	#1 Twister Switch – Red	XXX,-

#1 Switch Mounting 

QOP110560	#1 Switch Mounting, Link-It To The Hand Control Not With HD Allround Joystick Or Space Disk Mount	XXX,-
QOP110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) Armrest Rail Mount	XXX,-
QOP110561	#1 Switch Mounting, Link-It To The Headrest	XXX,-
QOP110563	#1 Switch Mounting, Link-It To The Seat Rail Operated By Knee, If Two Switches Selected, One Mounted On The Right Side, One On The Left	XXX,-
QOP110565	#1 Switch Mounting, Link-It To Chin Boom No Left / Right Switch Selection Required	XXX,-
QOP110569	#1 Omit Switch Mounting Switches Attached Lose	XXX,-

#1 Switch Mounting Location

QOP110572	#1 Switch Mounting Location, Right	XXX,-
QOP110571	#1 Switch Mounting Location, Left	XXX,-

#1 Satellite Switches Only For Allround Joystick, Maximum 2 Selectable

QOP110537	#1 Satellite Switch – Black	XXX,-
QOP110538	#1 Satellite Switch – Green	XXX,-
QOP110539	#1 Satellite Switch – Red	XXX,-

#2 Switch Selection

QOP110631	#2 Single Switch, Ø 64 mm, Black	XXX,-
QOP110632	#2 Single Switch, Ø 64 mm, Red	XXX,-
QOP110633	#2 Single Switch, Ø 64 mm, Green	XXX,-
QOP110088	#2 Egg Switch – Black	XXX,-
QOP110089	#2 MicroLite Switch – Black	XXX,-
QOP110671	#2 Single Switch, Ø 30 mm, Black	XXX,-
QOP110674	#2 Single Switch, Ø 30 mm, Blue	XXX,-
QOP110673	#2 Single Switch, Ø 30 mm, Green	XXX,-
QOP110672	#2 Single Switch, Ø 30 mm, Red	XXX,-
QOP110675	#2 Single Switch, Ø 30 mm, Yellow	XXX,-
QOP110095	#2 Twister Switch – Black	XXX,-
QOP110096	#2 Twister Switch – Green	XXX,-
QOP110097	#2 Twister Switch – Red	XXX,-

#2 Switch Mounting 

QOP110098	#2 Switch Mounting, Link-It To The Hand Control Not With HD Allround Joystick Or Space Disk Mount	XXX,-
QOP110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) Armrest Rail Mount	XXX,-
QOP110100	#2 Switch Mounting, Link-It To The Headrest	XXX,-
QOP110104	#2 Switch Mounting, Link-It To The Seat Rail Knee, If 2 Switches Selected Only Left And Right Mount	XXX,-
QOP110105	#2 Switch Mounting, Link-It To Chin Boom No Left / Right Switch Selection Required	XXX,-
QOP110106	#2 Omit Switch Mounting Switches Attached Lose	XXX,-

#2 Switch Mounting Location

QOP110108	#2 Switch Mounting Location, Right	XXX,-
QOP110107	#2 Switch Mounting Location, Left	XXX,-

#2 Satellite Switches Only For Allround Joystick, Maximum 2 Selectable

QOP110111	#2 Satellite Switch – Black	XXX,-
QOP110112	#2 Satellite Switch – Green	XXX,-
QOP110113	#2 Satellite Switch – Red	XXX,-

Multi-Function Switches Used In Combination With Assignable Button Programming

QOP110550	Mini 2 Button (Stereo Plug) Attached Lose, When Ordered. No Link-It Fixation Possible	XXX,-
QOP110551	Mini 2 Button (Dual Mono Plug) Attached Lose, When Ordered. No Link-It Fixation Possible	XXX,-
QOP110552	Mini 4 Button InLine (DB9 Connector) Attached Lose, When Ordered. No Link-It Fixation Possible	XXX,-
QOP110554	4-Way Toggle Switched Joystick	XXX,-

Accessories

QOP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	XXX,-
QOP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	XXX,-
QOP110582	Extension Cable DB9 (1 m)	XXX,-
QOP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	XXX,-

Safety Devices

€

QOP090147	Remote Stop Switch (9 m Activation Range, Fit Via OBC Port) Not With VR2	XXX,-
-----------	---	-------

Drive Wheels & Castor Wheels & Mudguards & Kerb Climber

◀ BACK



// QOP070069



// QOP080213



// QOP010007



// QOP010008



// QOP090038

LED Lights & Chargers

◀ BACK



// QOP110209



// QOP110209



// QOP160101



// QOP160104

SEDEO PRO Back Cushions

◀ BACK



// QOP030109



// QOP030113



// QOP030114



// QOP030115



// QOP030116



// QOP030117

JAY J3 Back

◀ **BACK**



// JAY030012



// JAY030013



// JAY030045



// JAY030014

Seat Width in mm	Shallow Contour	Mid Contour	Mid Deep Contour / Extra Deep	Deep Contour / Extra Deep
	JAY030012	JAY030013	JAY030045/JAY030050	JAY030014/JAY030051
Maximum User Trunk Width in mm				
360	330	295	280	280
410	380	345	330	330
460	430	395	380	380
510	480	445	430	430

JAY J3 Back, Spine Align Fit Kit & Positioning Boxes

◀ **BACK**



// JAY090060-JAY090034



// JAY090059-JAY090031



// JAY090035



// JAY090036

JAY J3 Back, Chest Strap

◀ **BACK**



// JAY090040 - JAY090042

Headrest Pads & Headrest Mounts & Push Handles

◀ **BACK**




// QOP030430 & QOP030431



// WHT030401-WHT030403



// WHT030404-WHT030406



// QOP030433



// WHT030407



// WHT030408



// WHT030409



// WHT030410



// WHT030411



// QOP030204

SEDEO Seat Cushions

◀ **BACK**



// QOP020073



// QOP020074



// QOP020075



// QOP020076



// QOP020077



// QOP020078

JAY Seat Cushions

◀ **BACK**



// JAY020005



// JAY020003



// JAY020006



// JAY020010



// JAY020011



// JAY020018



// JAY020019



// JAY020020

Armrest Types & Arm Pads

◀ **BACK**



// QOP040019



// QOP040021



// QOP040123




// QOP040027



// QOP040023



// QOP040204



// QOP040070

Fixed Centre Mount Legrest & Centre Mount Legrest

◀ **BACK**



// QOP050086



// QOP050140

SEDEO PRO Swing Away Legrest & Support Options & Power Centre Mount Legrest

◀ BACK



// QOP050016



// QOP050076

Remote Joystick– R-Net & Joystick Mounting

◀ BACK



// QOP110017



// QOP110013



// QOP110018



// QOP110056

Remote Joystick – VR2

◀ BACK



// QOP110036



// QOP110037

Joystick Knobs

◀ BACK



// QOP090505



// QOP090045



// QOP090507



// QOP090508



// QOP090509



// QOP090512



// QOP090511



// QOP090506

Electronic Programmers – R-Net

◀ **BACK**



// QOP090066

Power Modules – R-Net & Seating / Lighting Module Control – R-Net

◀ **BACK**



// QOP110004

Attendant Controls

◀ **BACK**



// QOP110021

#1 Ctrl+5- Switch Box Options & OMNI2 Input Control Module – R-Net

◀ **BACK**



// QOP110175



// QOP110176



// QOP090068



// QOP110116

Head Control

◀ BACK



// QOP110401 - QOP110406



// QOP110370 // QOP110371

Chin Control & Hand Control

◀ BACK



// QOP110470 // QOP110301



// QOP110471 // QOP110302



// QOP110472 // QOP110473 // QOP110304 // QOP110305



// QOP110306 - QOP110309 // QOP110474 - QOP110477



// QOP110475 // QOP110477 // QOP110307 // QOP110310

Table Top Controls

◀ BACK



// QOP110143

Switches

◀ BACK



// QOP110651 - QOP110655
// QOP110671 - QOP110675



// QOP110523



// QOP110524



// QOP110611 - QOP110613
// QOP110631 - QOP110633



// QOP110537 - QPP110539



// QOP110534 - QOP110536



// QOP110551



// QOP110552

Switch Mounting

◀ BACK



// QOP110526 // QOP110566



// QOP110566 // QOP110526
// QOP110529



// QOP110560 // QOP110572
// QOP110013 // QOP110056
// QOP110529 // QOP110526



// QOP110563

Hand Control Mounting

◀ **BACK**



// QOP110353 // QOP110302
// QOP110362



// QOP110353 // QOP110302
// QOP110362 // QOP110342



// QOP110302



// QOP110353 // QOP110361
// QOP110524 // QOP110526
// QOP110560



// QOP110353 // QOP110306
// QOP110362 // QOP110538
// QOP110539

Chin Control Mounting

◀ **BACK**



// QOP110480 // QOP110470
// QOP110471



// QOP110483 // QOP110484



// QOP110483 // QOP110484



// QOP110485 // QOP110486



// QOP110485 // QOP110474



// QOP110487 // QOP110471
// QOP110526 // QOP110529
// QOP110561 // QOP110572
// QOP110571



// QOP110487 // QOP110474

Comments For Customisations



Q500 M/F/R with SEDEO PRO

€

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year

**Sunrise Medical GmbH**

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu