



HILITE XTR



THE PERFECT START TO MADE TO MEASURE



HILITE XTR

Functional, lightweight and comfortable.

The Hilite XTR is the perfect Titanium chair to experience made to measure for the first time, whether you're a new user or an experienced one. With over 30 measurements to be tailored to the user, the clinical benefits and comfort levels are life changing.



TECHNOLOGY

The 'REAL' Made To Measure

Only RGK offers the Real Made to Measure with over 30 different measurements to personalise and create the perfect fit wheelchair for you.

Functional Options

A full range of active and assist options such as side guards, anti tips, push handles, wheel locks, armrests and much more to suit your individual requirements.

The highest quality materials

Using only the finest materials, the wheelchair frame is made from commercially pure Titanium grade 2, which provides a super smooth ride and corrosion resistance.

Completely Personalised

The Hilite XTR comes with a huge array of options to choose from, ensures superior functional performance in this classic design wheelchair.

KEY FEATURES

AIRTECH UPHOLSTERY

25mm SEAT RAIL TUBE

TITANIUM GRADE 2 FRAME

CRASH TESTED

OPEN OR CLOSED FRAME DESIGN

TUBE TO TUBE CASTOR ARM



OPTIONS AND ACCESSORIES

The Hilite XTR boasts a huge array of options to provide independence for any type of user. From posture and clinical improvements to ergonomic support to adjustment settings, this wheelchair has it all.

Castor Assembly



Underseat Scissor Wheel Lock



Screw in Sideguards



Closed Frame Design



Tubular Titanium Footrest



Height Adjustable Push Handles



Live Life

The possibilities are endless

Enabling you to live your life to the fullest



TECHNICAL STUFF

DESCRIPTION	MIN	MAX
User Weight	N/A	115kg
Overall length with leg rest (no wheels)	400mm	700mm
Overall width (seat width)	250mm	500mm
Folded height (only for chairs with folding backs)	500mm	750mm
Total mass (with wheels)	7.5kg	15kg*
Mass of the heaviest part (transport weight)	5.0kg	12kg*
Static stability downhill	N/A	5° of incline
Static stability uphill	N/A	5° of incline
Static stability sideways	N/A	5° of incline
Seat plan angle	0°	45° of tilt
Effective seat depth	250mm	500mm
Seat surface height at front edge	380mm	550mm
Backrest angle	-10° of tilt	20° of tilt
Backrest height	150mm	450mm
Footrest to seat distance	200mm	520mm
Leg to seat surface angle	70° of angle	90° of angle
Armrest to seat distance	150mm	300mm
Front location to armrest structure	N/A	N/A
Hand rim diameter	24"-in	26"-in
Horizontal location of axle	25mm	150mm
Minimum turning radius	270mm	N/A

NOTE: The above is a guide. For any requirements outside of this table may result in a design change.

*Based on additional options, excluding third party wheelchair add-ons.



rgkwheelchairs.com

HEADQUARTERS

RGK Wheelchairs Ltd
Unit 8 a/b/c Ring Road Zone 2,
Burntwood Business Park,
Burntwood, Staffordshire,
WS7 3JQ, UK

Phone: +44 (0) 1543 670077
Email: info@rgkwheelchairs.com

