







| | Till yptor | |
|--|---------------------------------------|------------------------------------|
| Deliver to: | Deliver to: | |
| Name: | Name: | |
| Street: | Street: | |
| Town: Postal code: | Town: | Postal code: |
| E-Mail: | E-Mail: | |
| Tel: Fax: | Tel: | Fax: |
| Order date: Purchase order: | Marked for: | |
| Number required: | Sunrise order no: | |
| | ● = Standard | □ = Cost option |
| □ Order □ Calculation □ Built 4 me (le | onger leadtime) ./. = Price reduction | = No cost option |
| Reimbursement code: | Max. user weight: | 125 kg |
| Art. no. IC : 770 00 001 | USER WEIGHT | kg |
| Art. No. DE: 770 00 000 | Mandatory information to ensure | |
| | | |

Base Price

KRYPTON R: STANDARD CONFIGURATION

KRYPTON R

- · Rigid wheelchair
- Carbon fibre frame
- · Adjustable back angle
- Foldable backrest
- Adjustable COG
- Rear seat height adjustable +/- 1cm
- · Carbon fibre axle-stems
- · Carbon fibre axle tube





| Kry | pton | R |
|-----|------|---|
| , | | |

= Standard

O = No cost option

☐ = Cost option

CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - in mm; measured from outside frame tube to outside frame tube

| 320 | 340 | 360 | 380 | 400 | 420 | 440 | 460 |
|-----|-----|-----|-----|-----|-----|-----|-----|
| 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |

SEAT DEPTH SELECTION - in mm

| 340 | 360 | 380 | 400 | 420 | 440 | 460 | 480 |
|-----|-----|-----|-----|-----|-----|-----|-----|
| 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |

FRONT SEAT HEIGHT SELECTION - in dependance of front castor size, fork, frame height version The high frame is 40 mm higher than the low frame.

Front seat height in mm

| | | | 1 1011 | t seat | | | | | | | | | | | |
|----|---------------|----------|--------|--------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| | Fork (length) | Frame: | 430 | 440 | 450 | 460 | 470 | 480 | 490 | 500 | 510 | 520 | 530 | 540 | 550 |
| | 111 mm | | 0 | 0 | 0 | | | | | | | | | | |
| | 115 mm | low | | | 0 | | | | | | | | | | |
| | Frog Leg | T IOW | | | | | | | | | | | | | |
| 3" | One Arm | | | | | | | | | | | | | | |
| ľ | 111 mm | | | | | | 0 | 0 | 0 | | | | | | |
| | 115 mm | high | | | | | | | 0 | 0 | 0 | | | | |
| | Frog Leg |] '''g'' | | | | | | | | | | | | | |
| | One Arm | | | | | | | | | | | | | | |
| | 111 mm | | | 0 | 0 | 0 | 0 | | | | | | | | |
| | 115 mm | low | | | | | | | | | | | | | |
| 4" | Frog Leg | T IOW | | | | | | | | | | | | | |
| | One Arm | 1 | | | | | | | | | | 0 | | | |
| | 111 mm | | | | | | | 0 | 0 | 0 | 0 | | | | |
| | 115 mm | high | | | | | | | | | | | | | |
| | Frog Leg |] '''g'' | | | | | | | | | | | | | |
| | One Arm | | | | | | | | | | | | | | |
| | 111 mm | | | | | 0 | 0 | 0 | | | | | | | |
| | 115 mm | low | | | | | | | | | | | | | |
| | Frog Leg | T IOW | | | | | | | | | | | | | |
| | One Arm | 1 | | | | | | | | | | | | | |
| 5" | 111 mm | | | | | | | | | 0 | 0 | 0 | | | |
| | 115 mm | high | | | | | | | | | | | | | |
| | Frog Leg | high | | | | | | | | | | | | | |
| | One Arm | <u> </u> | | | | | | | | | | | | | |
| 6" | 115 mm | high | | | | | | | | | | | | | |
| ٥ | Frog Leg | high | | | | | | | | | | | | | |

REAR SEAT HEIGHT SELECTION - in dependance of rear wheel size

The rear seat height is adjustable \pm 10 mm.

The rear seat height can only be between 0 cm and 110 mm below the selected front seat height.

Caution: with rear seat heights 460 mm and above, the width between the sideguards increases by 35 mm

Rear seat height in mm

| Rear Wheels | 380 | 390 | 400 | 410 | 420 | 430 | 440 | 450 | 460 | 470 | 480 | 490 | 500 |
|-------------|-----------|------------------|------------------|-----------|-----|-----|-----|-----|-----|-----|-----|-----------|-----------|
| 24" | $O^{(1)}$ | O ⁽²⁾ | O ⁽²⁾ | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | $O^{(3)}$ | |
| 25" | | $O^{(1)}$ | $O^{(2)}$ | $O^{(2)}$ | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | $O^{(3)}$ |

(1): No adjustment of rear seat-height possible

(2):

No downwards-adjustment of rear-seat-height possible. Only +1 cm adjustable.

(3): No upwards-adjustment of rear-seat-height possible. Only -1 cm adjustable.

For more information www.SunriseMedical.eu © Sunrise Medical 2024

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

a) For the user: Allow space between sideguard and hip on each side of the chair

(approx. 10 - 20 mm altogether)

b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

a) For the user: When sitting upright: leave space for 2 - 3 fingers between

front edge of the seat upholstery and back of knee

(approx. 30 - 60 mm). Care must be taken with adj. upholstery.

b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

a) For the user: When sitting with 90° kneeangle:

From back of knee to bottom of shoe sole.

b) On the wheelchair: From front of seat tube to the ground .

(min. 30 mm longer than user measurement) Remember to include the cushion height.

Rear seat height:

b) On the wheelchair: From back of seat tube to the ground.

Note: Max. seat height difference rear to front is 130 mm

Back height:

a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.

b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.

Note: The user's adility to sit upright has a major influence on the back height

and it can therefore differ considerably


Back angle: Important: The back height depends on the degree of disability.

 0° is the angle 90° to the seat. Negative angles are ${<}90^{\circ}$

and positive angles are >90°.

Low Seat Frame vs. High Seat Frame:

The high frame is 40 mm higher than the low frame. So you can choose shorter and lighter forks to minimize the weight of your wheelchair.



Technical Data Overview:

Max. safe slope:

Max. overall height: 950 mm Max. overall width: Camber 25" rear wheel 24" rear wheel Max. overall length: 930 mm For 0° camber: SW + 170 mm SW + 170 mm Max. turning circle: 700 mm For 3°camber: SW + 230 mm SW + 240 mm

Wheelchair weight: from 6,2 kg
Dismantled weight: from 3,7 kg
Weight of heaviest part: 2,1 kg

All changes in the configuration might change the seat heights.

Manufacturing tolerances are +/- 10 mm for all measurements.

© Sunrise Medical 2024

WEIGHT

We want to be very transparent on weight. The following information will be a guide to the lightest possible Krypton R: how to achieve a 6,2 kg wheelchair, and how to order the lightest possible configuration for your Krypton R

Krypton R: 6,2 kg configuration Transport weight: 3,7 kg

Frame: high frame, 92°, 30 mm inset, seat width: 360 mm, seat depth: 400 mm. Seat height: front 470 mm, rear 420 mm. Seat-sling Airflow superlight.

Wheel-lock: lightweight compact wheel-lock

Footrest: tubular footrest



Backrest: Carbon backrest, EXO superlight backrest-upholstery, without push-handles

Sideguards: Carbon fibre, slim-style

Rear wheels: 24" Proton wheels with Aluminium handrims (coated or anodized) and Schwalbe One

Forks and casters: 111 mm carbotecture fork, 3" caster

Please mind that there is a 5% tolerance on all weight specifications due to weight deviations on mounting parts (tyres, caster-wheels etc.)

LIGHTWEIGHT-GUIDE FOR ORDERING THE KRYPTON R

Sunrise Medical offers a wide variety of options for the Krypton R to suit every user's needs. However, most options add significant weight to the wheelchair. To make this more transparent when filling in the order-form, and to assist you in ordering the lightest possible wheelchair, we indicate the approximate weight-increase of the options on the order form:

- (L) indicates the lightest possible option within the range of options
- (+) indicates that this option adds more than 500 g to the wheelchair's weight compared to the lightest option within the range

LIGHTWEIGHT OPTIONS FOR THE KRYPTON R

We developed several options for the Krypton R to save weight on the chair without sacrifing daily-use suitability, comfort or durability



Carbon backrest (KRR030006): full carbon backrest with very close cross-tube (20mm crosstube to backrest-upholstery) (-100 g compared to Standard backrest)



Seat sling "airflow" superlight (KRR020019): fully adjustable seat-sling with composite buckles and lighter straps (max. user weight 85 kg) (-100 a compared to seat sling "airflow")



EXO superlight (KRR030005): lightweight version of the EXO backrest upholstery. EVA-foam padding, adjustable straps, wider cut of the upholstery to cover backrest (-160 g compared to Std. EXO

backrest-upholstery)





Axle-plua superliaht (KRR070049): (only 0° camber): glued into the axle-tube, no adjustability of the track-width possible (-160 g compared to widthadjustable axle-plug KRR070600 / 603)

| Kryptor | n R | | |
|-------------------------|--------|--|----------|
| | • | · · · · · · · · · · · · · · · · · · · | ice in € |
| 0 10" | | SEAT UPHOLSTERIES AND CUSHIONS | |
| Seat Sling KRR020012 | 0 | Cost aling "airflow" with 4 utility has | |
| KRR020012 KRR020017 | | Seat sling "airflow" with 1 utility bag Seat sling "airflow" with 2 utility bags | |
| KRR020017 | | Seat sling "airflow" with utility bag and catheter bag (from SW 340mm) | |
| KRR020019 | | Seat sling "airflow" superlight. Without bags, max user weight 85kg (1) (L) | |
| KRR020000 | 0 | Without seat sling | |
| | | | |
| Seat cushio | | | |
| Foam: (medi | um so | oft) Thickness: KRR020101 30 mm KRR020102 50 mm | |
| JAY Cushio | ns - f | or complete product range see Jay order forms | |
| JAY020002 | | JAY Basic | |
| JAY020003 | | JAY Soft Combi P | |
| JAY020001 | | JAY Easy Visco, Curved Base | |
| JAY020007 | | JAY Easy Fluid, Curved Base | |
| JAY020004 | | JAY Lite | |
| Seat access | | | |
| KRR090018 | | Positioning belt (only with Standard backrest) BACKREST | |
| KRR030001 | 0 | Standard backrest (aluminium, black) | |
| | | Choose back height factory setting (mm). 250 - 300 300 - 350 350 - 400 400 - 450 | |
| | | Height adjustments possible in the O 250 O 300 O 350 O 400 | |
| | | range 250-300; 300-400; 400-450 mm. O 275 O 325 O 375 O 425 | |
| | | O 300 O 350 O 400 O 450 | |
| KRR030006 | | Carbon backrest - fixed height, crosstube very close to the upholstery, "Krypton" rubber grip (2) (L) | |
| | | Choose back height (mm) | |
| | | no height adjustment possible. O 250 O 325 O 400 | |
| | | O 275 O 350 O 425 | |
| | | O 300 O 375 O 450 | |
| KRR030012 | 0 | Backrest, angle-adjustable (non foldable) (-31° to +15°) (only in black anodizing) | |
| KRR030015 | 0 | Backrest, angle-adjustable, fold-down (-31° to +15°) (only in black anodizing) | |
| KRR030102 | | Backrest, angle-adjustable, fold-down, locked when folded, colour: black (-31° to +15°) colour: Orange Silver Red Blue KRR030036 KRR030037 KRR030038 KRR030039 | |
| Duchhandla | a. ob | ages only one ention | |
| KRR030200 | 0 | oose only one option Omit push handles (L) | |
| KRR030201 | 0 | Push handles long | |
| KRR030203 | | Fold-down push handles (from 325mm backrest-height) | |
| KRR030204 | | Height adj. push handles (+) (only with Standard backrest; from 325mm backrest height) | |
| KRR030205 | | Height adj. push handles for Jay backs (+) (only with Standard backrest; from 325mm backrest height) | |
| Backrest up | holst | ery | |
| KRR030316 | 0 | EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version) | |
| KDD020247 | | Colour, binding tape: O Black O Orange O Silver O Red O Blue | |
| KRR030317 | | EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version) Colour, binding tape: O Black O Orange O Silver O Red O Blue | |
| KRR030005 | | EXO superlight, padded, adjustable, black (EXO Evo version) (L) | |
| KRR030300 | 0 | Colour, binding tape: O Black O Orange O Silver O Red O Blue Without backrest upholstery | |
| | | | |
| (1): | Max. | . user weight 85 kg, only available up to seat-width 400 mm For more informati | ion |
| (2): | Not | with height adj push handles, not with Jay J3 backs, not with positioning-belt www.SunriseMedica | al.eu |
| | (1)= | lightest option (+)= adds 500 g or more to the wheelchair-weight © Sunrise Medical | 2024 |

SEAT UPHOLSTERIES AND CUSHIONS



KRR020012 Seat sling, airflow



KRR020019 Seat sling, airflow, superlight

BACKREST



KRR030006 Carbon Backrest



KRR030201 Push handles long



KRR030316 EXO backrest sling (EXO Evo version)



KRR030102 Backrest, angle-adjustable, lock when folded



KRR030203 Push handles, fold-down



KRR030317 EXO PRO backrest sling (EXO Evo version)



KRR030204 Push handles , height-adjustable



KRR030005 EXO superlight backrest sling (EXO Evo version)

FOOTREST AND FOOTPLATES



KRR050029 Tubular footrest lightweight



KRR050025 Platform footrest LW composite plate



KRR050026 Platform footrest LW carbon plate



KRR050125 Performance footrest alu angle adj



KRR050127 Foot positioning plate

| Krypton | R | | |
|------------------------|----------|---|---|
| | <u>•</u> | · | rice in € |
| IAV I2 Back | Sh | JAY 3 BACKS AND ACCESSORIES - only with Standard backrest (+) | |
| Mid Thoracic: | | nallow Contour ☐ JAY030010 Medium - 340 mm ☐ JAY030033 Tall - 380 n | mm |
| Upper Thoracic. | | ☐ JAY030012 Short - 420 mm ☐ JAY030019 Medium - 460 mm ☐ JAY030025 Tall - 500 n | |
| Shoulder Heig | | ☐ JAY030015 Short - 530 mm ☐ JAY030022 Medium - 570 mm ☐ JAY030028 Tall - 610 n | |
| | | | |
| JAY J3 Back | | | 1 |
| Mid Thoracic: | | □ JAY030011 Medium - 340 mm □ JAY030018 Tall - 380 n □ JAY030013 Short - 420 mm □ JAY030020 Medium - 460 mm □ JAY030026 Tall - 500 n | |
| Upper Thorac | | □ JAY030013 Short - 420 mm □ JAY030020 Medium - 460 mm □ JAY030026 Tall - 500 n □ JAY030016 Short - 530 mm □ JAY030023 Medium - 570 mm □ JAY030029 Tall - 610 n | ***** |
| Shoulder Heig | | * | 1111 |
| | | ack - Shallow Contour | |
| Lower Thorac | | □ JAY060001 Short - 170 mm □ JAY060002 Tall - 240 n | |
| Mid Thoratic: | | □ JAY060003 Short - 300 mm □ JAY060004 Tall - 380 n | nm |
| Accessories | for J | JAY J3 Backs | |
| WHT030415 | | □ Whitmyer Headrest, COBRA XTRA, flip back (only with Jay 3 b | backs) |
| WHT030416 | | □ Whitmyer Headrest, AXYS, Standard (only with Jay 3 b | backs) |
| JAY090058 | | □ JAY Accessory bag (only with Jay 3 Carbon backs) | |
| | | OIDEOUADDO AND ADMOPOTO | |
| 1/20040404 | _ | SIDEGUARDS AND ARMRESTS Sideward earlier block | a |
| KRR040101 | | • | (L) |
| KRR040102 | | Sideguard, aluminium, cold weather protection (colour: black, matt) | 21 |
| KRR040103 | | Sideguard, aluminium, fender, cold weather protection (colour: black, matt) (min rear seat height 390 | |
| KRR040104 KRR040107 | | Sideguard, carbon fibre, fender, cold weather protection, black Sideguard, composite, removable, 21 cm high (31 cm long), black (min rear seat height 390) |) mm) |
| KRR040107 KRR040052 | | | |
| KRR040052 KRR040117 | | Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable Tubular armrest, Alu, swing away (sideways), padded, removable (+) | |
| KRR040117 KRR040000 | 0 | w/o Sideguard | ' |
| NRNU40000 | | W/O Sideguard | |
| | | FOOTREST AND FOOTPLATES | |
| KRR050029 | 0 | | (L) |
| KRR050030 | | Footrest, tubular, lightweight, aluminium, with cover | |
| KRR050025 | | Platform footplate lightweight angle & depth adjustable composite | |
| KRR050024 | | Platform footplate lightweight angle & depth adjustable aluminium (colour: black, matt) | |
| KRR050026 | | Platform footplate lightweight angle & depth adjustable carbon | |
| KRR050125 | | Performance footplate angle adjustable (3) | |
| KRR050128 | | Footboard, platform, plastic, angle and depth adjustable, flip-up (to the back) (4) | |
| KRR050127 | | Adjustable side protection, to create a side-stop on the footboard (5) | |
| | | LOWER LEG LENGTH (LLL) Attention LLL has to be 20 mm loss than cost beight front; all my | |
| | | LOWER LEG LENGTH (LLL) Attention: LLL has to be 30 mm less than seat height front; all me | asurements in mini |
| | 0 | 360 0 370 0 380 0 390 0 400 0 410 0 420 0 430 | |
| | 0 | 440 0 450 0 460 0 470 0 480 0 490 0 500 0 510 | |
| | | FRAME VERSIONS | |
| | | | |
| | | Frame length | |
| KRR010206 | 0 | Frame standard | |
| KRR010207 | 0 | Frame long (20 mm additional frame length) (not in combination with SD 480 mm) | |
| | | Frame angle | · · · |
| KRR010012 | 0 | 100° frame | |
| KRR010013 | | 92° Performance frame | |
| | | Frame inset | |
| KRR010101 | 0 | Inset: Small (10 mm) | |
| KRR010103 | | Inset: Performance (30 mm) (not available for seat width 320 mm) | |
| | | ` ; | |
| KRR010301 | 0 | Castor position (see explanation on next page) Castor house position wide | |
| KRR010301 KRR010302 | 0 | Castor house position narrow (not with seat-depths < 400mm, please see recommendations on next page) | - - |
| NAMOTOGOZ | | Castor flouse position flattow (not with seat-depuis > 400mm, please see recommendations on flext page) | |
| (1) | With | 24" rear wheel min. rear seat height 410 mm; with 25" rear wheel min. rear seat height 430 mm | |
| (2) | | with COG 56 mm | For more information |
| (3) | _ | with SW 320 mm and 30 mm inset | www.SunriseMedical.eu |
| (4) | | r weight limitation 75 kg, not with 4" castor | © Sunrise Medical 2024 |
| (1) (5) | | r for platform footplates lightweight angle & depth adjustable | o cumico modifica |
| (6) | _ | lightest option (+)= adds 500 g or more to the wheelchair-weight | |

SIDEGUARD & ARMREST



KRR040101 Slim style sideguard



KRR040104 Carbon fibre, cold weather protection, fender



KRR040107 Sideguard, composite, removable (black)




KRR040117 Tubular armrest, swing away, removable, padded




KRR040052 Single-post sideguard, armpad (33 cm)

FRAME VERSIONS: CASTOR HOUSE POSITIONS


Position wide:



Further to the front and wider position of the



Position narrow:





22 mm further rearwards and 24 mm per side more narrow position of the casters for a more narrow front-end and improved manoeuverability in tight spaces. Reduced front-stability compared to <u>position wide</u>, we recommend this option in combination with seat- height differences (front to rear) of 40mm or more or in combination with KRR010207 "frame, long".



KRR010301 Castor house position wide



KRR010302 Castor house position narrow



Castor house position wide vs. narrow

| Krypton | R | | |
|-----------------|------|---|------------------------|
| '' | • | = Standard O = No cost option □ = Cost option | Price in € |
| | | CENTRE OF GRAVITY (COG) - all measurements are in mm; measured from front of back post; choose of | nly one option |
| | | | |
| COG56 | 0 | 56 mm | |
| COG68 | 0 | 68 mm | |
| COG80 | 0 | 80 mm | |
| COG92 COG104 | 0 | 92 mm 104 mm | |
| COG 104 | | 104 IIIII | |
| | | FRAME ACCESSORIES | |
| KRR090031 | | Suspension up to 70 kg (1) (+) (chair will be built using aluminium axle-stems and axle-tube) | |
| KRR090032 | | Suspension up to 120 kg (1) (+) (chair will be built using aluminium axle-stems and axle-tube) | |
| | | Tip assist (step tube) (2) | |
| KRR090001 | | Left KRR090002 □ Right | |
| KRR090003 | | Pair (+) (not available with anti-tip) | |
| KRR090009 | | Crutch holder (only in combination with tip assist) | |
| | | Anti-tip tubes swing-away (2) | |
| KRR090004 | | Left KRR090005 □ Right (+) | |
| KRR090006 | | Pair (+) (not available with tip-assist) | |
| | | | <u> </u> |
| KRR090007 | | Quickie style anti-tip tubes pluggable (2) (+) (Pair) | |
| | | Anti tip flip up (user operated) (2) | |
| KRR090033 | | Left KRR090008 Right (+) (Please choose only one) | |
| KRR090010 | | Transit wheels (2) (+) | |
| KRF090016 | | Frame protection, black (textile, with velcro) | |
| KRR090577 | | Frame protection, Plastic, glued to the frame, black | |
| | | | <u> </u> |
| | | FORKS (compare with your choice on the first page) | |
| KRR080013 | 0 | Castor fork, Carbotecture ®, 111mm Colour: Black (L) | |
| KDD00000 | _ | Colour: Silver Orange Red (3) Blue (3) Castor fork. Aluminium. 115mm Colour: Black | |
| KRR080029 | 0 | Castor fork, Aluminium, 115mm Colour: Black Colour: | |
| KRR080008 | П | Castor fork, aluminium, one-arm Colour: Black (only with 4" & 5" Alu soft, Alu solid) | |
| | _ | Colour: ☐ Silver ☐ Orange ☐ Red ☐ Blue | |
| KRR080007 | | Castor fork, aluminium, suspension (Frog Legs ®), black (+) | |
| | | | |
| | | CASTORS (compare with your choice on the first p | page) |
| KRR080301 | 0 | Castor, solid (hard), rim: composite, black (3"-6", not with One arm fork) | |
| KRR080302 | 0 | Castor, solid (soft), rim: composite, black (only 5" and 6", not with One arm fork) | |
| KRR080303 | | Castor, solid (hard), light-up, rim: composite, black (not available in 6", not with One arm fork, only 4" ar | nd 5") |
| KRR080304 | | Castor, solid (hard), rim: aluminium Colour, rim: Silver (only 4" and 5") Castor, solid (soft), rim: aluminium Colour, rim: Silver (only 4" and 5") | |
| KRR080306 | | Colour, rim: Black | |
| | | Colour, IIII. Li black Li Grange Li Ned Li blue | |
| | | REAR WHEELS (compare with your choice on the first page; choose only one wheel | type) |
| | | NOTE: for users with higher user weight than 100 kg, we recommend cross spoked wheels | |
| KRR070000 | 0 | Without rear wheel | |
| KRR070001 | 0 | Rear wheel, Universal, 36 spokes (silver), cross-wise spokes (+) | |
| KRR070002 | 0 | Rear wheel, Design, 36 spokes (silver), straight spokes (only 24") | |
| KRR070003 | 0 | Rear wheel, Lightweight, 24 spokes (silver), straight spokes (Std.: black hub/rim; colour options last p | page) |
| KRR070004 | | Rear wheel, Proton, 24 spokes (black), straight spokes (L) weight limit 100kg) | |
| KRR070005 | | Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver | |
| KRR070007 | | Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black | |
| | | | |
| | | | |
| (1): | Not | possible with carbon axle-stems and axle: chair will be built with aluminium axle-stems and axle-tube | For more information |
| (2): | • | with KRR07049: camber adapter 0°, superlight. | www.SunriseMedical.eu |
| (3): | Red | and blue only in combination with soft castor | © Sunrise Medical 2024 |
| 1 | (L)= | lightest option (+)= adds 500 g or more to the wheelchair-weight | |

FRAME ACCESSORIES



KRR090032 Suspension



KRR090001 Tip assist



KRR090004 Anti-tip tubes SA



KRR090007 Quickie style anti tip tubes pluggable



KRR090033 Anti tip flip up (user operated)



KRR090010 Transit Wheels



KRR090010 Frame protection,Plastic, glued to the frame, black

FORKS



KRR080013 & KRR080306 Castor fork, Carbotecture & Castor solid (soft) aluminum



KRR080008 Castor fork, aluminum, one-arm, red

REAR WHEELS & HANDRIMS



KRR070001 Universal wheel



KRR070002 Design wheel



KRR070003 Lightweight wheel



KRR070004 Proton wheel



KRR070211 Handrims Surge®



KRR070208 Handrims Maxgrepp

TYRES



KRR070101 Schwalbe Right Run



KRR070107 Schwalbe One



KRR070104 Schwalbe Marathon Plus



KRR070102 Tyre, puncture proof

WHEEL LOCKS



KRR060001 Standard brake



KRR060002 Knee lever brake



KRR060023 Compact brake EVO



KRR060024 Lightweight compact brake EVO

| Kryptor | n R | | |
|------------------------|----------|---|------------------------|
| | • | | ice in € |
| KRR070101 | 0 | TYRES (choose only one tyre Tyre, pneumatic, fine profile, Schwalbe Right Run (not in 25") | іуре) |
| KRR070107 | | Tyre, pneumatic, slick, Schwalbe One (L) | |
| KRR070102 | <u> </u> | Tyre, puncture proof, tread (+) | |
| KRR070104 | | Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution (+) | |
| KRR090163 | | Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube (not with solid) | tyres) |
| | | HANDRIMS | |
| KRR070201 | 0 | Handrims, aluminium, powder-coated, clear (L) | |
| KRR070202 | 0 | Handrims, aluminium, anodised, silver (L) | |
| KRR070203 | | Handrims, aluminium, powder-coated, coloured (L) (select colour on last page) | |
| KRR070204 | 0 | Handrims, stainless steel, silver (+) | |
| KRR070205 | | Handrims, titanium, untreated | |
| KRR070206 | | Handrims, titanium, powder-coated, clear | |
| KRR070207 | | Handrims, Supergrip ® (+) | |
| KRR070208 | | Handrims, Maxgrepp ® (+) (only 24") | |
| KRR070212 | | Handrims, Ellipse 3R, black (+) (ergonomic handrims with rubber- | -strip) |
| KRR070250 | | Handrim cover, rubber (not with 25") | |
| | | Hand rim mounting | |
| KRR070307 | | Handrim mount, super narrow (7 mm) | |
| | (only v | ith stainless steel; not in combination with Spinergy, Proton and Lightweight Wheels and stainless steel Handrim, Not with Maxg | repp) |
| KRR070317 | 0 | Handrim mount, narrow (17 mm) (Not with Maxe | grepp) |
| KRR070327 | 0 | Handrim mount, wide (27 mm) | |
| | | WHEEL CAMBER (choose only one option) | |
| KRR070049 | 0 | 0° superlight, no width-adjustment (1) (L) | |
| KRR070600 | 0 | 0° width-adjustment +/-15mm | |
| KRR070603 | 0 | 3° width-adjustment +/-15mm | |
| | | SPOKE GUARDS (only available in 24") (+) | |
| KRR070831 | | Transparent | |
| KRR070823 | | Princess KRR070813 □ Galactica KRR07081 □ Jungle Fever | |
| KRR070816 | | Happy Lions KRR070821 □ Pirates KRR07081 □ Feather | |
| KRR070819 | | Monkeys KRR070814 □ Graffiti KRR07080 □ Chrome wheel | |
| KRR070820 | | Oriental KRR070802 □ Basketball KRR07080 □ Basketball 2 | |
| KRR070829 | | Superhero KRR070812 □ Flying water KRR07082 □ Space | |
| KRR070817 | | Hot wheels KRR070822 □ Pool KRR070827 □ Silk | |
| KRR070830 | | Tiger KRR070833 □ Zebra | |
| | | WHEEL LOCKS (Choose only one option) | |
| KRR060001 | 0 | Standard wheel lock | T |
| KRR060002 | | Knee lever wheel lock | |
| KRR060010 | | ☐ Extension for knee lever wheel lock, slip-on, swing-away (only with KRR06) | 0002) |
| KRR060023 | | Compact brake EVO | 3332) |
| KRR060024 | | Lightweight compact brake EVO (2) (L) | |
| KRR060006 | | One arm wheel lock (left hand actuation) (+) | |
| KRR060007 | | One arm wheel lock (right hand actuation) (+) | |
| | | ACCESSORIES | |
| KRR090000 | | ACCESSORIES Tool kit | |
| KRR090000 KRR090024 | | High pressure pump | |
| KRR090024 KRR090025 | | Air pump | |
| KRR090025 KRR090026 | | Mobile phone pocket | |
| KRR090030 | | Back pack | |
| | | | |
| (1): | Can | ber adapter glued into the axle-tube. No adjustment of the track-width possible. | For more information |
| (2): | Adji | stment of COG ± 2 cm on a chair in the field. | www.SunriseMedical.com |
| | (L)= | lightest option (+)= adds 500 g or more to the wheelchair-weight | © Sunrise Medical 2024 |

| Krypton | R | | | | | | |
|-------------------------------------|------|-----------------------|----------|--|------------------------------------|---------------------------------------|------------------------|
| - | • | = Standard | 0 | = No cost option | □ = Cost option | | |
| | | | | | HUB/ RIM OF THE STANDAR | | |
| | | F | landrir | n (2) Hub/rim (1) | | Handrim (2) | Hub/rim (1) |
| 02 Brilliant Silv | er | | | | 27 Graphitemetallic | | |
| 03 Silver Matt | 0. | | | | 29 Black | | |
| 06 Sun Yellow | | | | | 31 Matt Black | | |
| 14 Blue | | | | | 32 Signal White | | |
| 15 Brilliant Mar | rine | | | | 33 Pink | | |
| 16 Matt Marine | ; | | | | 34 Dark Chocolate Metallic | | |
| 19 Brilliant Pur | ple | | | | 35 Skyblue | | |
| 21 Red Ruby | | | | | 36 Cyan-Metallic | | |
| 22 Brilliant Red | t | | | | 37 Marsred (Orange) | | |
| 23 Matt Red | | | | | 40 Kiwi | | |
| 24 Tutti Frutti C | Orai | nge | | | 41 Flamingo-Glitter | | |
| | | | | <u>.</u> | 42 Matt-Anthracite | | |
| | | | | | 46 Mamba Green (new) | | |
| (1) Only with un | ivei | rsal cross spoked | rear wh | ieels | (2) only with KRR070203 | , handrim aluminium, po | wder coated |
| | | FRAME-COI | LOUF | RS | | | |
| KRR100139 | 0 | Natural Carbo | | | | | |
| KRR100137 | | Transparent, g | gloss | | | | |
| KRR100138 | | Transparent, r | natt | | | | |
| KRR100122 | | Black (gloss) | | | | | |
| KRR100123 | | Matte black | | | | | |
| | | EDAME CD | A DUI | CC "Elving Dival" | | : | !! -4: -l\ |
| KRR100590 | | Frame graphic | | | (scratch resistant, non-removeable | , includes the "Krypto | n" sticker) |
| | | Frame graphic | | | | | |
| | | Frame graphic | | | | | |
| | | - | - | | | | |
| | | | <u> </u> | ERS "SOPUR"/ "(| QUICKIE" and "Krypton" (3) | | |
| | 0 | black | | | | | |
| | 0 | silver | | | | | |
| | 0 | orange | | | | | |
| | 0 | red | | | | | |
| NRR 100554 | 0 | blue | | | | | |
| | | COLOURS F | OR I | LIGHTWEIGHT W | HEELS 24" (4) | | |
| KRR100336 | | Rear wheel, li | ghtwei | ght anodized, orange | (4) | | |
| | | | | ght anodized, blue (4 |) | | |
| | | | | ght anodized, red (4) | | | |
| | | | | "Flying Pixel" are chosei spokes. Only with light w | | | |
| 4) Allouiseu co | ioui | s for flub and filli, | DIACKS | spokes. Omy with light w | eight wheel in 24 | | |
| Comments | for | special buil | ld: | | | | |
| | | | | | | | i sa di Salana |
| | | | | | | le . | |
| | | | | | | DL | uilt <mark>4</mark> me |
| | | | | | | | |
| | | | | | | | |
| | | Base Price |) | | | | |
| Option price | | tronon and | | | | | |
| Packaging a Exclusive V <i>I</i> | | панѕроп | | | | | |
| VAT. | 11 | | | | | | |
| VAT. Total Pric | | | | | | | |
| | | Customer Serv | vice. | | | | For more information |
| Address: | ui ' | Castonier Gerv | 100 | | E-mail: | | vw.SunriseMedical.com |
| raaress. Tel. | | | | | E-commerce: | · · · · · · · · · · · · · · · · · · · | Sunrise Medical 2024 |
| Fax. | | | | | Website: | ı . | III WIJUI AVAT |

FRAME COLOURS



- Transparent, gloss Transparent, matt Black, gloss Black, matt

FRAME GRAPHIC "FLYING PIXELS"



From left to right:

- Frame graphics Flying Pixel, red Frame graphics Flying Pixel, silver Frame graphics Flying Pixel, blue

24 " LIGHT WEIGHT WHEEL COLOURS



KRR100336 Lightweight wheel, anodized orange



KRR100337 Lightweight wheel, anodized blue



KRR100338 Lightweight wheel, anodized red