S	IJ	N	T	2 1	S	E
						L.

Technical / News Bulletin

BULLETIN DETAILS						
Product Group:	Scooters					
Model Affected:	Elite2 XS and Elite2 Plus					
Effective Date: See Below	Serial No. Effective From:	N/A				
Bulletin Subject:	New Masterseat on Elite2 XS and Elite2 Plus (black versions only)					

As part of our continuous improvement programme, we are updating the seats on the Elite2 XS and Elite2 Plus.

The current seats will be replaced with the Masterseat. This will align the seats used on the Trophy and Elite2 family. The sliding mechanism under the seat has also been updated.

This update is only applicable for the black Elite2 XS and Elite2 Plus introduced in 2021.





Old Seat

New Masterseat

Features of the new Masterseat

- High backrest
- Easy adjustable lumbar supports

Comparable with the current seat

- Easy adjustable back angle
- Angle adjustable comfort armrests for setting the ideal height
- Sliding and turning mechanism for easy transfers

Spare Parts Elite2 XS / Plus

Most spare parts stay available:

Old Spare Part (will be available until January 2023) EB10284 – Eblo C1 Seat Two-tone

New Spare Part

1002861 – Masterseat with armrests
9012064 – Masterseat with sliders Elite2
(The masterseat itself on the Trophy and Elite will be the same, but the sliding mechanism will be Elite2 or Trophy specific caused by a difference in the seat post dimensions)

The new Masterseat will also fit on existing Elite2 Family Scooters; in this case, sliders need to be replaced as well.

Spare Parts Trophy

All spare parts stay available, except below:

Old Spare Part (will be available until October 2021) 1366 – Plastic cover for seat fixation strip

New Spare Part

9012384 – Plastic cover for seat fixation strip

For any queries on any of these updates please contact your local Area Sales Manager.

Customer Services	Customer Support	Technical Service Centre	E-mail	
Telephone:	+44 (0)845 605 66 88 Option 1-1	+44 (0)845 605 66 88 Option 1-2	enquiries@sunmed.co.uk	
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk	