

Combination Agreement

The below mentioned products are in principle appropriate to be combined with the below mentioned **LECKEY** products provided that the combination is made expertly and according to the restrictions of both manufacturers. The restrictions can be asked at the manufacturers.

The following Leckey seats and Sunrise bases have NOT been crash tested together as a combination but have been crash tested individually. It is therefore the responsibility of the end provider when issuing the equipment and any risk assessed and documented by them.

SUNRISE MEDICAL product	
Zippie TS (folding)	<p>Mygo Seat 1 - Leckey Interface A – 117-775 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction Seat User manual – Mygo LS151-07</p> <p>Mygo Seat 2 - Leckey Interface A – 117-775 Work Instruction – Mobility Interface 'A' – Generic fitting instruction Seat User Manual – Mygo LS151-07</p>
Zippie RS	<p>Mygo Seat 2 - Leckey Interface A – 117-775 Work Instruction – Mobility Interface 'A' – Generic fitting instruction Seat User Manual – Mygo LS151-07</p> <p>Mygo Seat 3 – Leckey Interface A – 117-775 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction Seat User Manual - Mygo LS151-06</p> <p>KIT Seat 1 – Leckey Interface 'A' – 133-790 Work Instruction internal assembly and customer – Mobility Interface 'L' – Generic fitting instruction Seat User Manual – Kit LS243-02</p>
Quickie TS	<p>Squiggles Seat - Leckey Interface A – 120-860 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction Seat User manual – Squiggles LS203-09</p> <p>Mygo Seat 1 - Leckey Interface A – 117-775 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction</p>

	<p>Seat User manual – Mygo LS151-07</p> <p>Mygo Seat 2 - Leckey Interface A – 117-775 Work Instruction – Mobility Interface 'A' – Generic fitting instruction Seat User Manual – Mygo LS151-07</p>
Quickie SR45	<p>Mygo Seat 1 - Leckey Interface A – 117-775 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction Seat User manual – Mygo LS151-07</p> <p>Mygo Seat 2 - Leckey Interface A – 117-775 Work Instruction – Mobility Interface 'A' – Generic fitting instruction Seat User Manual – Mygo LS151-07</p> <p>Mygo Seat 3 – Leckey Interface A – 117-775 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction Seat User Manual - Mygo LS151-06</p>
Quickie Iris	<p>Mygo Seat 3 – Leckey Interface K – 117-775 Work Instruction internal assembly and customer – Mobility Interface 'A' – Generic fitting instruction Seat User Manual - Mygo LS151-06</p>
Quickie Salsa Mini	<p>Mygo Seat 2 - Leckey Interface 'Salsa' – 137-860 Work Instruction – Generic fitting instruction Seat User Manual – Mygo LS151-07</p> <p>Mygo Seat 3 - Leckey Interface 'Salsa'– 137-860 Work Instruction – LWI-137-074 Seat User Manual – Mygo LS151-06</p> <p>Mygo Max seat - Leckey Interface 'Salsa' – 137-873 Work Instruction – Generic fitting instruction Seat User Manual – Mygo LS320-02</p> <p>KIT Seat 2 – Leckey Interface 'Salsa' – 137-873 Work Instruction – Generic fitting instruction Seat User Manual – Kit LS243-02</p>
Quickie Salsa M2	<p>Mygo Seat 2 - Leckey Interface 'Salsa' – 137-860 Work Instruction – Generic fitting instruction Seat User Manual – Mygo LS151-07</p> <p>Mygo Seat 3 - Leckey Interface 'Salsa'– 137-860 Work Instruction – LWI-137-074 Seat User Manual – Mygo LS151-06</p>

Mygo Max seat - Leckey Interface 'Salsa' – 137-873

Work Instruction – Generic fitting instruction
Seat User Manual – Mygo LS320-02

KIT Seat 1 – Leckey Interface 'Salsa' – 137-873

Work Instruction internal assembly and customer – Mobility
Interface 'L' – Generic fitting instruction
Seat User Manual – Kit LS243-02

KIT Seat 2 – Leckey Interface 'Salsa' – 137-873

Work Instruction internal assembly and customer – Mobility
Interface 'L' – Generic fitting instruction
Seat User Manual – Kit LS243-02

Please note that mounting has to be effected according to the **LECKEY** mounting instructions by authorized and skilled staff only. Restrictions from the **LECKEY** bracket list have to be observed.

The maximum torque specified in the mounting instructions of **LECKEY** has to be followed. If not specified the maximum torque of parts mounted directly on the frame should not exceed 7Nm.

The maximum values mentioned in the user manual of both manufacturers must be observed.

Machining, bending, welding or bracing on any safety relevant components is not allowed.

Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch
Germany
Malsch,

Senior Director European Mobility R&D

JAMES LECKEY DESIGN

19C Ballinderry Road
Lisburn
Northern Ireland
BT28 2SA



