



Technical / News Bulletin
UK

BULLETIN NO.:

TB UK180

DATE ISSUED

01/02/2022

BULLETIN DETAILS

Product Group:	QUICKIE	
Model Affected:	Q400 M Sedeo Lite, Q500 M Sedeo Pro & Q700 M Sedeo Pro	
Effective Date: See Below	Serial No. Effective From:	N/A
Bulletin Subject:	Q-Series Powered Seat Elevator & Tilt Module 50°	

At Sunrise Medical we are committed to our continuous improvement programme. As part of this program, we are pleased to introduce a new lift and tilt option that includes a 50° tilt to the Q400 M Sedeo Lite, Q500 M Sedeo Pro and Q700 M Sedeo Pro.

This option will provide the client with a higher angle for pressure relief and positioning, without having to compromise on the powered lift function.

With immediate effect, the powered lift & tilt module 50° is now available to order through the online configurator and on the paper order form. It can also be ordered as a spare part to be retro-fitted to an existing chair, please contact technical@sunmed.co.uk for further information.



QPL130006

Seat Elevator (300 mm)
(12 Inch) / Tilt Module
(Tilt 0 Degree - 50
Degree), Powered

NEW Order Form Option Code

xxx130006 Seat Elevator (300mm) / Tilt Module (Tilt 0° - 50°), Powered

NEW Spare Part Numbers

20012011-003 Lift / Tilt 50° Module, including mounting accessories

20010147 Tilt Actuator 50°

20010146 Cable Track

Existing Spare Parts

20010136 Lift Actuator

20010140 Cover

20010141 Micro Switch

For any queries on any of these updates please contact your local Area Sales Manager.

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 1-1	+44 (0)845 605 66 88 Option 1-2	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk