

# Pull Device StreetJet



Status:03/2022

## ORDER FORM

Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town / PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	Max. user weight: 120 kg		
Purchase Order No.:			User Weight: kg		
Art. No.: 47800001			Reimbursement Code: xxx		
Order Date:	Day	Month	Year		

Base Price

€ xxx,-

## Technical Data Overview

<b>Max. user weight</b>	120 kg	<b>Max. speed forward</b>	6 kph; optional 10 kph or 15 kph
<b>Front frame width</b>	220 – 425 mm <sup>(1)</sup>	<b>Wheel size</b>	12" pneumatic hybrid tyre / 16" pneumatic hybrid tyre
<b>Device IP Rating</b>	IPX4	<b>Motor power (max. power)</b>	12" 250 W (630 W) / 16" 250 W (750 W)
<b>Battery weight</b>	2,6 kg	<b>Battery type</b>	Li-Ion, 36V, 14 Ah, intelligent battery management system
<b>Total weight w/o battery</b>	from 12,9 kg <sup>(3)</sup>	<b>Brake type</b>	Twin-brake, system: mechanical disc-brake and separate proportional electric brake Optional Triple-brake, system: dual mechanical brake, and separate electric brake
<b>Total weight with battery</b>	from 15,5 kg <sup>(3)</sup>	<b>Operating Voltage</b>	36 Volt
<b>Charger</b>	36V, 2A	<b>Climbing power</b>	up to 11% (6°)
<b>Range of action</b>	up to 50km <sup>(2)</sup>		

(1) Measured inside tube to inside tube.

Use only the measured dimension and not the seat width to check compatibility.

(2) According to ISO 7176 - 4, depending on the ground surface, the weight of the user, driving behaviour, temperature and terrain.

(3) Depending on product configuration.

Base €		
STE070092	12 Inch wheel, pneumatic hybrid tyre, 250W/630W	STD
STE070093	16 Inch wheel, pneumatic hybrid tyre, 250W/750W	XXX,-
Ordering/Bundling €		
STE090213	Without wheelchair order	STD
STE090217	With wheelchair order, not mounted to the wheelchair <sup>(1)</sup>	–
Docking / One-Touch-Docking-Clamps €		
One-Touch-Docking-Clamps		
STE010053	One Touch Docking Clamps for round tube frames, adjustable diameter 19mm-24mm <sup>(2)(3)</sup>	–
STE010054	One Touch Docking Clamps for round tube frames, adjustable diameter 25mm-35mm <sup>(2)(3)</sup>	STD
STE010059	One Touch Docking Clamps for oval tubes (Quickie Helium, Quickie Nitrum, Quickie Xenon <sup>2</sup> FF) <sup>(2)</sup>	XXX,-
Adjustability €		
Front Frame width (measured inside to inside tube) <sup>(4)</sup>		
STE000160	220mm - 275mm	–
STE000161	275mm - 330mm	STD
STE000162	330mm - 425mm	–
Throttle €		
Throttle lever		
STE000082	Right hand mount (Proportional electrical brake on opposit side)	STD
STE000083	Left hand mount (Proportional electrical brake on opposit side)	–
Handlebar €		
STE000163	Narrow 480mm	STD
STE000164	Wide 540mm	–
Brake €		
	Proportional electrical brake (on opposite side to the throttle lever)	STD
	Mechanical disc brake Shimano, incl. parking brake function	STD
STE060026	Brake lever right hand mount	STD
STE060027	Brake lever left hand mount	–
STE060025	Mechanical double disc brake	XXX,-
Speed Setting €		
STE000025	Speed, max. 6 kp/h	STD
STE000027	Speed, max. 10 kp/h <sup>(5)</sup>	XXX,-
STE000028	Speed, max. 15 kp/h <sup>(5)</sup>	XXX,-
Display €		
	Display: Acceleration Setting, Speedometer, Average Speed, Maximum Speed, Trip Time, Trip Distance, Total Kilometers, State of Charge, Power Consumption	

(1) Wheelchair order on separate order form. Separate shipment of wheelchair and StreetJet.

(2) Always check installation space for the docking clamps before ordering. Depending on configuration and option, a particular configuration of a chair might not offer enough space.

(3) Always check combination agreements for non-Sunrise wheelchairs.

(4) Frame width measured from inside frame tube to inside frame tube at the mounting position.

Use only the measured dimension and not the seat width to check compatibility.

(5) Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via Fax with a signed waiver (next page). Max. speed-settings can only be altered by Sunrise Medical.

## Colours €

Frame Colour (Fork and main frame)		
RGK100200	Black matt 105	STD
RGK100171	White matt 101	XXX,-
RGK100202	Grey matt 131	XXX,-
RGK100203	Light grey matt 133	XXX,-
RGK100204	Navy blue matt 152	XXX,-
RGK100151	Black gloss 103	XXX,-
RGK100172	White gloss 100	XXX,-
RGK100152	Grey gloss 132	XXX,-
RGK100170	Cream gloss 135	XXX,-
RGK100166	Red gloss 108	XXX,-
RGK100169	Yellow gloss 120	XXX,-
RGK100163	Forest green gloss 143	XXX,-
RGK100161	Green gloss 107	XXX,-
RGK100159	Dark blue gloss 151	XXX,-
RGK100158	Royal blue gloss 154	XXX,-
RGK100156	Blue gloss 174	XXX,-
RGK100160	Duckegg blue gloss 141	XXX,-
RGK100168	Pink gloss 188	XXX,-
RGK100167	Orange gloss 180	XXX,-
RGK100157	Purple gloss 157	XXX,-

### Color Kit anodized (stand, central section, One-Touch-Clamps)

STE100146	Anodized, black	STD,-
STE100147	Anodized, red	XXX,-
STE100148	Anodized, blue	XXX,-
STE100149	Anodized, purple	XXX,-
STE100173	Anodized, gold	XXX,-

### Cable colour (brake cable, docking release cable)

STE060028	Cables, black	STD,-
STE060029	Cables, blue	XXX,-
STE060030	Cables, red	XXX,-

## Accessories €

	LED light	STD
	Carbon fender	STD
	Stand, with wheels	STD
STE090087	Additional Battery, Lithium-ion, 14 Ah - 36V	XXX,-
STE090189	Luggage rack (in frame colour)	XXX,-
	Front Frame Protector round tube (only in combination with STE010054 One Touch Docking Clamps for round tube frames, adjustable diameter 25mm-35mm) <sup>(3)/(4)</sup>	
STE090190	Frame protector, round tube 19mm	XXX,-
STE090191	Frame protector, round tube 25mm	XXX,-
STE090192	Frame protector, round tube 30mm	XXX,-
STE090194	Frame protector, round tube 32mm	XXX,-
STE090195	Frame protector, oval tube (Nitrum, Helium, Xenon) (only in combination with STE010059 One Touch Docking Clamps for oval tubes (Quickie Helium, Quickie Nitrum, Quickie Xenon² FF) <sup>(3)/(4)</sup>	XXX,-

## Empulse StreetJet



// STE070092 12 Inch wheel, pneumatic hybrid tyre, 250W/630W



// STE070093 16 Inch wheel, pneumatic hybrid tyre, 250W/750W

## Docking-clamps



// STE010059 One Touch Docking Clamps for oval tubes (Quickie Helium, Quickie Nitrum, Quickie Xenon<sup>2</sup> FF) <sup>(2)</sup>

## Accessories



LED light, Carbon fender, disc brake



Stand, with wheels



// STE090189 Luggage rack (in frame colour)



Front Frame Protector (example STE090192). Frame protector, (round tube 30mm)

## Brake



// STE060025 Mechanical double disc brake

## Colours

### Frame Colour (Fork and main frame)



### Color Kit anodized (stand, central section, One-Touch-Clamps)



### Cable colour (brake cable, docking release cable)



**Letter of Acknowledgement**

We herewith advise, that using the StreetJet with max. speed-settings higher than 6 kph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the StreetJet with speed-settings > 6 kph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the StreetJet speed settings higher than 6 kph (10/15 kph). When you're travelling with your StreetJet, always check the specific legislation in the respective country.

When placing an order for a max. speed-setting > 6 kph, please confirm with your signature below, that you have understood the above information and that you only intend the usage of your StreetJet on private property.

**I herewith confirm, that I am aware of the information above, that I will follow strictly the country specific legislation.**

**Date****Customer Signature**

Day

Month

Year

## Comments

## StreetJet

€

<b>Base Price</b>	<b>00,000.00</b>
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
<b>Total Price</b>	

## Date

## Signature

Day      Month      Year

**Sunrise Medical GmbH**

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu