

up! go!

C-ME

zippie®
Explore your world!

QUICKIE®
Live without limits.

Q Series
Powerchairs



Q300



Q400



Q500



Q700

YOUR CHAIR IS DESIGNED TO SUPPORT YOU, NOT TO DEFINE YOU

C-Me is all about seeing you, and just you.

Are you worried that you are most of the time out of sight of those around you? Do you mind that you cannot reach things independently? Would you like to go for a walk with friends and family and talk at eye level? With C-Me, you can achieve the independence you're looking for.

Elevate up to 300 mm whilst travelling 4.5 km/h and move shoulder-to-shoulder with everyone around you.

The C-Me function allows you to participate more actively in life without relying on others and decreases the risk of shoulder pain and injury. It allows you to reach areas that were out of reach, facilitates communication with other people and thus increases your level of independence and self-confidence.



FOR MORE INDEPENDENCE AND BETTER MOBILITY

Drive shoulder to shoulder

Drive at walking speed in full lift when you are out with your friends - C-Me takes away frustrating speed restrictions and you don't have to strain yourself when reaching up or constantly look up to communicate with people.

Ride in Comfort

Our suspension ensures a smooth and comfortable ride from start to finish. You can even lift or descend at the touch of a button with no need to stop while transitioning either!

Makes transferring easy

Together with the seat elevator the C-Me technology supports you while transferring on to surfaces which are higher.



Access new heights

More safety in road traffic, especially when crossing busy streets, due to a better overview and increased access to difficult to reach areas like cupboards, shelves and doorbells.



Scan the QR code with your iPhone's camera to get more information.

AVAILABILITY OF C-ME

C-Me is available on all our four drive wheel types and on five different seating systems (Kids, Teens, SEDEO LITE, SEDEO PRO and SEDEO ERGO) of the QUICKIE Q-series products from Q300 up to Q700.

CLINICAL NEED



Q50

SEATING
SLING SEAT

BASES
RWD



Q100

SEATING
SLING SEAT

BASES
RWD



Q200

SEATING
SLING SEAT

BASES
RWD



Q300

SEATING
SEDEO PRO
SEDEO LITE
KIDS SEAT
TEENS SEAT

BASES
MWD
RWD



Q400

SEATING
SEDEO PRO
SEDEO LITE

BASES
MWD
FWD
RWD



Q500

SEATING
SEDEO PRO
SEDEO LITE

BASES
MWD
FWD
RWD
HWD



Q700

SEATING
SEDEO PRO
SEDEO PRO ADV.
SEDEO PRO for HD

BASES
MWD
FWD
RWD



Q700-UP

SEATING
SEDEO PRO ADV. UP

BASES
MWD
FWD

C-ME AVAILABILITY

BASE PERFORMANCE

☐ RWD: Rear-Wheel Drive
 ☐ MWD: Mid-Wheel Drive
 ☐ FWD: Front-Wheel Drive
 ☐ HWD: Hybrid-Wheel Drive

Not all seats can be combined with all bases. Please refer to the order form for possible seat/base combinations.



Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.de
www.SunriseMedical.eu

Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines.

In line with our policy of continual product improvement, Sunrise Medical reserves the right to make changes at any time, and without notice. For the latest information, please visit www.sunrisemedical.co.uk

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu



Rev.A04/2023