



ATTITUDE Power



CE ORDER FORM

Rev: 05/2025 Valid:10/2025

| Billing Address Or Customer-ID | Shipping Address Identical To Billing Addres | | | |
|--|--|--|--|--|
| CID: | CID: | | | |
| Name: | Name: | | | |
| Street: | Street: | | | |
| Town/PC: | Town/PC: | | | |
| E-Mail: | E-Mail: | | | |
| Phone: | Phone: | | | |
| Order Information | User Weight | | | |
| Order Calculation Built-4-Me® | Max. user weight: 100 kg | | | |
| Purchase Order No.: | User Weight: kg | | | |
| Art. No.: 47500000 / 47500001 IC Reimbursement Code: xxx | | | | |
| Order Date: Day Month Year | Base Price € XXXX- | | | |

| Attitude - Mo | del Selection | € |
|---------------|--|-------|
| ATT010803 | ATTITUDE Power Speed 6 kph (standard); handlebar; disk brake; 9 Ah Lithium-ion battery front and rear lamp; docking system; speedometer integrated in control display; stand with wheels; parking brake. | XXX |
| ATT072111 | Hand throttle position right side | XXX,- |
| ATT072112 | Hand throttle position left side | - |

| Technical Data Overview | | | |
|--------------------------------------|----------------------------|-----------------|---------------------|
| Max. user weight | 100 kg | Engine | 250 W |
| Front frame width | 260 – 460 mm | Battery | 9 Ah oder 13,4 Ah |
| Weight of the attachment adapters | 85 kg | Handlebar width | 570 mm |
| Total weight | 15,3 kg | Brake | disc brake (180 mm) |
| Wheel size | 16" | Parking brake | yes |
| Maximum speed of the engine | 6 km/h (optional: 21 km/h) | Frame colors | 22 |

| nrise Wheelchairs | | Other Wheelchairs | |
|---|-------|---|----|
| Easy Life T (fixed front frame only) | - | Küschall K-serie; KSL Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s Fits Kueschall round-tube frames until build date 03/ 2020 | X> |
| Easy Life RT (fixed front frame only) | - | Küschall Advance, Champion Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's Fits Kueschall round-tube frames until build date 03/ 2020 | X |
| Easy Life (fixed front frame only) | - | Otto Bock Avantgarde CS/CLT Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s Only possible for frames without inset | X |
| Easy Life R (fixed front frame only) | - | Otto Bock Motus CS Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s Only possible for frames without inset | Х |
| Easy 200 | - | Meyra Avanti Pro Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | Х |
| Easy Max (fixed front frame only) | - | Meyra Hurricane Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | Х |
| Neon (fixed front frame only) | - | Meyra Nano Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | Х |
| Xenon/Xenon ² (fixed front frame only) | - | Meyra Smart F Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's | Х |
| Argon ² | - | PRO ACTIV Traveler Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | Х |
| Helium Helium with handbike axle is only approved with frame reinforecment kit for add-on bikes | - | PRO ACTIV Speedy Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | Х |
| Nitrum Nitrum with handbike axle is only approved with frame reinforecment kit for add-on bikes | - | Progeo Joker Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | Х |
| RGK Tiga Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo´s | XXX,- | TiLite TR, TRA | Х |
| QS5 X (fixed front frame only) | - | TiLite Aero T, Z, ZR, ZRA | X |
| | | Panthera S2, S3, U2, U3, Build 4 me - different lead-times and discount structures | X |
| | | Panthera U2 Light, U3 Light Build 4 me - different lead-times and discount structures | XX |

B4ME-Options



// To move the pedal position 2 cm closer or up to 6 cm further

| ATT090213 | ATTITUDE only - without a related wheelchair order | | |
|---|--|--|--|
| ATT090214 | ATTITUDE with wheelchair order - on a seperation | ate wheelchair order form | |
| FMG- MA03-3201- 01 | B4Me option: Adaptation link -2 to + 6 cm Build 4 me - different lead-times and discount structures | | |
| ieelchair | fornt frame width - measured from outsid | e front frame tube to outside front frame tube | |
| 260 mm | | 380 mm | |
| 280 mm | | 400 mm | |
| 300 mm | | 420 mm | |
| 320 mm | | 440 mm | |
| 340 mm | | 460mm | |
| 360 mm | | | |
| k O | | | |
| | | | |
| ATT080034 | Fork, Standard | | |
| | Fork, Standard | | |
| ATT080034 | Fork, Standard Tyre, pneumatic, tread, Mountainbike | | |
| ATT080034 | | | |
| ATT080034 es (C) ATT070110 ATT070111 | Tyre, pneumatic, tread, Mountainbike | able stem standard 🔘 | |
| ATT080034 es (C) ATT070110 ATT070111 | Tyre, pneumatic, tread, Mountainbike Tyre, pneumatic, tread, Big Apple | able stem standard 🔘 | |

| ATT060050 | Brake, manual actuation not with gear hub, 8 speed, coaster brake | | |
|-----------|---|-------------------------------------|---|
| ATT060053 | Performance brake lever (offers improved adjustability & modulation of the brake) | | |
| | ATT060102 | Brake lever mounted to right handle | - |
| | ATT060101 | Brake lever mounted to left handle | _ |

Tyres



// Tyre, pneumatic, tread, Mountainbike



// Tyre, pneumatic, tread, Big Apple



// S'em, multiadjustable



// Performance brake lever

| Speed only for | countries which allow a higher speed or for use on private ground only. This option can only be ordered 🔟 | € |
|-----------------|---|----------------|
| ATT000025 | Speed, max. 6 kph | - |
| ATT000027 | Speed, max. 10 kph | - |
| ATT000028 | Speed, max. 15 kph | - |
| Battery - charg | er always included 🔘 | € |
| ATT160021 | Battery, Lithium-ion, 9 Ah | XXX,- |
| ATT160022 | Battery, Lithium-ion, 13,4 Ah | XXX,- |
| ATT090086 | Additional battery, Lithium-ion, 9 Ah | XXX,- |
| ATT090087 | Additional battery, Lithium-ion, 13,4 Ah | XXX,- |
| ATT090089 | Additional charger, Power | XXX,- |
| Accessories | ۵ | € |
| ATT090074 | Front- & rear lamp | - |
| ATT090080 | Speedometer | - |
| ATT090095 | Mudguard | XXX,- |
| ATT090067 | Rear view mirror | XXX,- |
| ATT090081 | Bottle holder | XXX,- |
| ATT090082 | Low rider carrier for saddle bags | XXX,- |
| ATT090024 | High pressure pump (0 - 11 bar) | XXX,- |
| | | |
| ATT090084 | Additional weight, removable, 5 kg, to improve traction | XXX,- |
| | Additional weight, removable, 5 kg, to improve traction Bell | XXX,- XXX,- |
| ATT090084 | | , |

Chain guard Colours You can choose separate colours for the frame, the cranks and the chain guard. Frame Chain guard Crank Frame Chain guard Crank 01 Matte Black 12 Royal Blue 02 Gloss Black 13 Brilliant Ocean 03 Metallic **14 Scarlet Noir** Shadow 04 Starlight White **15 Red Pepper** 05 Silver Ore **16 Cherry Lustre** 06 Storm Gray 17 Turbo Green 07 Desert Tan **18 Pink Power 08 Desert Sage 19 Satin Rose** 20 Mystic Purple **09 Salute Green 10 Velvet Sapphire** 21 Solar Orange 11 Midnight Blue 22 Velocity Yellow

| Sticker colours | | |
|-----------------|--------|---|
| ATT100550 | Black | - |
| ATT100551 | Silver | - |

Accessories



// Mudguard, bolt on, black



// Multitool



// Rear view mirror



// Bottle holder





// Ccnnect Display ccnnects with your smartphone (download of neodrives app necessary) Records trip data such as average speed, maximum speed, trip distance, and total distance.



// Additional weight, removable, 5 kg, to increase traction



// High pressure pump (0 - 11 bar)



// Bell



// Low rider carrier for saddle bags

Letter of Acknowledgement

We herewith advise, that using the ATTITUDE Power with max. speed-settings higher than 6 kph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the ATTITUDE Power with speed-settings > 6 kph may not be driven on public roads, bicycle-lanes, or pavement. You need to inform yourself about the country specific legislation regarding usage and your legal obligations regarding the usage of the ATTITUDE Power with speed settings higher than 6 kph (10/15 kph). When you're travelling with your ATTITUDE Power, always check the specific legislation in the respective country.

When placing an order for a max. speed-setting > 6 kph, please confirm with your signature below, that you have understood the above information and that you only intend the usage of your ATTITUDE Power on private property.

I herewith confirm, that I am aware of the information above, and that I will strictly follow the country specific legislation.

| Date | | | Customer Signature |
|------|-------|------|--------------------|
| Day | Month | Year | |

Comments

Please find some examples for already existing B4me options below:

FMG-MJ02-3187-02 - Tetra operation (ATTITUDE Hybrid) FMG-MA03-3201-01 Adaptation link -2 to + 6 cm (ATTITUDE Manual, Hybrid & Power)



| ATTITUDE Manual | € |
|-------------------------------|-----------|
| Base Price | 00,000.00 |
| Total Of Additional Equipment | |
| Packaging And Transport | |
| Exclusive VAT | |
| VAT | |
| Total Price | |

| Date | | | Signature |
|------|-------|------|-----------|
| Day | Month | Year | |



Sunrise Medical GmbH

Kahlbachring 2-4 D-69254 Malsch/HD Phone: +49 (0) 72 53 / 9 80-0 Fax: +49 (0) 72 53 / 9 80-222 E-Mail: export@sunrisemedical.de www.SunriseMedical.eu