

Q300[®] SERIES

Compact in size. BIG on performance.



Q300 M Mini Kids



Q300 M Mini Teens



Q300 M Mini with SEDEO LITE



Q300 M Mini with SEDEO PRO

Q300M Mini

One of the narrowest
TRUE mid-wheel drive
power wheelchairs.

ULTRA-RELIABLE,
HIGH TORQUE 4-POLE MOTORS
WITH SPEEDS UP TO 10 KM/H

OCCUPIED TRANSIT
TESTED TO ISO-7176



SOLID SHEET METAL CHASSIS

SIX WHEEL DRIVE SUSPENSION

CLIMB CURBS
UP TO 70 MM

RELIABLE

With a solid sheet metal chassis made, the industry best motors and end-of-line testing, every Q300 M Mini was built to keep you moving.

TINY TURNING

Thanks to its TRUE mid-wheel drive technology, the Q300 M Mini turns on the spot for an ultra-small 1120 mm turning circle. Other compact MWD chairs can have offset drive wheels, increasing their overall turning radius.

ULTRA-SLIM

Ridiculously narrow at 520 mm wide, the Q300 M Mini goes where others can't. Whether a complicated living room, a cluttered kitchen, or a tight doorway, it'll fit into your life, not the other way around.

ELEVATE AND DRIVE WITH C-ME®

C-Me enables you to elevate up to 300 mm while travelling up to 4.8 km/h! You can lift or descend at the touch of a button with no need to stop while transitioning. Our suspension doesn't lock out when you drive with C-Me, giving you a smooth and comfortable ride from start to finish.





SEDEO® LITE™

Truly modular in design, this seating system provides a new level of functionality and support for mild to moderate postural needs. The integrated seat rail can mount additional positioning elements for further postural needs.



SEDEO® PRO™

A truly versatile seating system, SEDEO PRO combines a highly adjustable seat frame with an excellent choice of seating to accommodate a wide range of moderate to complex needs. SEDEO PRO can be easily and quickly adapted to the changing needs of an individual. The integrated seat rail can mount additional positioning elements for further postural needs.



KIDS SEAT

This seating system provides a power chair that can keep up with growing children, offering 280 - 400 mm width and 260 - 400 mm depth adjustment. The contoured SEDEO LITE backrest provides comfort and stability. As the lateral support is situated inside the backrest, it won't shift or slip over time. Great for kids on the move!



TEENS SEAT

Packed with all the growth adjustability and features of the Kids version, the Q300 M Mini Teens stretches to an extended seat depth range of 360 - 500 mm. Ideal for growing teenagers who require a narrow seat width but larger adjustments to seat depth over time.

Q300 M MINI SEAT FRAME COMPARISON

	SEDEO LITE	SEDEO PRO	KIDS SEAT	TEENS SEAT
Seat Width (mm)	385 - 510	380 - 560	280 - 400	300 - 400
Seat Depth (mm)	410 - 560	400 - 560	260 - 400	360 - 500
Seat height (mm)	410 - 490	410 - 490	410 - 490	410 - 490
Fixed angle tilt	0 - 9°	0 - 9°	0 - 9°	0 - 9°
Backrest height (mm)	460 - 510	530 - 720	380 - 485	460 - 510
Backrest angle	85 - 120°	90 - 150°	86° - 120°	90 - 120°
Maximum user weight	136 kg	136 kg	76 kg	120 kg
Elevate with C-Me (mm)	300	300	300	300
SEDEO PRO Flexible Back Cushion	●	●		
JAY & WHITMYER	●	●	●	●



Additional information, videos and selection guide on the seating options.

Q300 M Mini

See order form as specifications may vary based upon configuration

Total width	From 520 mm
Total length (with centre-mount)	From 1060 mm
Turning circle	1120 mm
Maximum safe slope	8°
Maximum kerb climbing	50 mm (with 12" drive wheels) 60 mm (with 13" drive wheels) 70 mm (with 14" drive wheels)/h
Maximum range (ISO 7176-4)	Up to 25 km (with 41 Ah batteries) Up to 36 km (with 56 Ah batteries)
Product weight (with batteries)	From 103 kg
Drive wheels	12" / 13" / 14"
Max. speed	6 km/h / 8 km/h / 10 km/h
Battery size	41 Ah / 50 Ah / 56 Ah
Electronics	VR2, R-NET
Transport without user	Fix the wheelchair safe in the vehicle (e.g. car or plane).
Intended use & environment	For people with limited mobility, indoor and outdoor use. Wheelchair type class B
Transit Brackets	ISO 7176-19:2008 Transit Approved

* Exact measurements may vary depending on options.



Order Form



Web Page

COLOR OPTIONS

AVAILABLE SHROUD COLORS



AVAILABLE ACCENT COLORS



Important Note—
Print processes cannot accurately represent painted colors. Actual colors may vary.



Sunrise Medical GmbH
Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.de
www.SunriseMedical.eu

Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines.

In line with our policy of continual product improvement, Sunrise Medical reserves the right to make changes at any time, and without notice. For the latest information, please visit www.sunrisemedical.eu
If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu