















GEMINO 30 GEMINO 30 Carbon GEMINO 30 Comfort

GEMINO 60 GEMINO 30 Walker

r GEMINO 60 Walker

GEMINO 30 SpeedControl GEMINO 30 Parkinson

Built to perform indoors and out, GEMINO lightweight rollators set a new standard in comfort, safety and design. Made to stand out from the crowd, every rollator is infused with a style and elegance to suit every individual taste. GEMINO is a symbol of independence, and will be with you for every journey.

Lightweight Rollators For An Active Lifestyle

Are you looking for a rollator that fits your active lifestyle and allows you to walk comfortably and confidently, indoors and out? Then a GEMINO rollator is the perfect choice for you. Our awardwinning lightweight rollators set new standards in quality, comfort, safety and design.

The GEMINO comes in various models, sizes and colours. Naturally, there's also a special/medical GEMINO rollator for you, if you need something extra to walk safely. So there's always a GEMINO that fits your needs and personal style.





Creating Enjoyable Walking Experiences

GEMINO rollators have been designed with your comfort in mind. At the heart of this is the lightweight aluminium frame that offers maximum stability, light braking and easy movability, ensuring a great walking experience indoors and out.

Walk With Confidence

Wherever you go, you'll feel confident with a GEMINO rollator. The sturdy, stable frame design ensures maximum safety and security in all settings.





Easy Braking

Ergonomically designed, you'll find that operating the brake handles requires very little effort. The smart brake design is easy on the wheels and ensures optimal grip.






Memory Lock Push Handles

The ergonomic handles are always in the right position thanks to the easy-to-use height adjustment system and memory lock function, so once the handles are at the right height for you, they stay like that.

Don't Let Obstacles Get In Your Way

Kick tubes make it easy to climb kerbs or thresholds, while the protective side bumpers help you avoid obstacles.









Designed for comfortable walking indoors-

and-outdoors, even on rough terrain. Large,

soft wheels allow overcoming of obstacles

Gemino 30

23

GEMINO30 **GEMINO**30 Comfort

Award-winning design combining high indoor manoeuvrability with outdoor performance. Lots of colour options and accessories for easy personalisation.

Everyday use

150 kg 130 kg (M)

125 kg (S)

150-200 cm 135-170 cm (M)

125-165 cm (S) Also suitable for children

Carbon

with ease. Use indoors and out, mainly in the city

Indoor and outdoor use, also

on uneven terrain

Everday use / Active use

150 kg 130 kg (M)

> 150-200 cm 135-170 cm (M)

Large, soft wheels ensure ultimate comfort and easy overcoming of obstacles

Easy personalisation with 7 frame colours and 8 bag colours. For GEMINO 30 Carbon only Carbon black available.

Small turning radius for great

indoor manoeuvrability


Super lightweight 5.4 (M)/5.6 kg Carbon

version available for easy transportation and lifting

Extra strong bag (10 kg load capacity)



GEMINO 60



GEMINO30

GEMINO30



GEMINO30

Specially designed for people with Parkinon's. Featuring a slowdown brake promoting a steady pace, and a reverse braking system and laser

GEMINO 30 Enjoy Life To The Fullest

Looking for a lightweight rollator that stands out? The GEMINO 30 range of rollators exceed all expectations when it comes to comfort, safety and design. The modern look of the GEMINO 30 is complemented with the highest quality materials to create a stylish rollator that will fit any lifestyle. Depending on your height, we also offer the GEMINO 30 M and GEMINO 30 S so you can walk with confidence.

Easy Folding With Self-Locking

The GEMINO 30's unique 'fold-click system' makes it very easy to fold to a compact size for storing. It is practical for small spaces and portable for transportation and, you can even safely manoeuvre the GEMINO 30 in tight places when folded.

Intergrated Net Bag

Got a lot to carry? The net bag can be easily opened with one hand, offers plenty of space for purchases and carries up to 5kg.

Anatomical Shaped Handgrips

Ergonomically designed hand grips provide extra comfort and support while walking.

Make A Statement

Looking for individuality? Express your personality with a selection of 7 anodised frame colours.





GEMINO 30 Carbon The Premium Carbon Rollator For An Active Lifestyle

Looking for something lighter? The GEMINO 30 Carbon and GEMINO 30 M Carbon, have all the advantages of the GEMINO 30 but with a super lightweight and strong carbon fibre frame, offering support and making it effortless use as part of any active lifestyle.

Super Lightweight

The GEMINO 30 Carbon includes the same benefits as the GEMINO 30 with a weight from 5.4 kg making the GEMINO 30 Carbon easier to lift, transport and store, while also being super manoeuvrable in tight spaces.



GEMINO 30 Comfort High Indoor Manoeuvrability With Outdoor Performance

Are you looking for a lightweight rollator that allows you to walk comfortably and confidently, indoors-and-outdoors? Then the GEMINO 30 Comfort rollator is the perfect choice for you. With large soft wheels and superb rolling performance for extra comfort.

Wheels That Make The Difference:

The large, soft, shock absorbing wheels allow you to walk comfortably on any surface.

Indoor And Outdoor Use: The compact frame gives you high indoor manouevrability and the large soft wheels allow you to overcome obstacles with ease whether it is indoor tresholds or outdoor kerbs.

Handle Extra Weight: Lots of things to carry? GEMINO 30 Comfort can hold upto 10 kg in its basket.



GEMINO 60 Ultimate Comfort & Safety On Any Terrain

Whether strolling the highstreet or enjoying a walk in countryside, the GEMINO 60 is the perfect companion for your active lifestyle. With its extra-large soft wheels, this rollator allows you to walk comfortably and securely, even on rough terrain.

Designed For Superior Stability:

The GEMINO 60's sturdy, extra wide frame combined with the large wheels ensures excellent stability and safety, even on rough terrain.

Enjoy More Space: Designed for maximum stability, whilst offering extra walking and seat space.

Heavy-Duty:

All GEMINO 60 models are suitable for heavy-duty use and take a maximum user weight of 150 kg.



GEMINO 30 Walker Extra Support & Stability

If you need extra support when walking, the GEMINO 30 Walker is the perfect choice for you. Combining the benefits of the GEMINO 30 with enhanced support and stability makes the GEMINO 30 Walker ideal for gait training.

Forearm Supports For Enhanced Stability And Safety

You can easily adjust the height, depth and angle of the forearm supports to suit your needs. The brakes and push handles are fully integrated, easy to reach and light to operate, making walking safe and comfortable.

Advanced Braking Options

Looking for even more reassurance? Complement your GEMINO 30 Walker with a slowdown brake or innovative SpeedControl braking system. With SpeedControl your rollator will brake automatically when you start to accelerate.



GEMINO 60 Walker Optimal Support & Stability

With the GEMINO 60 Walker get the best of both worlds. You get the comfortable rolling performace and stability of the GEMINO 60 with additional support and safety from a walker. The GEMINO 60 Walker is great for gait training.

Comfortable And Safe Outdoor Use:

The GEMINO 60 Walker includes the same benefits as the GEMINO 30 Walker combined with the sturdy frame design with extra-large soft wheels, allowing you to go out in total comfort and safety.



GEMINO 30 SpeedControl The Latest Innovation

With the GEMINO 30 innovative SpeedControl dynamic braking system you'll feel completely secure while walking. SpeedControl helps prevent falls and is also ideal for gait training.

Dynamic Braking For Maximum Security

SpeedControl promotes a steady pace, ensuring that you walk comfortably. When you start to accelerate the brakes in the rear wheels will automatically activate, preventing the rollator from 'running away' from you.

Adjustable To Your Needs

The SpeedControl system features 3 levels of braking that can be easily adjusted to your needs.





GEMINO 30 Parkinson Specially Designed To Support Those With Parkinson's Disease

Our GEMINO 30 Parkinson rollator has been developed people with Parkinson's Disease in mind and offers features that are fully tailored for particular needs.

Slowdown Brake Promotes A Steady Pace

The GEMINO 30 Parkinson is fitted with an adjustable slowdown brake which provides additional resistance and helps promote walking at a steady pace.

Increased Safety With Reverse Brake

The reverse braking system immediately stops moving once the brake is released. This increases safety if you experience 'freezing' while walking.

Laser Device Helps Overcome Freezing

If 'freezing' occurs you can activate a laser line with the push of a button. This can help you get moving again by concentrating on stepping over the laser line as it appears on the ground.



Accessories

GEMINO rollators are not just easily adjustable, they also come with a full range of attractive accessories to make everyday life easier. A stylish bag, a practical crutch holder or a light for extra safety. Simply choose the options that fit your needs and lifestyle.



Umbrella^{*}

With this height-adjustable umbrella you can stay dry while keeping both hands on the push handles.



LED Light

Stay visible with this LED light that attaches to your frame.



Table Tray

A practical tray with cup holder, allows you to transport food and beverages safely with ease.



Crutch Holder

Keep your crutch or walking cane securely fastened when walking. Easy to mount, easy to access.



Universal Bracket

Suitable for fitting your choice of personal accessories.



Adjustable Back Rest**

Provides support to your back to allow you to sit more comfortably while resting.



Slowdown Brake*

Provides additional resistance to help you walk at a steady pace. Easy to adjust to the desired level – and easy to turn on/off. Not available with SpeedControl.



SpeedControl***


Brake automatically when you start to accelerate (such as going down hill), ensuring smooth walking and a comfortable, upright position at all times. Not available with slowdown brake. Available on all GEMINO 30 models.

* not available for GEMINO 30 Carbon / 30 M Carbon ** not available for GEMINO 30 / 30 M Walker *** check availability per model





Graphite Grey






Pink



Green



Alcantara Black







Mint Green


Colourful Advanced Bag

Enhance your GEMINO with a stylish, closable bag in your favourite colour. With handles and shoulder strap for easy carrying.

Blue



Technical Data








GEMINO30 **GEMINO**30 GEMINO 30 Carbon Comfort BOM 305 30 30**M** 30 not Carbon Maximum user weigh 125 kg 150 150 Recommended user height cm 150-200 135-170 125-165 150-200 135-17 Maximum load capacity in basket kα Total length 🗛 cm Total width **B** cm 60 55 62 62 7.0 (7.3) 7.1 (7.4) Total weight without basket (with basket) 7.8 (7.5 kg 5.6 (5.9) Carbon 5.4 (5.7) Carbon 6.6 (6.8) 8 (7.7) Width between centre of push handles C cm 42 47 Seat height D cm 62 55 6/ 69 - 88 Total height 🕒 cm 73 - 88 Carbon 78 -65 - 77 80. 71 - 9 Total height folded cm Total lenght folded 68 cm Total width folded 23 23 23 23 cm Diameter front wheels 20 x 3.6 cm 20 x 3.6 20 x 3.6 23 x 4.2 23 x 4.2 Diameter rear wheels **G** 20 x 3.6 23 x 4.2 23 x 4.2 cm 20 x 3.6 20 x 3.6 PERFORMANCE Turning radius cm 8/ Black Silver grey √ not Carbon √ not Carbor Champagne √ not Carbon √ not Carbo Sapphire blu √ not Carbon √ not Carbo Pink √ not Carbon Apple green √ not Carbon √ not Carbo Scarlett red √ not Carbon CE EN 12182 ISO 11199-2 ISO 11199-3

GEMINO Walker rollators: EN 12182, ISO 11199-3 - All other GEMINO rollators: EN 12182, ISO 11199-2













		Gemino 60		GEMINO 30 Walker		GEMINO 50 Walker		GEMINO 30 SpeedControl			GEMINO 30 Parkinson	
		50	60 M	IO	BOM	50	60 M	IO	EOM	eos	30	Bom
DIMENSIONS												
Maximum user weight	kg	150	150	150	130	130	130	100	100	100	130	130
Recommended user height	cm	150-200	135-170	150-200	135-170	150-200	135-170	150-200	135-170	125-165	150-200	135-170
Maximum load capacity in basket	kg	10	10	5	5	10	10	5	5	5	5	5
Total length \Lambda	cm	74	74	65	65	74	74	65	65	68	65	65
Total width	cm	64	64	60	60	64	64	60	60	55	60	60
Total weight without basket (with basket)	kg	8.5 (8.8)	8.2 (8.5)	10.6 (10.9)	10.4 (10.7)	11.8 (12.1)	11.6 (11.9)	8.7 (8.9)	8.6 (8.8)	8.2 (8.4)	7.5 (7.8)	7.4 (7.7)
Width between centre of push handles 🕑	cm	51	51	41	41	45	45	47	47	42	47	47
Seat height D	cm	62	53	62	55	62	53	62	55	48	62	55
Total height 🕒	cm	76-95	73 - 83	98 - 115	84 - 100	98 - 115	89 - 98	78 - 100	69 - 88	65 - 77	78 - 100	69 - 88
Total height folded	cm	84	80	97	89	109	101	80	76	69	80	76
Total lenght folded	cm	74	74	65	65	74	74	65	65	68	65	65
Total width folded	cm	23	23	35.5	35.5	35.5	35.5	23	23	23	23	23
Diameter front wheels 🕫	cm	28 x 4.2	28 x 4.2	20 x 3.6	20 x 3.6	28 x 4.2	28 x 4.2	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6
Diameter rear wheels G	cm	23 x 4.2	23 x 4.2	20 x 3.6	20 x 3.6	23 x 4.2	23 x 4.2	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6
PERFORMANCE												
Turning radius	cm	95	95	84	84	95	95	84	84	79	84	84
COLOUR OPTIONS												
Black												
Silver grey		\checkmark		\checkmark	\checkmark	\checkmark	\checkmark	\checkmark		\checkmark	\checkmark	
Champagne												
Sapphire blue										\checkmark	\checkmark	
Pink										\checkmark		
Apple green												
Scarlett red												
CERTIFICATES			,		1						1	,
CE		√	√	√	V	√ ,	√	√	√	√	√	V
EN 12182		√	√		\checkmark		\checkmark	√	√	√	√	√
ISO 11199-2		V						V		V	V	V
ISO 11199-3												

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

For more information on any of the GEMINO products in this cataloge, vist our website by scanning the QR code with your phone.



Gemino



IMPROVING PEOPLE'S LIVES® Thorns Road, Brierley Hill, West Midlands, DY5 2LD t: +44 (0) 845 605 6688 | e: marketing@sunmed.co.uk sunrisemedical.co.uk

Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines. In line with our policy of continual product improvement, Sunrise Medical reserves the right to make changes at any time, and without notice. For the latest information, please visit www.sunrisemedical.co.uk. If you are visually impaired, this document can be viewed in PDF format at www.sunrisemedical.co.uk

Gemino

771111