

### **TWO IN ONE!** VERSATILITY IS A BREEZE WITH THE P100

The BREEZY P100 is a folding powered wheelchair suitable for indoor and outdoor use. It offers quick assembly and disassembly, and is easily converted to a manual chair with the simple addition of self propelling 24" wheels included in the package.

## FROM POWER.....





**POWERED FUNCTION** Simply attach the portable motor unit under the base of the chair and you're ready to zoom around town.....

# .....TO MANUAL



MANUAL FUNCTION ....OR remove the powered unit and simply attach the 24" Quick release PU Wheel provided.

# **IN TWO EASY STEPS!**





- HYBRID both a power & manual chair in one
- AUSTRALIAN STANDARDS tested
- ELECTRONICS Easy to use VR2 electronics
- LITHIUM ION BATTERIES for longer life expectancy
- COLOURS available in Red
- VERSATILE for indoor & outdoor use
- FOLDING convertable power wheelchair
- BACK POUCH for convenient storage
- ADJUSTABLE seat height

COMPACT & FOLDABLE





**FEATURES** 





DIMENSIONS

							<u>H</u>	
CODE	Seat Width	Frame Width	Seat Depth	Folding Width	Armrest Height	Back Height	Frame Depth	Adj Seat Height
BRP100-16	400mm (16")	520mm (20")	398mm (16")	300mm (12")	210mm (8")	435mm (17")	1040mm (41")	880 – 930mm (34.5" - 36")
BRP100-18	450mm (18")	570mm (22")	398mm (16")	300mm (12")	210mm (8")	435mm (17")	1040mm (41")	880– 930mm (34 5″ - 36″)

#### **SPECIFICATIONS**

Front wheels	200*50 PU Tyre	Max. Groove Crossing Ability	1000mm	
Rear wheels	2.50 – 6 pneumatic tyre / Max. 50PSI (3.5 BAR)	Max. Speed Forward	7.8km / hr	
Seat Width	450mm (18"), 400mm (16")	Max. Speed in Low Speed Mode	2km / hr	
Motor Output	200W*DC24V*2PCS	Max Weight Capacity	120kg	
Brakes	Magnetic Brake	Weight (without batteries)	33.4kg	
Battery	20AH * 24V * 1PCS (Lithium Battery)	Front Frame Weight	15.4kg	
Range	25km (on the flat road, 120kg load)	Rear Frame Weight	18kg	
Slope Grade Durability	10° max.	Battery Weight	5.4kg	
Min Turning Radius	Approx. 870mm	Anti-tip Wheel	Yes (on power unit only)	
Max. Kerb Climbing Ability	40mm	Colour	Red	



©Sunrise Medical. Breezy P100 Rev A.AU

Sunrise Medical Australia Pty Ltd 11 Daniel Street Wetherill Park Phone: (02) 9678 6600 Email: enquiries@sunrisemedical.com.au www.sunrisemedical.com.au