

CERTIFICATE AIR TRANSPORTATION Li-ION BATTERIES 288 Wh FOR Empulse StreetJet, F55, and Attitude Hybrid USERS

- **Sunrise Medical GmbH** hereby certifies that our battery 288 Wh (8,0 Ah) Lithium Ion battery has successfully passed the United Nations Manual of Tests and Criteria, Part III, subsection 38.3 test.
- Therefore, **Sunrise Medical GmbH** certifies that our battery 288 Wh (8,0 Ah) Lithium Ion battery has a rated energy of 288 Wh, so the battery may be carried-on passenger aircraft without restriction according to IATA 2.3.2.3. and 2.3.2.4. limit of 300 Wh for powered mobility aids.

BATTERY DESCRIPTION

Battery Name	E-bike Battery Li-ion	Model	ANS-AEI-DHF1U-PS102B.40T
Voltage	36 V	Rated Energy	288 Wh
Rated Capacity	8,0 Ah	Dimension L x W x H	312,5mm x 91,9mm x 93,3mm
Battery weight	2,1 kg	Cell Weight	70 g

FIRST AID MEASURES

If materials from a ruptured or otherwise damaged cell or battery contact skin, flush immediately with water, wash affected area with soap and water. For eye contact, flush with copious amounts of water for 15 min. and see physician at once. Avoid inhaling any vented gases. If irritation persists, seek medical assistance.

REGULATORY INFORMATION

Regulations specifically applicable to the product:

- ADR (road transportation)
- IATA-DGR (air transportation)
- IMDG-Code (sea transportation)
- US Department of Transportation 49CFR - Code of Federal Regulations [USA]

TRANSPORT INFORMATION

The battery is specifically designed to be removed and so the user must follow this specific procedure:

- The battery must be removed from the StreetJet/F55/Attitude Hybrid and carried in the passenger cabin.
- In order to remove the battery properly, switch the battery off (on/off switch), turn the key to the OPEN position, grab the battery and gently pull to rider's right and up until the battery comes fully out of the holder.
- Once the battery is removed, turn the key to the OFF position and take off the key.
- The battery must be stored in such a way that prevents damage or switching on by the movement of the luggage or other cargo.
- The battery poles must be individually protected to prevent short circuits.
- The pilot-in-command must be informed of the location of the battery when removed and carried in the cabin.

This certificate does not apply to damaged or defective batteries. This certificate is valid meanwhile the DGR 65th ed. 2024 IATA regulation is in force. Due to possible amendments or modifications of the regulation it is crucial to check our website www.sunrisemedical.com where information about these modifications will be published when applicable.

- Before you travel, please contact the appropriate carrier. The travel operator will be able to supply details of any special requirements/instructions.



Daniel Haüter
Quality Director Europe

28-April-2024
Date

Sunrise Medical GmbH

Kahlbachring 2-4 · 69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: kundenservice@sunrisemedical.de
www.SunriseMedical.de



SUNRISE MEDICAL
Sunrise Medical GmbH
Kahlbachring 2-4
D-69254 Malsch/Heidelberg
Telefon +49 7253 980 0
Telefax +49 7253 980 111
www.SunriseMedical.de
Sitz: Malsch/Heidelberg
Handelsregister Mannheim HRB-Nr.: 729025
Gerichtsstand: Heidelberg
Geschäftsführer: Thomas Hardt, Kai Schyktanz, Udo Loosen

Bankverbindung:
HSBC Trinkhaus & Burkhardt AG Düsseldorf
IBAN: DE 03300308800670033007
BIC-Code: TUBDEDD