

Quickie® All Court Ti



All COURT Ti

The NEW Quickie All Court Ti basketball chair is extremely versatile and can be used for any position on any court. The new titanium frame reduces the overall weight of the chair while improving performance and turning responsiveness. The All Court Ti provides adjustments for specific types of disabilities, player positions and performance levels.

Options:

- Toe Clips
- Half Solid Seat Pan
- “Stand-Up” option for front bumper
- Quick-Release Anti-tip Wheel

Quickie All Court Specifications

Chair Weight	21.5 lbs (9.8 kg)
Seat Width	12" to 20" 1" increments
Seat Depth	10" to 20" 1" increments
Front Seat Height	16" to 21" 1/2" increments
Rear Seat Height	13" to 21" 1/2" increments
Backrest Heights	Fixed backrests 9", 10", 11", 12", 13", 14", 15", 16", 17" and 18" in 1" increments
Backrest Upholstery	Standard, sports adjustable
Casters	3" Micro, (New improved high performance) 3" Micro Lighted
Rear Wheels	24" X/S Spoke (540mm) 24" X/S Heat Treated (559mm) 26" X/S Heat Treated (590mm) 24" Spinerger (540mm) 24" Spinerger (559mm) 26" Spinerger (590mm)
Camber	16°, 20°
Hand Rim Options	Anodized, painted black, plastic coat, & titanium
Anti-tip Wheel Options	Quick-release, lighted
Maximum User Weight	250 lbs (120 kg)
Color Options	23 different colors
Side Guards	Plastic & aluminum
Positioning Belts	Hook & loop closure, nylon buckle closure, snowboard ratchet-style pos. belt (Lg & XLg.)
Chest Straps	2" & 4"
5th Wheel Support	Anti-tip brace (reduces rearward seat flex) Performs like a rigid chair, while maintaining adjustability

Specifications subject to change without notice.

New Features



Special anti-tip brace (optional)



Redesigned front bumper



Ratchet-style positioning belt (optional)



Extra-hard caster wheels



www.sunrisemedical.com

Sunrise Medical
7477 East Dry Creek Parkway
Longmont, CO 80503
303.218.4600
800.333.4000
800.263.3390 in Canada

Quickie® Match Point Ti



MATCH POINT Ti

The NEW Quickie® Match Point Ti Tennis chair offers a dual front-end caster design for a wider footprint and unsurpassed stability. The integral fifth wheel enhances maneuverability and stability. As with other Quickie sports chairs, the Quickie Match Point Ti features the patented center-of-mass adjustment that eliminates the need to re-square the casters and/or adjust the wheel heights when moving the center of gravity.

Options:

- Toe Clips
- Quick-release anti-tip wheel

Quickie Match Point Specifications

Chair Weight	20 lbs (9.1 kg)	
Seat Width	12" to 20"	1" increments,
Seat Depth	10" to 18"	1" increments
Front Seat Height	16" to 21" 1/2" increments	
Rear Seat Height	13" to 21" 1/2" increments	
Backrest Heights	Fixed backrests 9", 10", 11", 12", 13", 14", 15", 16", 17" and 18" in 1" increments	
Backrest Upholstery	Standard, sports adjustable	
Casters	3" Micro, (New improved high performance) 3" Micro Lighted,	
Rear Wheels	24" X/S Spoke (540mm) 24" X/S Heat Treated (559mm) 26" X/S Heat Treated (590mm) 24" Spinergy (540mm) 24" Spinergy (559mm) 26" Spinergy (590mm)	
Camber	16°, 20°	
Hand Rim Options	Anodized, painted black, plastic coat, & titanium	
Anti-tip style	Integral, quick-release, lighted	
Maximum User Weight	250 lbs (120 kg)	
Color Options	23 different colors to choose from	
Side Guards	Plastic & aluminum	
Positioning Belts	Hook & loop closure, nylon buckle closure, snowboard ratchet style pos. belt (Lg & XLg.)	
Chest Straps	2" & 4"	
5th Wheel Support	Anti-tip brace (reduces rearward seat flex) Performs like a rigid chair, while maintaining adjustability	

Specifications subject to change without notice.

New Features



Special anti-tip brace (optional)



Integral 5th wheel (anti-tip)



Ratchet-style positioning belt (optional)



Extra-hard caster wheels