


zippie®

GS™

The ideal system for  
the growing pediatric user




 **SUNRISE**  
MEDICAL®

improving people's lives

The Zippie GS, named for its built-in Growth System, meets the client's needs in both folding and rigid frames as well as standard and reverse wheel configurations.

Ability to grow the frame as the child's needs change. **The Built-in Growth System** easily "grows". The dual seat rails allows for the seating system and accessories to both be mounted without interference in the same small space.



A fully adjustable, universal axle plate allows the GS to be set up in standard or **reverse wheel configurations** and provides half-inch vertical and horizontal adjustments for increased wheel access, without sacrificing growth potential.

**Optional One-Arm Drive with Dual Handrims** is a proven design that's easy to maneuver and offers trouble-free durability.

#### ZIPPIE GS SPECIFICATIONS

	Folding	Rigid
Overall Weight:	11 kgs w/o footrest	13 kgs w/o footrest
Weight Capacity:	74 kgs	74 kgs
Overall Product Width:	19" -27"	19" -27"
Transit Option:	Yes	Yes
Front Seat Height:	14" - 20.25"	14" - 20.25"
Rear Seat Height:	12.5" - 20"	12.5" - 20"
Hanger Options:	60°, 70°, 80°, 90°(Swing-in/out & HD), ELR, ALR	60°, 70°, 80°, 90°(Swing-in/out & HD), ELR, ALR
True Fit Growth System:	Yes	Yes
Warranty:	Lifetime on frame & cross brace, 1 year on parts	Lifetime on frame & cross brace, 1 year on parts
TGA numbers:	100167	100167

\* Actual measurements may vary according to the wheelchair's configuration and accessories