EMINEOS



T: +31 (0)492 593 888 F: +31 (0)492 546 497

customerservice@sunrisemedical.nl

www.SunriseMedical.nl



01/11/2023 Version 1.0 Valid from Until 01/11/2024 Date Dealer Reference No. Quantity Delivery address

Basic configuration exclusive VAT € 1.399

- Seat width 42 cm
- Seat depth 37 cm (possible seat depth adjustments 34 till 40 cm)
- Seat height 39,7 cm
- Frame colour red
- Brake : Pull to lock, foldable
- Armrest : Height adjustable with side guards
- Velcro back angle adjustable

Angle adjustable +7 (97°)

Adjustable from 0° to +20° (90° - 110°)

Back height 55 cm, stepless adjustable 50 - 60 cm

- Push handles
- Front wheel : 5"" solid • Front Fork: 120 mm • Rear wheels 22"", 25 mm
- Anti tip RHS

The standard configuration below is indicated with 'basic'. If an alternative option is required please mark clearly by ticking the associated box.

THE CODES WHICH ARE MENTIONED BEHIND THE OPTIONS, ARE NOT ARTICLE NUMBERS. ALWAYS USE THESE CODES IN COMBINATION WITH THE PRICE LIST-ID NUMBER MENTIONED ABOVE.

CHASSIS		option ID	Price
1	1 Short frame, seat width 42 cm	Y01384	Basic
250			

WHEELS			option ID	Price
1 2 3	1	Front wheel solid	Y00829	Basic
12 (Sept.)	1	Front wheel 5""	Y00827	Basic
	2	Rear wheel 20""	Y00824	0
	2	Rear wheel 22""	Y00825	Basic
	2	Rear wheels high pressure 25 mm	Y00880	Basic
	3	Silicon cover for pushrim, set	Y00021	141

Pricelist Emineo S

BACKRESTS			option ID	Price
1 77 2 0 3	1	Velcro back, angle adj. +7 (97°) Adjustable from 0° to +20° (90°-110°) Ba	ack heigh Y00625	C
		Backrest adjustment with gas spring	Y00387	115
	2	Push bar	Y01115	C
	3	Height adjustable push handles, set	Y00124	77
ARMRESTS			option ID	Price
	∐1	Armrest height adjustable 23 - 33 cm, right	Y00075	Basic
	1	Armrest height adjustable 23-33 cm, left	Y00062	Basic
LEGRESTS			option ID	Price
LEGRESTS	1	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, right	Y00741	Basic
LEGRESTS 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	<u></u> 1	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left	Y00741 Y00740	Basi Basi
LEGRESTS	1 2	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right	Y00741 Y00740 Y00737	Basi Basi 6
LEGRESTS 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	1 2 2	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left	Y00741 Y00740	Basio Basio 62 62
LEGRESTS 2 3 6 7 6 7 7 7 7 7 7 7 7 7 7	1 2	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right	Y00741 Y00740 Y00737 Y00736	Basi Basi 62 2
LEGRESTS 2 2 3 5 6 6 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7	1 2 2 2 3	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right	Y00741 Y00740 Y00737 Y00736 Y00641	Basic
LEGRESTS	1 2 2 3 3 4	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable without heelstrap, left	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889	Basi Basi 6: 6: 2: 2: Basi Basi
LEGRESTS	1 2 2 3 3 4 4	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable without heelstrap, left Footplate angle adjustable with heelstrap, right	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892	Basi Basi 6: 6: 2: 2: Basi Basi
LEGRESTS S S 6	1 2 2 3 3 4 4 5 5	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable with heelstrap, right Footplate angle adjustable with heelstrap, right Footplate angle adjustable with heelstrap, left	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892 Y00891	Basi Basi 6: 6: 2 2 Basi Basi 4:
LEGRESTS 1 2 3 6 6 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7	1 2 2 3 3 4 4 5 5	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable without heelstrap, right Footplate angle adjustable with heelstrap, right Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle and depth adjustable without heelstrap, right	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892 Y00891 Y00894	Basic Basic 66 67 2° Basic Basic 44 44
LEGRESTS	1 2 2 3 3 4 4 5 5 6	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle and depth adjustable without heelstrap, right Footplate angle and depth adjustable without heelstrap, right Footplate angle and depth adjustable without heelstrap, left	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892 Y00891 Y00894 Y00893	Basi Basi 6: 6: 2: 2 Basi Basi 4: 4:
LEGRESTS	1 2 2 3 3 4 4 5 5	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable without heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle and depth adjustable without heelstrap, right Footplate angle and depth adjustable without heelstrap, left Footplate angle and depth adjustable without heelstrap, left Footplate angle and depth adjustable without heelstrap, left	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892 Y00891 Y00894 Y00893 Y00836	Basi Basi 6 6 2 2 Basi Basi 4 4
LEGRESTS 1 2 5 6 6 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7	1 2 2 2 3 3 3 4 4 5 5 6 6 6 7	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle and depth adjustable without heelstrap, right Footplate angle and depth adjustable without heelstrap, right Footplate angle and depth adjustable without heelstrap, left	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892 Y00891 Y00894 Y00893	Basi Basi 6 6 2 2 Basi Basi 4 4 8 8
LEGRESTS 1 2 4 5 6 7 ACCESSORIES	1 2 2 2 3 3 3 4 4 5 5 6 6 6 7	Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left Angle adjustable legrest hanger, right Angle adjustable legrest hanger, left Calf support SW 42-51 right Calf support SW 42-51 left Footplate angle adjustable without heelstrap, right Footplate angle adjustable without heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle adjustable with heelstrap, left Footplate angle and depth adjustable without heelstrap, right Footplate angle and depth adjustable without heelstrap, left Footplate angle and depth adjustable without heelstrap, left Footplate angle and depth adjustable without heelstrap, left	Y00741 Y00740 Y00737 Y00736 Y00641 Y00637 Y00890 Y00889 Y00892 Y00891 Y00894 Y00893 Y00836	Basi Basi 6 6 2 2 Basi Basi 4 4 8 8

ACCESSORIES			option ID	Price
1 2 3	1	Anti-tip stabiliser 24"" wheels, right, 20"" and 22"" for Emineo S	Y00866	Basic
	1	Anti-tip stabiliser 24"" wheels, extra (in addition), left, 20"" and 22"" for Emineo S	Y00865	35
M 18	2	Anti-tip stabiliser with Step-tip for 20"" and 22"" rear wheels, right	Y00868	73
(2) (2) (2)	2	Anti-tip stabiliser with Step-tip for 20"" and 22"" rear wheels, left	Y00867	73
	3	Side support , right	Y01434	71
4 5 6	3	Side support , left	Y01435	71
78 20 The second	4	Headrest Basic	Y00665	98
	5	Headrest Advanced	Y00667	260
	∐6	Cover headrest Basic pad	Y00648	55
	6	Cover headrest Advanced pad Seat cushion	Y00649 Y00476	84 71
	_			
REMARKS				

Sunrise Medical HCM can tailor solutions to meet your requirements. Please contact your Sunrise Medical HCM contact person for all additional wishes that are not described above.