




ASSEMBLY INSTRUCTION

Legrest hanger






Description	Article number	
Angle adjustable, round tube RHS	☐	3013501
Angle adjustable, round tube LHS	☐	3013511
Angle adjustable, square tube RHS	☐	3013502
Angle adjustable, square tube LHS	☐	3013512
Fixed angle 70°, 80°, 90°	☐	3013200





1







Sunrise Medical HCM B.V.
Vossenbeemd 104
5705 CL Helmond
The Netherlands

T: +31 (0)492 593 888
customerservice@sunrisemedical.nl
www.SunriseMedical.eu