

Technical / News Bulletin

BULLETIN DETAILS					
Product Group:	Zippie, Quickie				
Model Affected:	Q300 M Mini Kids, Q300 M Mini Teens				
Effective Date: See Below	Serial No. Effective From:	N/A			
Bulletin Subject:	Introduction of Q300 M Mini Kids & Teens				

As part of the continued expansion of the Q-Series Family, we are delighted to introduce the Q300 M Mini Kids and Q300 M Mini Teens to the portfolio.

Base Overview

The Kids and Teens seating systems are mounted to the Q300 M Mini base, offering all the benefits of indoor manoeuvrability due to a 52cm wide chassis and true mid-wheel drive technology, as well as the outdoor improvements such as greater range (36km) and increased performance via higher torque motors. Use of this chassis also means more shared components and accessories across the Q-Series fleet, reducing stock holdings and maintaining familiarity with the product.

Seating Overview

The Kids and Teens seating is based on the Kids and Teens seating from the Salsa range, meaning it is a tried and tested seating system in terms of quality and durability. However, we have introduced a number of Q-Series options and accessories to improve the clinical application of the product, such as:

- Smaller Seat Sizes
- NEW 14" Drive Wheel Option
- NEW Sedeo Lite Backrest
- Q-Series Armpads
- NEW Z-Finity Footplates
- NEW Tall Sideguards
- NEW Remote Stop Switch
- NEW Capsulated Joystick Mount





The Q300 M Mini Kids and Teens seating systems allow you to go shorter and narrower than the Salsa M2 Mini version, offering more flexibility and a wider clinical application for smaller clients

- Kids Seating: 11"-16" width, 10"-16" depth
- Teens Seating: 11"-16" width, 14"-20" depth

NEW: 14" Drive Wheels	Now available on the Q300 M Mini Kids and Teens is the option for 14" drive wheels.		
	This option significantly improves the outdoor performance of the base in terms of kerb climbing, terrain handling, ramp transitions and general comfort and stability.		
	Even better, this is a no-cost option that can be ordered straight from the configurator.		
NEW: Sedeo Lite Kids & Teens Backrest	Dedicated Sedeo Lite Deep Contour backrests in paediatric		
	sizes offer increased comfort via softer foam, integrated laterals and comfort covers.		
	As well as Introducing Paediatric sizes of both Sedeo and Whitmyer headrests directly onto the order form.		
NEW: Q-Series Armpads	We have also introduced the standard Q-Series arm pads on both the Kids and Teens seating.		
	These arm pads container softer and more comfortable foam for increased comfort as well and furthering integration with the rest of the Q-Series family		
NEW Option: Z-Finity Footplates			
	As children grow and their needs change, the Z-Finity System can be adjusted to provide proper Support of the legs and feet at all times. Its rotating hardware allows infinite adjustment of the footrests' angle and placement to achieve optimal positioning of the feet, ankles, and knees		
NEW Option: Tall Sideguards			
	We received feedback that a lot of clinicians like to utilise the sideguard for mild positioning support, not just to hold the cushion in place, which was not possible the standard Salsa sideguards.		
	For the Q300 M Mini, we have introduced a tall sideguard option to the order form which allows the sideguard to protrude above even the tallest cushions, allowing for mild positioning support.		
NEW Option: Remote Stop Switch			
	 We have introduced a new feature to the order form that was previously only available via Built4Me, the remote stop switch. If the child is driving somewhere they shouldn't, or somewhere potentially dangerous, then the switch will cut the power to the chair and not allow them to travel any further. This cannot be overridden by the child until the button is released, and as it is fitted inside the base it is also tamper proof 		



Now Available on the Configurator

Both products are also available to order directly from the configurator, rather than via a paper order for or via B4M like the Salsa version. This reduces order processing time, shows no-go's and pricing in real time, and allows you to save quotes, orders and templates.

	Sunris Medica		A Welcome Samuell	Select the model's display cuit	ture: English (UK) V	
	BASE CHAIR) MODEL (VER			CENTRE MOUNT LEGRESTS	Outstanding Selections Seat Width Seat Depth	A
					Back Frame Size Seat Cushion Options	
* Model	Q200 M	Seat Width	380 mm (38 cm)		Seat Cushion Type Seat Rail	
* Seat Depth		SW40 SW42	400 mm (40 cm) 420 mm (42 cm)		Seat Heights Seating Option	
5041 50435 5046	410 mm (41 cm) 435 mm (43.5 cm) 460 mm (46 cm)	5W44 5W46 5W43 5W50	440 mm (44 cm) 460 mm (46 cm) 480 mm (48 cm) 500 mm (50 cm)		pering upper Backness Frame Type Backness Cushion Type Backness Spile - See Jay section for (3 Backness op	
SD485 5061 50535	485 mm (48.5 cm) 510 mm (51 cm) 535 mm (53.5 cm)	* Country		Ð	Backness Height Backness Width jay Backness	
Back Frame	560 mm (56 cm)	QML170111 QML170101 QML170114	International United Kingdom Republic of Ireland		jay backnest Arminest Type Arm Pad Right	
QML030154	Backrest, Width-Adjustable, 305 - 405 mm Backrest, Width-Adjustable, 460 - 510				Arm Pad Left Centre Mount Legrest Legrest Combination	
QML020155	mm		Déll		Swine Augu Formianes	
					_	

Next Steps

To arrange a demonstration of these products please contact your Area Sales Manager.

Alternative, if you have a client in mind that would benefit from the features mentioned above, again please contact your Area Sales Manager who will be able to support an assessment in the field

For more information or for any Demo Stock requests please contact your local Area Sales Manager.

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 1-1	+44 (0)845 605 66 88 Option 1-2	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk